



Solid Phase Extraction (SPE)

Basic principles of SPE	2
CHROMABOND® phases for SPE	
Summary	7
Description of individual phases and products	10
SPE accessories	
Vacuum manifolds and accessories	50
Empty columns and accessories	52
High-throughput SPE	
Products for automated SPE	53
CHROMABOND® MULTI 96	54
Liquid-liquid extraction	
CHROMABOND® XTR	56
Phase separation	
CHROMABOND® PTS and PTL	58
Flash chromatography	
Glass columns for low pressure flash chromatography	59
CHROMABOND® Flash RS cartridges	60
Adsorbents for CHROMABOND® Flash cartridges	64

Sample clarification

CHROMAFIL® syringe filters	66
CHROMAFIL® filtration cartridges	75
CHROMAFIL® MULTI 96 filter plates	76

Vials, closures, and accessories

General information on vials and caps	77
Vials, closures, and accessories	78
Autosampler compatibility	103



Basic principles of SPE




Solid phase extraction (SPE) is a powerful method for sample preparation and is used by most chromatographers today.

About 25 years ago MACHEREY-NAGEL designed and introduced CHROMABOND® SPE cartridges containing silica-based adsorbents. Since then we developed the widest range of phases and products for SPE based on silica and polymeric materials.

SPE has capabilities in a broad range of applications:

- ✦ Environmental analyses
- ✦ Pharmaceutical and biochemical analyses
- ✦ Organic chemistry
- ✦ Food analysis

Solid Phase Extraction




SPE is a form of digital (step-wise) chromatography designed to extract, partition, and/or adsorb one or more components from a liquid phase (sample) onto a stationary phase (adsorbent or resin). An adsorbed substance can be removed from the adsorbent by step-wise increase of elution strength of the eluent (step gradient technique). SPE extends a chromatographic system's lifetime, improves qualitative and quantitative analysis, and the demand placed on an analytical instrument is considerably lessened.


In general, SPE is used for three important purposes in state-of-the-art analyses:

- ✦ Concentration of the analyte (up to factor 10.000 – increase of chromatographic sensibility and improved limits of detection)
- ✦ Removal of interfering compounds (protection of subsequent analyses like HPLC, GC, TLC, UV or IR spectroscopy, ...)
- ✦ Changing an analyte's environment to a simpler matrix more suitable for subsequent analyses

Advantages of SPE compared to classical liquid-liquid extraction:

- ✦ Lower consumption of solvents
- ✦ Faster – enormous time savings
- ✦ Lower costs per sample
- ✦ Potential for automation
- ✦ High consistency in individual sample handling
- ✦ More specific selectivity because of the broad range of adsorbents and different retention mechanisms
- ✦ Optimization of extraction by variation or adjusting of the solid phase and chromatographic conditions





Since analytes can be either adsorbed on the SPE packing material or directly flow through while the interfering substances are retained, two general separation procedures are possible – both cases are shown in the figure above.

Main steps of the SPE procedure

1. Conditioning of the adsorbent

Conditioning of the adsorbent is necessary in order to ensure reproducible interaction with the analyte. Conditioning, also called solvation, results in a wetting of the adsorbent and thus produces an environment, which is suitable for adsorption of the analyte. Non-polar adsorbents are usually conditioned with 2–3 column volumes of a solvent, which is miscible with water (methanol, THF, 2-propanol etc.), followed by the solvent in which the analyte is dissolved (pure matrix, e.g., water, buffer). Polar adsorbents are conditioned with nonpolar solvents.

After the conditioning step the adsorbent bed **must not run dry**, because otherwise solvation is destroyed (de-conditioning).

2. Sample application (adsorption)

Sample application can be performed with positive or negative pressure with a flow rate of ~3 mL/min. Sample volumes vary from a few mL up to liters.

3. Washing of the adsorbent

Washing of the adsorbent is usually achieved with a special wash solution; however, in some cases it may not be necessary. If the polarity difference between wash solution and eluent is very large, or if both are not miscible, drying of the adsorbent bed after washing is recommended to improve elution and recovery.

4. Elution

Elution with a suitable eluent should not be too fast. The elution speed depends on the column or cartridge dimension and the quantity of adsorbent (about 1 mL/min).



Basic principles of SPE

Molecular interactions in SPE


SPE adsorbents are most commonly categorized by the nature of their primary interaction mechanism with the analyte of interest. The three most common extraction

mechanisms used in SPE are reversed phase (RP), normal phase (NP) and ion exchange.

Typical extraction mechanisms


- ◆ **Reversed Phase** Extraction of hydrophobic or polar organic analytes from aqueous matrix
- ◆ **Normal Phase** Extraction of polar analytes from non-polar organic solvents
- ◆ **Ion Exchange** Extraction of charged analytes from aqueous or non-polar organic samples

Types of retention mechanisms:




Nonpolar interactions

Silica-based: C₁₈ ec, C₁₈, C₁₈ Hydra, C₈
Polymer-based: HR-X, HR-P, Easy, PS-RP
Interactions: hydrophobic
Sample: mostly aqueous
Elution: solvents with lower polarity (compared to water) CH₃OH, CH₂Cl₂, CHCl₃, hexane




Polar interactions

Silica-based: SiOH, CN, NH₂, OH (diol), C₆H₅
Other: Alox, Florisil®
Interactions: hydrogen bonds, dipole-dipole and π-π interactions
Sample: mostly organic
Elution: polar solvents (compared to sample solvent), e.g., (nonprotic) ethers, ketones (MTBE, THF, acetone) CH₂Cl₂, CHCl₃



Cation exchangers

Silica-based: SA (SCX), PCA (WCX), PSA,
Polymer-based: HR-XC, HR-XCW, PS-H⁺
Interaction: between charged analytes and functional group of cation exchanger
Sample: aqueous (pH 3-5)
Elution: acidic: pH 2 (e.g., HCl, or 20% AcOH in CH₃OH - CH₃CN)
 basic: pH 8-9 (e.g., 5% NH₃ in CH₃OH - CH₃CN)
 solvents or buffers with higher ionic strength and counter ions with high selectivity (e.g., Ca²⁺)



Anion exchangers

Silica-based: SB (SAX), NH₂, DMA
Polymer-based: HR-XA, HR-XAW, PS-OH⁻
Interaction: between charged analytes and functional group of anion exchanger
Sample: aqueous (pH 8-9)
Elution: basic: pH 10 (e.g., 20% NH₃ in CH₃OH - CH₃CN)
 acidic: pH 4-5 (e.g., HCl, or 5% AcOH in CH₃OH - CH₃CN)
 solvents or buffers with higher ionic strength and counter ions with high selectivity (e.g., citrate)

It should be noted, that in SPE the interactions described above are not found in pure form, but in combination. For example, modified silicas, unless they have been subjected to endcapping (silanization of residual silanol groups with short-chain silanes), still possess free silanol groups, which can enter into secondary interactions.



Sample pretreatment

For direct extraction with adsorbents the sample matrix (sample environment) has to fulfill three conditions:

- The matrix has to be liquid, if possible with low viscosity.
- Solids should be removed from the liquid matrix.
- The matrix (sample environment) should be suitable for retention of the analyte.

For solid samples there are different methods to convert the sample into a suitable matrix:

- Dissolution of the solid sample in a suitable solvent
- Lyophilization of the sample and dissolution in a suitable solvent
- Extraction of the solid sample with a suitable solvent
- Homogenization of the sample in a suitable solvent

In order to find the suitable solvent, one has to consider all desired sample components. Also, the suitable solvent should enhance retention of the analyte. For example, samples with large contents of solids are often homogenized in nonpolar solvents like hexane, while for samples with high water content dissolution in acids, bases, buffers or very polar solvents such as methanol is recommended.

Additionally, SPE allows to alter the properties of the sample matrix. If, for example, natural products are extracted with methanol or acetone, the polarity of the extracts can be increased by dilution with water, in order to enhance nonpolar solid phase extraction on the C₁₈ material.

Our CHROMABOND® QC policy

- **Highest production standard**
our facilities are EN ISO 9001:2008 certified
- All of our bonded phases and SPE products are vigorously tested for perfect **reproducibility** from lot-to-lot and within every single batch
Careful attention to particle size distribution and pore diameters assures consistent column flow
Chemical reproducibility is guaranteed by strict quality control throughout manufacturing
- All products are individually tested to meet our **strict quality specifications**, ensuring our outstanding product reproducibility, reliability and performance
- Each product is supplied with a **certificate of analysis** stating the results of internal examinations and quality control






Basic principles of SPE


CHROMABOND®

Sample reservoir
polypropylene or glass




CHROMAFIX®

Luer fitting, female




CHROMABOND® LV

Sample reservoir
polypropylene, 15 mL



CHROMABOND® MULTI 96



Design of columns, cartridges and 96-well plates

All CHROMABOND® columns, cartridges and 96-well plates are manufactured from polypropylene (PP) with lowest content of extractables (plasticizers, stabilizers, ...) offering blank value free results when using most common solvents. The high quality CHROMABOND® adsorbents are kept in place by chemically very inert polyethylene filter elements (PE, standard pore size 20 µm).

CHROMABOND® polypropylene columns

- PP columns with PE filter elements
- Different sizes from 1, 3, 6 up to 150 mL
- Adsorbent weights from 20 mg to 50 g
- Male Luer tip as exit
- Compatible with most robots (e.g., Gilson ASPEC™, Caliper AutoTrace®)

CHROMABOND® glass columns

- Glass columns with chemically very inert glass fiber filter elements (nominal pore size 1 µm)
- Two different sizes: 3 and 6 mL
- Available with all CHROMABOND® phases
- Excludes any influence from the column material (e.g., plasticizers)

CHROMAFIX® cartridges

- PP cartridges with PE filter elements
- Three different sizes with different adsorbent weights: Small (0.4 mL), Medium (0.8 mL), Large (1.8 mL)
- Female Luer tip at the inlet, male Luer tip as exit
- Offers alternative way of handling using positive pressure by syringes or peristaltic pumps
- Especially suited for convenient solid phase extraction of small sample volumes

CHROMABOND® LV columns

- Large volume PP columns with PE filter elements
- Three different adsorbent weights (100, 200 and 500 mg)
- Funnel-shaped reservoir with 15 mL volume
- Especially for clinical samples – the whole sample (e.g., urine, serum, blood) can be applied to the column in one step
- Can be directly used in the Zymate® lab robots of Zymark

CHROMABOND® MULTI 96 · SPE in 96-well format

- 96-well PP plates with PE filter elements
- Cavity volume 1.5 mL
- Adsorbent weights 10, 25, 50 and 100 mg
- Supplied with any CHROMABOND® SPE adsorbents
- For simultaneous preparation of 96 samples
- Easy method transfer from CHROMABOND® columns or CHROMAFIX® cartridges to CHROMABOND® MULTI 96
- Readily adaptable to all common automated / robotic handling systems (for details see page 54)

Online-SPE (see page 53)

- Online columns and cartridges
- SPE columns with caps and needles for the Gerstel MultiPurposeSampler (MPS)
- Columns for Gilson ASPEC™ systems (ASP)



For the development kits as well as for all individual CHROMABOND®, CHROMABOND® LV and CHROMAFIX® types columns are sealed in units of five columns each to prevent adsorption of contaminants from the environment, e.g., laboratory air.



Ordering information

Designation	Contents of the kit	REF
Investigating the best separation mechanism for a clean-up procedure		
CHROMABOND® HR-Xpert development kit I	columns with 3 mL, 60 mg each (particle size 45 µm): 10 columns with HR-X; 5 columns each with HR-XC, HR-XA, HR-XCW, HR-XAW	730723
CHROMABOND® HR-Xpert development kit II	columns with 3 mL, 200 mg each (particle size 85 µm): 10 columns with HR-X; 5 columns each with HR-XC, HR-XA, HR-XCW, HR-XAW	730726
CHROMABOND® polymer development kit	5 columns each with 3 mL, 200 mg: HR-X, HR-XC (MCX), HR-XA (MAX), HR-P, Easy, PS-H ⁺ , PS-OH ⁻	730288
CHROMABOND® standard development kit	5 columns each with 3 mL, 500 mg: C ₁₈ , C ₁₈ ec, C ₈ , C ₆ H ₅ , NH ₂ , DMA, OH, CN, SiOH, SA (SCX), SB (SAX)	730496
Selecting the optimum RP phase for a clean-up procedure		
CHROMABOND® RP development kit I	10 columns each with 3 mL, 500 mg: C ₁₈ , C ₁₈ ec, C ₈ , C ₄ and 10 columns each with 3 mL, 200 mg HR-P, HR-X	730197
CHROMABOND® RP development kit II	10 columns each with 1 mL, 100 mg: C ₁₈ , C ₁₈ ec, C ₈ , C ₄ , HR-P, HR-X	730207
CHROMAFIX® RP development kit I	10 cartridges each CHROMAFIX® S: C ₁₈ , C ₁₈ ec, C ₈ , C ₄ , HR-P, HR-X	731883
CHROMABOND® RP development kit III	10 columns each with 3 mL, 500 mg: C ₁₈ , C ₁₈ ec, C ₁₈ Hydra, C ₈ and 10 columns each with 3 mL, 200 mg HR-P, HR-X	730490
CHROMABOND® RP development kit IV	10 columns each with 1 mL, 100 mg: C ₁₈ , C ₁₈ ec, C ₁₈ Hydra, C ₈ , HR-P, HR-X	730491
CHROMAFIX® RP development kit II	10 cartridges each CHROMAFIX® S: C ₁₈ , C ₁₈ ec, C ₁₈ Hydra, C ₈ , HR-P, HR-X	731886
CHROMABOND® RP development kit V	10 columns each with 3 mL, 500 mg: C ₆ H ₅ , NO ₂ , C ₆ H ₁₁ ec, C ₄ , C ₂	730492
CHROMABOND® RP development kit VI	10 columns each with 1 mL, 100 mg: C ₆ H ₅ , NO ₂ , C ₆ H ₁₁ ec, C ₄ , C ₂	730493
CHROMAFIX® RP development kit III	10 cartridges each CHROMAFIX® S: C ₆ H ₅ , NO ₂ , C ₆ H ₁₁ ec, C ₄ , C ₂	731887
Selecting the optimum polar phase for a clean-up procedure		
CHROMABOND® polar development kit I	10 columns each with 3 mL, 500 mg: SiOH, Florisil®, NH ₂ , CN, OH	730199
CHROMABOND® polar development kit II	10 columns each with 1 mL, 100 mg: SiOH, Florisil®, NH ₂ , CN, OH	730208
CHROMAFIX® polar development kit	10 cartridges each CHROMAFIX® S: SiOH, Florisil®, NH ₂ , CN, OH	731884
Selecting the optimum ion exchanger for a clean-up procedure		
CHROMABOND® ion exchange development kit I	10 columns each with 3 mL, 500 mg: SA (SCX), SB (SAX), HR-XC (MCX), HR-XA (MAX), PS-OH ⁻ , PS-H ⁺ , DMA	730206
CHROMABOND® ion exchange development kit II	10 columns each with 1 mL, 100 mg: SA (SCX), SB (SAX), HR-XC (MCX), HR-XA (MAX), PS-OH ⁻ , PS-H ⁺ , DMA	730209
CHROMAFIX® ion exchange development kit I	10 cartridges each CHROMAFIX® S: SA (SCX), SB (SAX), HR-XC (MCX), HR-XA (MAX), PS-OH ⁻ , PS-H ⁺ , DMA	731885
CHROMABOND® cation exchange development kit I	10 columns each with 3 mL, 500 mg: SA (SCX), PSA, PCA, HR-XC (MCX), HR-XCW (WCX), PS-H ⁺	730494
CHROMAFIX® cation exchange development kit	10 cartridges each CHROMAFIX® S: SA (SCX), PSA, PCA, HR-XC (MCX), HR-XCW (WCX), PS-H ⁺	731888
Phase selection for clean-up procedures for environmental samples		
CHROMABOND® kit I for environmental sample preparation	10 columns each with 3 mL, 200 mg HR-P, 6 mL, 1000 mg C ₁₈ ec, 6 mL, 2000 mg C ₁₈ PAH, 6 mL, 500/1000 mg CN/SiOH, 3 mL, 500/500 mg SA/SiOH	730205
CHROMABOND® kit II for environmental sample preparation	5 columns each with 3 mL, 500/500 mg SiOH-H ₂ SO ₄ /SA, 3 mL, 500 mg SiOH, 6 mL, 1000 mg Florisil, 3 mL, 500/500 mg SA/SiOH, 6 mL, 700/2000/700 mg NAN	730349



Summary of MN phases for SPE

Code	Matrix	Modification / Application	Similar phases*	Page
Reversed phases				
HR-X	PS/DVB		ENVI-Chrom P · Strata™-X · Oasis® HLB · Nexus	14
Easy	PS/DVB	polar, bifunctional	Strata™-X · Oasis® HLB · Porapak™ RDX · Nexus, Bond Elut® PPL, Focus™ · Styre Screen® DVB Bakerbond™ H ₂ O-philic DVB · Isolute® ENV+	20
HR-P	PS/DVB		Strata™ SDB-L · Bond Elut® ENV, Bond Elut® LMS · DCS-PS/DVB, ENV PS-DVB · Bakerbond™ H ₂ O-phobic DVB · Isolute® 101 · LiChrolut® EN	21
PS-RP	PS/DVB	removal of organic components	like HR-P	22
C ₁₈ ec	silica	octadecyl, endcapped	Strata™ C18-E · Sep-Pak® tC18 · Bond Elut® C18 · DSC-18(Lt), ENVI-18, LC-18 · CLEAN-UP® C18, Bakerbond® Octadecyl · Isolute® C18(EC), LiChrolut® RP-18 E	23
C ₁₈ ec f	silica	as above, fast flow		23
C ₁₈	silica	octadecyl, not endcapped	Strata™ C18-U · AccuBond® C18 · Bakerbond™ PolarPlus · Isolute® C18 · LiChrolut® RP-18	24
C ₁₈ f	silica	as above, fast flow		24
C ₁₈ PAH	silica	special octadecyl phase, for enrichment of PAHs from water	Bakerbond™ Octadecyl Lightload	42
C ₁₈ Hydra	silica	octadecyl, not endcapped, for polar analytes		25
C ₈	silica	octyl	Strata™ C8 · Sep-Pak® C8 · Bond Elut® C8 · DSC-8, ENVI-8, LC-8 · CLEAN-UP® C8 · AccuBond® C8 · Bakerbond™ Octyl · Isolute® C8(EC)	26
C ₄	silica	butyl		27
C ₂	silica	dimethyl	Bond Elut® C2	27
C ₆ H ₁₁ ec	silica	cyclohexyl, endcapped	Bond Elut® CH	28
C ₆ H ₅	silica	phenyl	Strata™ PH · Bond Elut® PH · DSC-Ph · CLEAN-UP® Phenyl · AccuBond® Phenyl · Bakerbond™ Phenyl · Isolute PH(EC)	29
Normal phases				
SiOH	silica	unmodified	Strata™ Si-1 · Bond Elut® silica · DSC-Si, LC-Si · CLEAN-UP® silica · Accubond® silica, Bakerbond™ silica gel · Isolute® silica · LiChrolut® Si	32
NH ₂	silica	aminopropyl	Strata™ NH ₂ · Sep-Pak® NH ₂ · Bond Elut NH ₂ · DSC-NH ₂ , LC-NH ₂ · CLEAN-UP® aminopropyl · Accubond® NH ₂ · Bakerbond™ amino · Isolute® NH ₂ · LiChrolut® NH ₂	31
OH (Diol)	silica	diol	DSC-Diol, LC-Diol · Accubond® Diol (OH)	30
CN	silica	cyano	Strata™ CN · Sep-Pak® CN · Bond Elut® CN-U · DSC-CN, LC-CN · CLEAN-UP® CN · Accubond® CN · Bakerbond™ cyano · Isolute® CN · LiChrolut® CN	30
Alox A	aluminium oxide acidic		LC-Alumina-A · Accubond® aluminium oxide A	33
Alox N	aluminium oxide neutral		LC-Alumina-N · Accubond® aluminium oxide N	33
Alox B	aluminium oxide basic		LC-Alumina-B · Accubond® aluminium oxide B	33
Florisil®	magnesium silicate		Strata™ FL-PR · Sep-Pak® Florisil® · Bond Elut® Florisil® · ENVI-Florisil®, LC-Florisil® · CLEAN-UP® Florisil® · Accubond® Florisil® · Bakerbond™ Florisil® · Isolute® FL · LiChrolut® Florisil®	34
PA	polyamide 6		DPA-6S	34
Ion exchangers				
SB	silica	quaternary ammonium anion exchanger (SAX)	Strata™ SAX, Sep-Pak® SAX, Bond Elut® SAX · DSC-SAX, LC-SAX · CLEAN-UP® Quaternary Amine · Accubond® SAX · Bakerbond™ Quaternary Amine · Isolute® SAX · LiChrolut® SAX	37
* Phases which provide a similar selectivity based on chemical or physical properties (list not complete)				



Code	Matrix	Modification / Application	Similar phases*	Page
SA	silica	benzenesulfonic acid cation exchanger (SCX)	Strata™ SCX · Bond Elut® SCX · DSC-SCX, LC-SCX · CLEAN-UP® Benzenesulfonic Acid · Accubond® SCX · Bakerbond™ Aromatic Sulfonic Acid · Isolute® SCX · LiChrolut® SCX	36
PCA	silica	propylcarboxylic acid cation exchanger (WCX)	Strata™ WCX · Bond Elut® CBA · DSC-WCX, LC-WCX · CLEAN-UP® Carboxylic Acid · Bakerbond™ Carboxylic Acid · Isolute® CBA	35
PSA	silica	propylsulfonic acid cation exchanger		35
HR-XC	PS/DVB	strong mixed mode cation exchanger for basic analytes (MCX)	Oasis® MCX · Strata™-X-C · HyperSep™ Retain™-CX · Styre Screen® DBX	16
HR-XA	PS/DVB	strong mixed mode anion exchanger for acidic analytes (MAX)	Oasis® MAX · Strata™-X-A · HyperSep™ Retain™-AX · Styre Screen® QAX	17
HR-XCW	PS/DVB	weak mixed mode cation exchanger for basic analytes (WCX)	Oasis® WCX · Strata™ X-CW	18
HR-XAW	PS/DVB	weak mixed mode anion exchanger for acidic analytes (WAX)	Oasis® WAX · Strata™ X-AW	19
PS-OH ⁻	PS/DVB	strong anion exchanger in OH ⁻ form		22
PS-H ⁺	PS/DVB	strong cation exchanger in H ⁺ form		22
PS-Mix	PS/DVB	mixture of PS-OH ⁻ and PS-H ⁺		22
PS-Ag ⁺	PS/DVB	strong cation exchanger in Ag ⁺ form		22
PS-Ba ²⁺	PS/DVB	strong cation exchanger in Ba ²⁺ form		22
Phases for special applications				
Dry	Na ₂ SO ₄	for drying organic samples		47
Drug	silica	bifunctional C ₈ /SA, for enrichment of drugs from urine	Strata™ Screen-C · Bond Elut® Certify I · DSC-MCAX · Clean Screen® DAU · Accubond® Evidex · Bakerbond™ Narc-2 · Isolute® HCX · LiChrolut® TSC · HyperSep™ Verify CX	38
Drug II	silica	bifunctional C ₈ /SB, for extraction of THC and derivatives and of acidic analytes from biological fluids	Strata™ Screen-A · Bond Elut Certify II · Clean Screen® THC · Bakerbond® Narc-1 · Isolute® HAX · HyperSep™ Verify AX	39
Crosslinks	cellulose	for enrichment of collagen crosslinks		40
Tetracycline	silica	special octadecyl phase, for enrichment of tetracyclines		40
HR-P-AOX	PS/DVB	for extraction of AOX from water (DIN 38409 - H22)		41
CN/SiOH	silica	combination phase for enrichment of PAHs from soil		44
NH ₂ /C ₁₈	silica	combination phase for enrichment of PAHs from water		42
Na ₂ SO ₄ /Florisil®		combination phase for extraction of hydrocarbons from water (DIN H-53 / ISO DIS 9377-4)		43
SA/SiOH	silica	combination phase for enrichment of PCB from waste oil	Bakerbond™ PCB-N	45
SiOH-H ₂ SO ₄ /SA	silica	combination phase, used together with SiOH for enrichment of PCB from oil		46
NAN	silica / AgNO ₃ + Na ₂ SO ₄	combination phase for enrichment of PCB from sludge		44
ABC18	silica	octadecyl, with ion exchange functions, for acrylamide analysis	Isolute® M-M	47
Diamino	silica	primary and secondary amine functions (PSA), for determination of pesticides in food samples (QuEChERS method)	Supelclean™ PSA, Bond Elut® PSA	48
Phase separation		CHROMABOND® PTL/PTS		58
Liquid-liquid extraction		CHROMABOND® XTR	EXTrelut® · Chem Elut™ · Hydromatrix™	56

* Phases which provide a similar selectivity based on chemical or physical properties (list not complete)





The professional concept of innovative SPE phases

The **CHROMABOND® HR-Xpert** family comprises 5 polymer-based RP and mixed-mode ion exchange phases:

- **CHROMABOND® HR-X** hydrophobic PS/DVB copolymer
- **CHROMABOND® HR-XC** strong mixed-mode cation exchanger
- **CHROMABOND® HR-XA** strong mixed-mode anion exchanger
- **CHROMABOND® HR-XCW** weak mixed-mode cation exchanger
- **CHROMABOND® HR-XAW** weak mixed-mode anion exchanger

These innovative SPE phases offer

- **State-of-the-art spherical polymer**
 - Two particle sizes (45 µm and 85 µm) adequate for different sample volumes and matrices
 - Broad spectrum of application with special suitability for enrichment of pharmaceuticals from biological matrices
 - Ideal flow properties due to low content of particulate matter
- **Optimized pore structure and high specific surface**
 - High loadability and outstanding elution properties
 - Low solvent consumption
 - Rapid, economical analyses
- **High-purity adsorbent material**
 - Allows highest reproducibility with extremely low blind values
 - Reliable analyses at ultra trace level
 - No method adaptation for new batches necessary




The HR-Xpert concept guarantees:

- RP and mixed-mode SPE phases with distinct ion exchange and reversed phase properties: excellent enrichment of neutral, acidic and basic compounds
- Modern, spherical support polymer with optimized pore structure and high surface: good reproducibility, reliable and cost-efficient analysis
- Possibility for more aggressive washing procedures for matrix removal: cleaner samples and protection of your HPLC and GC instruments
- Quantification of analytes also from heavily contaminated samples: lower limits of detection also for critical matrices

CHROMABOND® HR-Xpert is the perfect combination for all tasks in sample preparation



Chemical structures of the phases:



Similar phases:


CHROMABOND® HR-X:	Oasis® HLB, Strata™-X, Nexus, ENVI-Chrom P
CHROMABOND® HR-XC:	Oasis® MCX, Strata™-X-C, HyperSep™ Retain™-CX, StyreScreen® DBX
CHROMABOND® HR-XA:	Oasis® MAX, Strata™-X-A, HyperSep™ Retain™-AX, StyreScreen® QAX
CHROMABOND® HR-XCW:	Oasis® WCX, Strata™-X-CW
CHROMABOND® HR-XAW:	Oasis® WAX, Strata™-X-AW



CHROMABOND® HR-Xpert


2 particle sizes – 1 goal: HR-Xpert for optimized sample preparation

For different application requirements the particle sizes complement each other perfectly.



Ideal for:

- + smaller sample volumes
- + smaller adsorbent weights
- + lower elution volumes



Recommended for:


- + large volume or viscous samples, heavy matrix load
- + operation without vacuum possible (e.g., for volatile analytes)
- + higher adsorbent weight without increase in back pressure

Features of 45 µm particles

- About half the radius results in 8fold particle number per volume for approx. equal adsorbent weight
- Same specific surface for both particle sizes:
→ considerably larger freely accessible external surface for 45 µm particles
- Denser adsorbent packing:
→ enhanced interaction of the analyte with the adsorbent, better extraction results


Ideal elution characteristics

Method: 1 mL column with 30 mg CHROMABOND® HR-X, 1 mL standard solution (1 mg/mL hexobarbital), dry-ing, elution in portions of 100 µL with methanol (see application 305490 at www.mn-net.com/apps)



Advantages of 45 µm particles:

- + faster elution
- + lower elution volumes required



Breakthrough behavior in enrichment

Method: 1 mL column with 15 mg CHROMABOND® HR-X, apply portions of 1 mL standard solution (250 µg/mL hexobarbital in water), collect eluates (see application 305480 at www.mn-net.com)


45 µm (red)

The analyte is completely retained up to fraction 10.

85 µm (blue)

Small amounts even break through with fraction 4.

45 µm particles provide better enrichment and breakthrough behavior for small adsorbent weights. When using larger adsorbent weights this effect is less pronounced, since then analytes have sufficient contact with the 85 µm adsorbent particles as well.




45 µm particles are ideal for small sample and elution volumes, while for large amounts of sample and adsorbent 85 µm particles show advantages due to better flow properties.



The CHROMABOND® HR-Xpert concept for neutral, acidic and basic analytes

3 paths – 1 goal: cleaner samples

Depending on the character of the analytes HR-Xpert offers suitable adsorbents and optimal methods for sample preparation, cleaning and concentration.



* Under organic washing and elution conditions the following compounds will be also eluted:

HR-X:	polar compounds such as organic acids and bases
HR-XC, HR-XCW:	acidic components and impurities
HR-XA, HR-XAW:	basic components and impurities



Polymer-based reversed phases for SPE

HR-X spherical, hydrophobic polystyrene-divinylbenzene adsorbent resin

- Hydrophobic polystyrene-divinylbenzene copolymer
- pH stability 1–14

High-purity material with highest reproducibility and lowest blank values due to an optimized manufacturing process

Spherical particles, size 45 µm and 85 µm (standard)
pore size 55–60 Å; very high surface 1000 m²/g
capacity 390 mg/g (caffeine in water)

Excellent recovery rates especially for the enrichment of pharmaceuticals and active ingredients due to the spherical structure of the particles, very homogeneous surface, and optimized pore structure

- Recommended application:

Pharmaceuticals / active ingredients from tablets, creams and water / waste water

Drugs and pharmaceuticals from urine, blood, serum and plasma

Trace analysis of pesticides, herbicides, phenols, PAHs and PCBs from water

Drugs from water

Column type:
CHROMABOND® HR-X, 3 mL, 200 mg
REF 730931

Sample: 1 µg/mL each in water

Column conditioning: 5 mL methanol, 5 mL dist. water

Sample application: slowly aspirate 500 mL water (pH 3) through the column

Column washing: 5 mL water

Elution: after drying 3 x 2 mL acetonitrile

Further analysis: HPLC on NUCLEODUR® C₁₈ Gravity, 5 µm;
see MN Appl. No. 121690

Recovery rate [%]

Compound	HR-X	Strata™ X
Ketoprofen	98	92
Ibuprofen	91	93
Pentobarbital	99	95
Meclofenamic acid	92	93
Protriptyline	63	45
Nortriptyline	53	39

MN Appl. No. 304240

Pesticides from water

Column type:
CHROMABOND® HR-X, 3 mL, 200 mg
REF 730931

Sample pretreatment: samples are spiked with 500 ng of each pesticide in 1000 mL water, adjusted to pH 2 with HCl or pH 7

Column conditioning: 10 mL methanol, 10 mL dist. water

Sample application:

slowly pass 1000 mL spiked water sample through the column with the aid of a tubing adapter (REF 730243)

Elution: after drying 5 mL methanol – THF (1:1, v/v)

Further analysis: HPLC

Recovery rates [%]

Compound	HR-X pH 2	Compound	HR-X pH 7
Metamitron	86	Desisopropylatrazine	90
Quinmerac	90	2,4-Dichlorobenzamide	95
Chloridazon	93	Desethylatrazine	89
Picloram	83	Hexazinone	95
Metribuzin	84	Bromacil	103
Cyanazine	83	Simazine	91
Metabenzthiazuron	94	Desethylterbutylazine	89
Chlortoluron	91	Atrazine	88
Isoproturon	89	Metaxyl	97
Diuron	91	Metazachlor	93
Dimethenamid-P	89	Propazine	88
Linuron	94	Terbutylazine	86
Epoxyconazole	85	Metolachlor	97
Penconazole	90		
Alachlor	93		
Propiconazole-1	89		
Flufenacet	91		
Diffenencam	58		
Triallate	42		

MN Appl. No. 304250/304260



Options for online-SPE and automated SPE see page 53




Highest reproducibility

- ✓ within each batch ✓ from batch to batch

Compounds:

- A phenobarbital
B pentobarbital
C hexobarbital



Barbiturates from serum

Column type:
CHROMABOND® HR-X, 3 mL, 200 mg
REF 730931

Sample: 100 ng/mL each in serum

Column conditioning: 5 mL methanol, 5 mL dist. water

Sample application: 1 mL spiked serum

Column washing: 5 mL water

Elution: after drying 3 x 2 mL methanol

Further analysis: HPLC on NUCLEODUR® 100-5 C₁₈ ec, see MN Appl. No. 117820

MN Appl. No. 304290

Standard protocol for CHROMABOND® HR-X

Column type:
CHROMABOND® HR-X, 3 mL, 200 mg
REF 730931

Sample pretreatment: if necessary, adjust pH value

Column conditioning: 5 mL methanol

Equilibration: 5 mL water

Sample application:

slowly aspirate the sample through the column

Column washing: 5 mL water – methanol (95:5, v/v)





Elution: after drying 3 x 2 mL methanol

Further analysis: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC



MN Appl. No. 304310

Ordering information

	Volume	Adsorbent weight					Pack of	
	CHROMABOND® HR-X polypropylene columns (85 µm)							
		30 mg	60 mg	100 mg	200 mg	500 mg	1 g	
	1 mL	730934		730935			30	
	3 mL		730936		730931	730937	30	
	6 mL				730938	730939	30	
	15 mL					730940	730941	20
	CHROMABOND® HR-X polypropylene columns (85 µm) • BIGpacks							
				200 mg	500 mg			
3 mL				730931.250		250		
6 mL				730938.250	730939.250	250		
CHROMABOND® HR-X polypropylene columns (45 µm) • NEW!								
	30 mg	60 mg	100 mg	200 mg				
1 mL	730934P45		730935P45			30		
3 mL		730936P45		730931P45		30		
	CHROMABOND® LV-HR-X (85 µm)							
		30 mg	60 mg		200 mg			
	15 mL	732130	732131		732132		30	
	CHROMABOND® MULTI 96 HR-X							
		96 x 10 mg (45 µm)	96 x 25 mg (45 µm)		96 x 50 mg (85 µm)	96 x 100 mg (85 µm)		
		738530.010M	738530.025M		738530.050M	738530.100M	1	
	CHROMABOND® HR-X adsorbent (85 µm)							
						730663	20 g	



Polymer-based ion exchangers for SPE

HR-XC

- Strong acidic benzenesulfonic acid cation exchanger
exchange capacity 1.0 meq/g, $pK_a < 1$
Base material polystyrene-divinylbenzene copolymer
pH stability 1–14; high purity material, highest reproducibility and lowest blank values due to an optimized production process
Spherical particles, size 45 μm and 85 μm (standard); pore size 65–75 \AA
very large specific surface 800 m^2/g ; pore volume 1.4 cm^3/g
RP capacity 300 mg/g (caffeine in water)
Outstanding recovery rates especially for the enrichment of basic analytes

strong cation exchanger

- Recommended application:
Basic active ingredients from heavily matrix-contaminated samples like, e.g., urine, plasma, serum
Fungicides from food, melamine from milk
Basic analytes like, e.g., amines; bases with pK_a 2–10

Solid Phase Extraction



Standard protocol for CHROMABOND® HR-XC

- Column type:**
CHROMABOND® HR-XC, 3 mL, 200 mg
REF 730952
- Sample pretreatment:** individual sample preparation with reference to analytes and matrix
- Column conditioning:** 5 mL methanol
- Equilibration:** 5 mL water
- Sample application:** slowly aspirate sample through the column
- Washing 1:** 2 mL 0.1 mol/L HCl in water
- Washing 2/Elution 1:** 2 mL methanol (neutral and acidic compounds); if necessary, further washing steps
- Elution 2:** after drying 5 mL methanol – 5% NH_3 (basic compounds)
- Further analysis: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC

MN Appl. No. 304740

Fractionation of acidic, neutral and basic




- Column type:**
CHROMABOND® HR-XC, 3 mL, 200 mg
REF 730952
- Sample:** 1 mL spiked matrix, acidified with 200 μL 2% H_3PO_4
- Column conditioning:** 5 mL methanol, then 5 mL water
- Sample application:** slowly aspirate sample through the column
- Washing:** 2 mL 0.1 mol/L HCl
- Elution:** 2.5 mL methanol (fraction A: neutral and acidic analytes); then 5 mL methanol – NH_3 90:10, v/v (fraction B: basic analytes)
- Further analysis for fraction A: HPLC, e.g., on NUCLEODUR® C₁₈ Gravity, see MN Appl. No. 122230; for fraction B: HPLC on NUCLEODUR® C₈ Gravity, see MN Appl. No. 118520

Recovery rates [%]

Fraction A: neutral and acidic analytes		Fraction B: basic analytes		
Compound	HR-XC	Compound	HR-XC	Oasis® MCX Strata™ X-C
Suprofen	108	1. Doxepin	101	68 82
Naproxen	85	2. Imipramine	95	71 85
Tolmetin	73	3. Amitriptyline	94	72 78
Phenobarbital	108	4. Trimipramine	92	70 81
Indomethacin	33			
Hexobarbital	80			

MN Appl. No. 304780

Ordering information

Volume	Adsorbent weight						Pack of	
	CHROMABOND® HR-XC polypropylene columns (85 µm)							
		30 mg	60 mg	100 mg	150 mg	200 mg	500 mg	
	1 mL	730969		730049				30
	3 mL		730956			730952	730953	30
	6 mL				730957		730955	30
CHROMABOND® HR-XC polypropylene columns (45 µm) · <i>NEW!</i>								
	1 mL	730969P45		730049P45				30
	3 mL		730956P45			730952P45		30
	CHROMAFIX® HR-XC cartridges (85 µm)							
	Size	S		M		L		
	adsorbent weight Ø	155 mg		240 mg		500 mg		
		731755		731756		731757		30
	CHROMABOND® HR-XC adsorbent (85 µm)							
							730664	100 g




HR-XA

- Strong basic quaternary ammonium anion exchanger
exchange capacity 0.25 meq/g, pKa ~ 18
- Base material polystyrene-divinylbenzene copolymer
pH stability 1–14; high purity material with highest reproducibility and lowest blank values due to an optimized production process
- Spherical particles, size 45 µm and 85 µm (standard); pore size 55–65 Å
very large specific surface 850 m²/g; pore volume 1.4 cm³/g
- RP capacity 350 mg/g (caffeine in water)
- Outstanding recovery rates especially for the enrichment of acidic analytes

strong anion exchanger

- Recommended application:
Acidic active ingredients from heavily matrix-contaminated samples like, e.g., urine, plasma, serum
- Phenolic acids, acidic herbicides
- Weak/medium-strength acids with pKa 2–8

analytes from serum



Standard protocol for CHROMABOND® HR-XA

Column type:
CHROMABOND® HR-XA, 3 mL, 200 mg
REF 730951

Sample pretreatment: individual sample preparation with reference to analytes and matrix

Conditioning: 5 mL methanol

Equilibration: 5 mL water

Sample application: slowly aspirate sample through the column

Washing 1: 2 mL 0.1 mol/L NaOH in water

Washing 2 / Elution 1: 2 mL methanol (neutral and basic compounds), if necessary, further washing steps

Elution 2: after drying 5 mL methanol – 1 to 10 % formic acid (acidic compounds)

Further analyses: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC




MN Appl. No. 304970

For further applications on
CHROMABOND® polymer phases see
our online application database at

www.mn-net.com/apps

Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight						Pack of
	CHROMABOND® HR-XA polypropylene columns (85 µm)							
		30 mg	60 mg	100 mg	150 mg	200 mg	500 mg	
	1 mL	730968		730727				30
	3 mL		730950			730951	730954	30
	6 mL				730958		730966	30
	CHROMABOND® HR-XA polypropylene columns (45 µm) · NEW!							
	1 mL	730968P45		730727P45				30
	3 mL		730950P45			730951P45		30
	CHROMAFIX® HR-XA cartridges (85 µm)							
	Size	S		M		L		
	adsorbent weight Ø	155 mg		240 mg		500 mg		
		731768		731769		731770		50
	CHROMABOND® HR-XA adsorbent (85 µm)							
							730671	100 g



Polymer-based ion exchangers for SPE

HR-XCW

- Weak acidic carboxylic acid cation exchanger
exchange capacity >0.7 meq/g, pKa ~ 5
- Base material spherical PS/DVB copolymer, pH stability 1–14
- high purity material, highest reproducibility and lowest blank values due to an optimized production process
- Spherical particles, size 45 µm and 85 µm (standard); pore size 50–60 Å
- very large specific surface 850 m²/g; pore volume 1.2–1.4 cm³/g
- RP capacity 350 mg/g (caffeine in water)
- Outstanding recovery rates especially for enrichment of strongly basic analytes

weak cation exchanger

- Recommended application:
Basic compounds like quaternary amines
- Active ingredients from heavily matrix-contaminated samples like, e.g., urine, plasma, serum
- Strong bases with pKa > 10

Standard protocol for CHROMABOND® HR-XCW

Column type:
CHROMABOND® HR-XCW, 3 mL, 200 mg
REF 730739

Sample pretreatment: individual sample preparation with reference to analytes and matrix

Column conditioning: 5 mL methanol

Equilibration: 5 mL acidified water

Sample application: slowly aspirate sample through the column

Washing 1: 2 mL acidified water

Washing 2/Elution 1: 2 mL methanol (neutral and acidic compounds), if necessary, further washing steps

Elution 2: after drying 2 x 2 mL methanol – 1 to 5 % formic acid (strongly basic compounds)

Further analysis: if necessary, evaporate and redissolve in a suitable solvent; HPLC or GC

MN Appl. No. 305300



Analysis of perfluorinated


Column: 125 x 2 mm NUCLEODUR® Sphinx RP, 3 µm

Eluent: A) 10 mmol/L NH₄Ac in water – methanol (75:25, v/v); B) 10 mmol/L NH₄Ac in acetonitrile – methanol (75:25, v/v)
10–30 % B in 3 min, 30–55 % B in 8 min, 55–10 % B in 4 min

Flow rate: 0.30 mL/min, temperature 50 °C

Injection: 2.5 µL (5 mg/L each after SPE enrichment)





Detection: MS, ESI negative



MN Appl. No. 123340

Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight					Pack of	
	CHROMABOND® HR-XCW polypropylene columns (85 µm)							
		30 mg	60 mg	100 mg	150 mg	200 mg	500 mg	
	1 mL	730731		730733				30
	3 mL		730735			730739	730741	30
	6 mL				730737		730743	30
	CHROMABOND® HR-XCW polypropylene columns (45 µm) · NEW!							
	1 mL	730731P45		730733P45				30
	3 mL		730735P45			730739P45		30
	CHROMAFIX® HR-XCW cartridges (85 µm)							
	Size	S		M		L		
	adsorbent weight Ø	155 mg		240 mg		500 mg		
		731774		731775		731776		50
	CHROMABOND® HR-XCW adsorbent (85 µm)							
							730674	100 g



HR-XAW

- Weak basic secondary and tertiary ammonium anion exchanger
exchange capacity >0.5 meq/g, pKa ~ 6
- Base material spherical PS/DVB copolymer, pH stability 1–14
high purity material with highest reproducibility and
lowest blank values due to an optimized production process
- Spherical particles, size 45 µm and 85 µm (standard); pore size 55–65 Å
very large specific surface 850 m²/g; pore volume 1.2–1.4 cm³/g
RP capacity 350 mg/g (caffeine in water)
- Outstanding recovery rates especially for enrichment of acidic analytes

weak anion exchanger

- Recommended application:
Perfluorinated surfactants
Acidic compounds like sul-
fonates
Active ingredients from
heavily matrix-contaminat-
ed samples like, e.g., urine,
plasma, serum
Strong acids with pKa < 1

surfactants from water

Application in accordance with DIN 38407-42

Column type:
CHROMABOND® HR-XAW, 3 mL, 60 mg
REF 730747

Sample: 500 mL water, spiked with 1 mL standard solution
(20 µg/L of each compound)

Conditioning: 2 mL methanol + 5 % ammonia, then 2 mL metha-
nol, finally 2 mL water

Sample application: slowly aspirate sample through the column

Washing: 2 mL water, then 2 mL acetone – acetonitrile – formic
acid (50:50:1, v/v/v), finally 2 mL methanol

Elution: 2 mL methanol with 5 % ammonia

Further analysis: evaporate to dryness in a stream of nitrogen
under slight heating, and redissolve in a suitable solvent for
HPLC

Recovery rates [%]:

Compound	Recovery
1 Perfluoropropionic acid (PFPrA)	103
2 Perfluoropentanoic acid (PFPeA)	94
3 Perfluorohexanoic acid (PFHxA)	94
4 Perfluorooctanoic acid (PFOA)	95
5 Perfluorooctane sulfonate K salt (PFOS)	81
6 Perfluorododecanoic acid (PFDoDA)	82

MN Appl. No. 305140



impregnated with fluorosurfactants?

Standard protocol for CHROMABOND® HR-XAW

Column type:
CHROMABOND® HR-XAW, 3 mL, 200 mg
REF 730748

Sample pretreatment: individual sample preparation with refer-
ence to analytes and matrix

Conditioning: 5 mL methanol

Equilibration: 5 mL water

Sample application: slowly aspirate sample through the column

Washing 1: 25 mmol/L ammonium acetate





Washing 2/Elution 1: 2 mL methanol (neutral and basic com-
pounds), if necessary, further washing steps

Elution 2: after drying 2 x 2 mL methanol – 1 to 5 % ammonia
(strongly acidic compounds)

Further analyses: if necessary, evaporate and redissolve in a
suitable solvent; HPLC or GC

MN Appl. No. 305200

Ordering information

	Volume	Adsorbent weight						Pack of
	CHROMABOND® HR-XAW polypropylene columns (85 µm)							
		30 mg	60 mg	100 mg	150 mg	200 mg	500 mg	
	1 mL	730728		730729				30
	3 mL		730747			730748	730744	30
	6 mL				730749		730745	30
	CHROMABOND® HR-XAW polypropylene columns (45 µm) · NEW!							
	1 mL	730728P45		730729P45				30
	3 mL		730747P45			730748P45		30
	CHROMAFIX® HR-XAW cartridges (85 µm)							
	Size	S		M		L		
	adsorbent weight Ø	155 mg		240 mg		500 mg		
		731771		731772		731773		50
	CHROMABOND® HR-XAW adsorbent (85 µm)							
							730673	100 g



Polymer-based reversed phases for SPE

Easy polar, bifunctionally modified polystyrene–divinylbenzene copolymer

- Polar modified polystyrene–divinylbenzene copolymer with a weak anion exchanger

Specific surface 650–700 m²/g, particle size 80 µm, pore size 50 Å, pH stability 1–14

The Easy effect:

- Without preconditioning
- Due to bifunctional modification much more hydrophilic than conventional polystyrene–divinylbenzene polymers
- Easily wettable with water

- Recommended application:

polar herbicides and pesticides from water (acidic, neutral, basic)
polar phenols from water
polyaromatic compounds
polychlorinated biphenyls
drug analysis from urine, blood, serum, plasma, pharmaceuticals and active ingredients from tablets, creams

Recovery of pesticides

Private communication: Mr. Kühn, GUB, Waldshut Tiengen, Germany

Column type:
CHROMABOND® Easy, 3 mL, 200 mg
REF 730754

Column conditioning:
1 mL water, 3 mL methanol, 1 mL water

Sample application:
aspirate the sample through the column

Elution:
3 x 1 mL acetone

Further analysis:
HPLC with NUCLEOSIL® 120-5 C₁₈

MN Appl. No. 303220

Recovery rates [%]:

Compound	Recovery	Compound	Recovery
Desisopropylatrazine	90	Metalaxyl	96
2,6-Dichlorobenzamide	93	Isoproturon	94
Desethylatrazine	93	Diuron	94
Hexazinone	69	Metazachlor	97
Terbacil	65	Propazine	95
Simazine	81	Terbutylazine	93
Cyanazine	93	Linuron	96
Desethylterbutylazine	91	Metolachlor	97
Methabenzthiazuron	94	Triallate	61
Chlortoluron	91	Standard	64
Atrazine	92		

Ordering information

Volume	Adsorbent weight						Pack of
	CHROMABOND® Easy polypropylene columns						
		30 mg	60 mg	100 mg	200 mg	500 mg	1 g
	1 mL	730751		730752			30
	3 mL		730753		730754	730759	30
	6 mL				730755	730756	30
	15 mL					730757	730758
	CHROMABOND® Easy polypropylene columns • BIGpacks						
					200 mg		
	3 mL				730754.250		250
	6 mL				730755.250		250
	CHROMABOND® LV-Easy						
					200 mg		
	15 mL				732472		30
	CHROMABOND® MULTI 96 Easy						
		96 x 25 mg		96 x 50 mg		96 x 100 mg	
		738520.025M		738520.050M		738520.100M	1
	CHROMABOND® Easy adsorbent						
						730661	20 g

Glass columns on request.



HR-P

polystyrene–divinylbenzene adsorbent resin

- Highly porous polystyrene–divinylbenzene copolymer
- specific surface 1200 m²/g
- particle size 50–100 µm
- very high binding capacity, up to 30% of adsorbent weight (for comparison: silica adsorbents about 3%)

- Recommended application:
- aromatic compounds
- phenols from water
- nitroaromatics from water
- pesticides from water
- PAHs from oil

Aromatic amines from water samples

Private communication M. Leß, T.C. Schmidt, Department of Chemistry, University Marburg, 1997

Compounds investigated: aromatic amines

Column type:

CHROMABOND® HR-P, 3 mL, 200 mg
REF 730108

Sample pretreatment: adjust to pH 9 using 10 mol/L NaOH

Column conditioning: 2 mL each of methanol, acetonitrile and 10⁻⁵ mol/L sodium hydroxide

Sample application:

aspirate sample through the column with about 10 mL/min

Column washing:







wash with 2 mL dist. water, dry 5 min under vacuum

Elution: 3 x 1 mL methanol – acetonitrile (1:1, v/v)

For recovery rates of numerous aromatic amines please see application 301810 at www.mn-net.com/apps.

MN Appl. No. 301810

Ordering information

	Volume	Adsorbent weight				Pack of
	CHROMABOND® HR-P polypropylene columns					
		100 mg	200 mg	500 mg	1 g	
	1 mL	730280				30
	3 mL		730108	730117		30
	6 mL		730119	730111	730118	30
	CHROMABOND® HR-P polypropylene columns - BIGpack					
		200 mg				
	3 mL	730108.250				250
	CHROMABOND® HR-P glass columns					
			200 mg	500 mg	1 g	
	3 mL		730108G			30
	6 mL			730111G	730118G	30
	CHROMABOND® LV-HR-P					
	15 mL		200 mg	732108		
	CHROMAFIX® HR-P cartridges					
	Size	S	M	L		
	Adsorbent weight Ø	200 mg	330 mg	680 mg		
		731839	731840	731841	50	
	CHROMABOND® MULTI 96 HR-P					
			96 x 100 mg			
			738111.100M			1
	CHROMABOND® HR-P adsorbent					
			730615			20 g



Polymer-based phases for SPE

PS-RP / PS-OH⁻ / PS-H⁺ PS-Mix / PS-Ag⁺ / PS-Ba²⁺

phases for RP and ion chromatography

- Base material high purity polystyrene-divinylbenzene copolymers (PS/DVB), pore size 100 Å, particle size 100 µm
- Very low degree of swelling, thus very well suited for chromatography
- Reliable function over the whole pH range from 0-14
- Different modifications for different applications from elimination of nonpolar compounds up to the removal of specific polar components

- Recommended application:
Removal of interfering compounds
 - Improves chromatographic separation, if the interfering components overlap with the analyte in the chromatogram
 - Improves lifetime of the chromatographic column, since interfering components can irreversibly block the column packing
- Enrichment of the analytes

Properties of the individual modifications:

PS-RP	hydrophobic PS/DVB copolymer	removal of organic interfering components from water
PS-OH ⁻	strong PS/DVB anion exchanger, OH ⁻ form capacity 0.6 meq/g	removal or concentration of anions from water increasing the pH value in acidic samples
PS-H ⁺	strong PS/DVB cation exchanger, H ⁺ form capacity 2.9 meq/g	removal or concentration of cations from water decreasing the pH value of basic samples
PS-Mix	mixture of PS-OH ⁻ and PS-H ⁺	desalting of water
PS-Ag ⁺	strong PS/DVB cation exchanger, Ag ⁺ form	removal of halide ions from water
PS-Ba ²⁺	strong PS/DVB cation exchanger, Ba ²⁺ form	removal of sulfate ions from water

Application 301930/302750: removal of halides from aqueous samples shown for the trace analysis of nitrate besides an excess of chloride or bromide

Compounds investigated: 20 ppm nitrate besides 2500 ppm chloride or 500 ppm bromide, respectively



Column type:
CHROMAFIX® PS-Ag⁺ (M)
0.8 mL, Ø 480 mg, REF 731865
Column conditioning: 1 mL dist. water

Sample application and elution:

apply 4 x 1 mL sample fractions to the cartridge, discard 1st mL, collect 2nd, 3rd and 4th mL separately

Further analysis: HPLC with column 250 x 4 mm NUCLEOSIL® Anion II; eluent 2 mmol/L K H phthalate pH 6, 2 mL/min; detection: indirect UV, 280 nm (see applications 110440 and 110450 at www.mn-net.com/apps)

Ordering information

Phase		Volume / Adsorbent weight				Pack of		
	CHROMABOND® PS polypropylene columns							
		3 mL 200 mg	3 mL 500 mg	6 mL 500 mg	6 mL 900 mg			
	PS–RP	730765	730692	730693			30	
	PS–OH [–]	730396	730344	730378			30	
	PS–H ⁺	730690	730376	730377			30	
	PS–Mix		730394		730310		30	
	CHROMAFIX® PS cartridges							
		Size S	Adsorbent weight Ø	Size M	Adsorbent weight Ø	Size L	Adsorbent weight Ø	
	PS–RP	731877	200 mg	731875	320 mg			50
	PS–OH [–]	731868	200 mg	731860	380 mg	731862	800 mg	50
	PS–H ⁺	731867	230 mg	731861	430 mg	731863	900 mg	50
	PS–Mix	731909	230 mg					50
	PS–Ag ⁺	731866	240 mg	731865	480 mg			50
	PS–Ba ²⁺	731871	280 mg	731870	550 mg			50












C₁₈ ec / C₁₈ ec f (f = fast flow)

- Base material silica, pore size 60 Å, particle size 45 µm for C₁₈ ec, 100 µm for C₁₈ ec f (for fast flow), specific surface 500 m²/g, pH stability 2–8
- Octadecyl phases, endcapped, carbon content 14%
- Very nonpolar, hydrophobic interactions with a wide variety of organic compounds
- Advantageous for clean-up of samples with large structural variations (polarity differences)

octadecyl silica, endcapped

- Recommended application:
nonpolar compounds
aflatoxins, amphetamines, antibiotics, antiepileptics, barbiturates, caffeine, drugs, preservatives, fatty acids, nicotine, PAHs, pesticides, PCBs, heavy metals, vitamins
very well suited for desalting of samples
- C₁₈ ec f for viscous samples

Ordering information

	Volume	Adsorbent weight						Pack of
	CHROMABOND® C ₁₈ ec polypropylene columns							
		100 mg	200 mg	500 mg	1 g	2 g	5 g	10 g
	1 mL	730011						100
	3 mL		730012	730013				50
	6 mL			730014	730015	730141		30
	15 mL					730404		20
	45 mL						730405	20
	70 mL							730259
	CHROMABOND® C ₁₈ ec polypropylene columns · BIGpacks							
			500 mg	1 g				
	3 mL		730013.250					250
	6 mL		730014.250	730015.250				250
	CHROMABOND® C ₁₈ ec glass columns							
		200 mg	500 mg	1 g				
	3 mL	730012G	730013G					50
	6 mL		730014G	730015G				30
	CHROMABOND® LV-C ₁₈ ec							
		200 mg	500 mg					
	15 mL	732012	732013					30
	CHROMAFIX® C ₁₈ ec cartridges							
	Size	S	M	L				
	Adsorbent weight Ø	270 mg	530 mg	950 mg				
		731804	731805	731806				50
	CHROMABOND® MULTI 96 C ₁₈ ec							
		96 x 25 mg	96 x 50 mg	96 x 100 mg				
		738011.025M	738011.050M	738011.100M				1
	CHROMABOND® C ₁₈ ec adsorbent							
							730611	100 g
	CHROMABOND® C ₁₈ ec f polypropylene columns (fast flow)							
		200 mg	500 mg	1 g				
	3 mL	730269	730018					50
	6 mL		730016	730010				30
		CHROMABOND® C ₁₈ ec f adsorbent (fast flow)						
							730613	100 g



Silica-based reversed phases for SPE













C₁₈ / C₁₈ f (f = fast flow)

- Base material silica, pore size 60 Å, particle size 45 µm for C₁₈, 100 µm for C₁₈ f (for fast flow), specific surface 500 m²/g, pH stability 2–8
- Octadecyl phases, not endcapped, carbon content 14%
- Similar to C₁₈ ec, however possesses more free silanols (SiOH), which allow secondary interactions with polar groups of the analytes

octadecyl silica

- Recommended application: nonpolar compounds, pesticides
- C₁₈ f for viscous samples

Ordering information

	Volume	Adsorbent weight						Pack of	
	CHROMABOND® C ₁₈ polypropylene columns								
		100 mg	200 mg	500 mg	1 g	2 g	5 g	10 g	
	1 mL	730001						100	
	3 mL		730002	730003				50	
	6 mL			730004	730005	730130		30	
	15 mL					730028		20	
	45 mL						730400	20	
	70 mL							730261	10
	CHROMABOND® C ₁₈ polypropylene columns · BIGpacks								
			500 mg	1 g					
	3 mL		730003.250					250	
	6 mL		730004.250 730005.250					250	
	CHROMABOND® C ₁₈ glass columns								
			500 mg	1 g					
	3 mL		730003G					50	
	6 mL		730004G 730005G					30	
		CHROMABOND® LV-C ₁₈							
			200 mg						
15 mL			732002					30	
	CHROMAFIX® C ₁₈ cartridges								
	Size	S		M		L			
	Adsorbent weight Ø	240 mg		480 mg		950 mg			
		731801		731802		731803		50	
		CHROMABOND® MULTI 96 C ₁₈							
			96 x 25 mg				96 x 100 mg		
		738001.025M				738001.100M		1	
	CHROMABOND® C ₁₈ adsorbent								
								730602	100 g
		CHROMABOND® C ₁₈ f polypropylene columns (fast flow)							
		200 mg	500 mg	1 g					
3 mL			730402	730008				50	
	6 mL		730403	730009			30		
		CHROMABOND® C ₁₈ f adsorbent (fast flow)							
									730612



C₁₈ Hydra

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Special octadecyl phase for polar analytes, not end-capped, carbon content 15 %

octadecyl silica for polar analytes

- Recommended application: more polar compounds like pesticides and their polar degradation products, phenols, phenoxycarboxylic acids, nitroaromatics, pharmaceuticals

Pesticides from water

Compounds investigated: triazines and carboxylic amides

Column type:
CHROMABOND® C₁₈ Hydra, 6 mL, 2 g
REF 730301

Sample pretreatment: adjust 1000 mL water to pH 7–8 with diluted NH₃ and add 100 µL of the internal standards (1 µg/L).

Column conditioning: 2 x 5 mL methanol, then 2 x 5 mL dist. water

Sample application: force or aspirate the sample through the column. Then dry for 2 h with 2 bar N₂.

Elution: slowly aspirate 10 mL methanol through the column. Evaporate the eluate to dryness in a tapered flask with a rotation evaporator at 30 °C and store in a refrigerator for ~ 15 min. Redissolve the residue in 200 µL cold, fresh *n*-hexane and transfer the solution to a conic HPLC vial (e.g., REF 702891). Store the solution in a refrigerator until chromatography.







Recovery rates: between 95 and 100 %

Further analysis: GC with OPTIMA® δ-3 or OPTIMA® δ-6 (e.g., application 250420) or HPLC in accordance with EN ISO 11369: 1997 on NUCLEOSIL® 120-3 C₁₈ (application 110880)

MN Appl. No. 302060



Ordering information

	Volume	Adsorbent weight						Pack of
	CHROMABOND® C ₁₈ Hydra polypropylene columns							
		50 mg	100 mg	200 mg	500 mg	1 g	2 g	3 g
	1 mL	730294	730295					100
	3 mL			730296	730297	730298		50
	6 mL				730299	730300	730301	730302
	CHROMABOND® C ₁₈ Hydra glass columns							
			200 mg	500 mg	1 g			
	3 mL		730296G	730297G	730298G			50
	6 mL			730299G	730300G			30
	CHROMABOND® LV-C ₁₈ Hydra							
	15 mL		200 mg					30
	CHROMAFIX® C ₁₈ Hydra cartridges							
	Size	S		M		L		
	Adsorbent weight Ø	270 mg		530 mg		950 mg		
		731730		731731		731732		50
	CHROMABOND® MULTI 96 C ₁₈ Hydra							
						96 x 100 mg		
						738294.100M		1
	CHROMABOND® C ₁₈ Hydra adsorbent							
							730628	100 g



Silica-based reversed phases for SPE






C₈

octyl silica

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Octyl phase, not endcapped, carbon content 8%
- Similar to C₁₈, however slightly more polar
- secondary interactions with polar compounds are more pronounced due to shorter alkyl chains

- Recommended application: pesticides, PCB

Ordering information

	Volume	Adsorbent weight				Pack of
	CHROMABOND® C ₈ polypropylene columns					
		100 mg	200 mg	500 mg	1 g	
	1 mL	730021				100
	3 mL		730022	730023		50
	6 mL			730024	730134	30
	CHROMABOND® C ₈ glass columns					
				500 mg		
	6 mL			730024G		30
	CHROMABOND® LV-C ₈					
				500 mg		
	15 mL			732023		30
	CHROMAFIX® C ₈ cartridges					
	Size	M				
	Adsorbent weight Ø			520 mg		
	CHROMABOND® MULTI 96 C ₈					
					96 x 100 mg	
					738021.100M	1
	CHROMABOND® C ₈ adsorbent					
					730601	100 g






C₄

butyl silica

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Butyl phase, not endcapped, carbon content 7%
- Slightly more polar than C₁₈ or C₈, due to shorter alkyl chains the silica surface is not completely shielded

- Recommended application: compounds, which are too strongly retained on C₁₈ or C₈ e.g., analgetics from blood

Ordering information

	Volume	Adsorbent weight		Pack of
	CHROMABOND® C ₄ polypropylene columns			
		100 mg	500 mg	
	1 mL	730225		100
	3 mL	730227		50
	CHROMAFIX® C ₄ cartridges			
	Size	S	M	
	Adsorbent weight Ø	220 mg	440 mg	
		731740	731741	50
	CHROMABOND® C ₄ adsorbent			
			730651	100 g

Glass columns, LV columns and MULTI 96 on request.



C₂

dimethyl silica

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Dimethyl phase, not endcapped, carbon content 4%
- Similar to C₄

- Recommended application: e.g., antiepileptics from plasma

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® C ₂ polypropylene columns				
		100 mg	500 mg	1 g	
	1 mL	730169			100
	3 mL	730221			50
	6 mL	730409	730410	30	
	CHROMABOND® C ₂ adsorbent				
			730652	100 g	

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.



Silica-based reversed phases for SPE

C₆H₁₁ ec

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Cyclohexyl phase, endcapped, carbon content 9%
- Alternative phase for the midpolar range

cyclohexyl silica, endcapped

- Recommended application:
phenols from water
chloroanilines from waste water
anthelmintics from tissue

Comparison of different phases for phenol analysis

Compounds investigated:

phenol, 2,4-dinitrophenol, pentachlorophenol

Column types:

CHROMABOND® C₁₈, 6 mL, 2000 mg, REF 730130

CHROMABOND® C₆H₁₁ ec, 6 mL, 2000 mg, REF 730469

Column conditioning:


10 mL acetone, 10 mL methanol, and 10 mL dist. water (pH 2)

Sample application:

aspirate the sample through the column.

Elution:



10 mL methanol



MN Appl. No. 302150

phenol 2,4-dinitrophenol pentachlorophenol

Ordering information

	Volume	Adsorbent weight		Pack of
	CHROMABOND® C ₆ H ₁₁ ec polypropylene columns			
		500 mg	1 g	
	3 mL	730442		50
	6 mL	730443	730444	30
	CHROMABOND® C ₆ H ₁₁ ec adsorbent			
			730631	100 g

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.

For further applications on CHROMABOND® phases
see our online application database at

www.mn-net.com/apps



C₆H₅

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Phenyl phase, carbon content 8%
- Polarity similar to C₈
- In addition to hydrophobic interactions more selective adsorption is possible by π–π interactions due to the electron density of the phenyl ring.

phenyl silica

- Recommended application:
 - aflatoxins
 - caffeine
 - phenols

Flavor compounds from brandy

Compounds investigated: asarone, quinine, coumarin, quassin

Column type:
CHROMABOND® C₆H₅, 6 mL, 1000 mg
REF 730412

Sample pretreatment:
mix 10 mL sample with 90 mL water and 10 g sodium chloride and adjust to pH 7 with 0.1 mol/L sodium hydroxide solution

Column conditioning:
10 mL methanol, then 10 mL dist. water

Sample application:
slowly force or aspirate the sample through the column

Column washing:
2.5 mL water, then 2.5 mL pentane



- Elution:*
- 2 x 2.5 mL pentane – diethyl ether (7:3, v/v):
asarone, coumarin
 - 10 mL 1 mol/L basic methanol – diethyl ether (9:1, v/v): quinine
 - 5 mL chloroform: quassin

MN Appl. No. 300170



Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® C ₆ H ₅ polypropylene columns				
		100 mg	200 mg	500 mg	
	1 mL	730083			100
	3 mL		730411	730084	50
	CHROMABOND® C ₆ H ₅ adsorbent				
				730606	100 g

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.



Silica-based normal phases for SPE



CN

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Cyanopropyl phase, carbon content 5.5 %
- Polar to midpolar
- In addition to weak hydrophobic interactions selective interactions are possible due to the high electron density of the CN group.

cyanopropyl silica

- Recommended application:
cyclosporins
carbohydrates

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® CN polypropylene columns				
		100 mg	200 mg	500 mg	
	1 mL	730061			100
	3 mL		730420	730063	50
	6 mL			730421	30
	CHROMABOND® CN adsorbent				
				730607	100 g

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.



OH (Diol)

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Diol phase, carbon content 5.5 %
- Polar
- Properties similar to SiOH

diol silica

- Recommended application:
antibiotics
prostaglandins

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® OH (Diol) polypropylene columns				
		100 mg	200 mg	500 mg	
	1 mL	730051			100
	3 mL		730417	730053	50
	6 mL			730418	30
	CHROMABOND® OH (Diol) adsorbent				
				730605	100 g

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.



NH₂

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Aminopropyl phase, carbon content 3.5 %
- Polar, weak anion exchanger

aminopropyl silica

- Recommended application:
trace elements
lipids

Metals: trace elements from water

Compounds investigated: Al, Be, Cu, Cr(VI), Mo(VI), V(V)

Column type:

CHROMABOND® NH₂, 3 mL, 500 mg

REF 730033

Sample pretreatment:

mix 100 mL water sample with 5 mL 0.001 % alizarinsulfonic acid solution and adjust to pH 5.5 with acetic acid or sodium acetate

Column conditioning:

2 column volumes 1 mol/L nitric acid, then 2 column volumes dist. water

Sample application:

force or aspirate sample through the column with 3–4 mL/min

Column washing:

2 mL dist. water; dry column under vacuum for 4 min






Elution:

2 column volumes 2 mol/L nitric acid

MN Appl. No. 301910



Ordering information

	Volume	Adsorbent weight				Pack of
	CHROMABOND® NH ₂ polypropylene columns					
		100 mg	200 mg	500 mg	1 g	
	1 mL	730031				100
	3 mL		730413	730033		50
	6 mL			730180	730626	30
	CHROMABOND® NH ₂ polypropylene columns · BIGpack					
			500 mg			
3 mL		730033.250			250	
	CHROMABOND® NH ₂ glass columns					
				500 mg	1 g	
	3 mL			730033G		50
	6 mL			730180G	730626G	30
	CHROMABOND® LV-NH ₂					
				500 mg		
15 mL		732033			30	
	CHROMAFIX® NH ₂ cartridges					
	Size	S				
	Adsorbent weight Ø	220 mg				
		731813				50
	CHROMABOND® MULTI 96 NH ₂					
				96 x 100 mg		
				738031.100M		1
	CHROMABOND® NH ₂ adsorbent					
				730603		100 g



Silica-based normal phases for SPE

SiOH

- Unmodified, weakly acidic silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8

Very polar






Adsorbs humidity from air, for this reason it should be kept well closed and if necessary dried before use

Due to its high affinity for polar compounds it should not be conditioned with polar (e.g., methanol) or water-containing solvents.

unmodified silica

- Recommended application:
aflatoxins
chloramphenicol
pesticides
steroids
vitamins

Ordering information

Volume	Adsorbent weight							Pack of	
	CHROMABOND® SiOH polypropylene columns								
	100 mg	200 mg	500 mg	1 g	2 g	5 g	10 g	50 g	
	1 mL	730071						100	
	3 mL		730214	730073				50	
	6 mL			730070	730075	730107		30	
	15 mL					730217		20	
	45 mL						730406	20	
	70 mL							730072	10
150 mL								730473	10
CHROMABOND® SiOH polypropylene columns • BIGpacks									
		500 mg	1 g	2 g					
3 mL		730073.250						250	
6 mL			730075.250	730107.250				250	
CHROMABOND® SiOH glass columns									
		200 mg	500 mg	1 g	2 g				
3 mL		730214G	730073G					50	
6 mL			730070G	730075G	730107G			30	
	CHROMABOND® LV-SiOH								
		200 mg	500 mg						
	15 mL		732072	732073				30	
	CHROMAFIX® SiOH cartridges								
	Size	S		M		L			
	Adsorb. weight Ø	230 mg		420 mg		880 mg			
		731828		731829		731830		50	
	CHROMABOND® MULTI 96 SiOH								
						96 x 100 mg			
						738071.100M		1	
	CHROMABOND® SiOH adsorbent								
						730608		100 g	



Alox A / Alox N / Alox B






aluminium oxide, acidic, neutral, basic

Aluminium oxide, high purity, pore volume 0.90 mL/g, particle size 60–150 µm, specific surface 150 m²/g

Properties of the individual modifications:

Alox A:	aluminium oxide, acidic	pH value 4 ± 0.5
Alox N:	aluminium oxide, neutral	pH value 7 ± 0.5
Alox B:	aluminium oxide, basic	pH value 9.5 ± 0.5

Ordering information

	Phase	Volume	Adsorbent weight			Pack of
	CHROMABOND® Alox polypropylene columns					
			500 mg	1 g	4 g	
	Alox A	3 mL	730452			50
	Alox A	6 mL	730453	730017		30
	Alox A	45 mL			730455	20
	Alox N	3 mL	730446			50
	Alox N	6 mL	730447	730139		30
	Alox N	45 mL			730250	20
	Alox B	3 mL	730429			50
	Alox B	6 mL	730466	730020		30
	Alox B	45 mL			730467	20
	CHROMABOND® Alox glass columns					
				1 g		
	Alox N	6 mL		730139G		30
	Alox B	6 mL		730020G		30
	CHROMABOND® LV-Alox					
				1 g		
	Alox A	15 mL		732210		30
	Alox N	15 mL		732091		30
	Alox B	15 mL		732205		30
	CHROMAFIX® Alox cartridges					
		Size	M	L		
		Adsorb. weight Ø	850 mg	1700 mg		
	Alox N		731844	731845		50
	CHROMABOND® MULTI 96 Alox					
				96 x 100 mg		
	Alox A			738253.100M		1
	Alox N			738251.100M		1
	Alox B			738252.100M		1
	CHROMABOND® Alox adsorbents					
	Alox A				730642	100 g
	Alox N				730641	100 g
	Alox B				730640	100 g



Normal phases for SPE




Florisol®

Matrix magnesium silicate (MgO – SiOH 15:85), high purity, particle size 150–250 µm

magnesium silicate

Recommended application:
organic tin compounds,
aliphatic carboxylic acids,
PCBs, PAHs

Ordering information

	Volume	Adsorbent weight				Pack of
	CHROMABOND® Florisil® polypropylene columns					
		200 mg	500 mg	1 g	2 g	
	3 mL	730457	730081			50
	6 mL		730238	730082	730239	30
	CHROMABOND® Florisil® polypropylene columns · BIGpack					
				1 g		
	6 mL			730082.250		250
	CHROMABOND® Florisil® glass columns					
			500 mg	1 g	2 g	
	6 mL		730238G	730082G	730239G	30
	CHROMAFIX® Florisil® cartridges					
	Size				L	
	Adsorbent weight Ø				990 mg	
					731848	50
	CHROMABOND® Florisil® adsorbent					
					730622	100 g

LV columns and MULTI 96 on request




PA

Matrix polyamide 6, unmodified, high purity, particle size 40–80 µm

polyamide 6

Recommended application:
flavonoids, PAHs

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® PA polypropylene columns				
		200 mg	500 mg	1 g	
	3 mL	730384	730126		50
	6 mL		730007	730127	30
	CHROMAFIX® PA cartridges				
	Size	S		L	
	Adsorbent weight Ø	170 mg		620 mg	
		731849		731851	50
	CHROMABOND® PA adsorbent				
				730660	100 g

Glass columns, LV columns and MULTI 96 on request






PCA

propylcarboxylic acid cation exchanger based on silica

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Propylcarboxylic acid modified silica
- Weakly acidic cation exchanger (WCX)

- Recommended application: strong cations

Ordering information

	Volume	Adsorbent weight		Pack of
	CHROMABOND® PCA polypropylene columns			
		500 mg	1 g	
	3 mL	730482		50
	6 mL	730483	730484	30
	CHROMABOND® LV-PCA			
		500 mg		
	15 mL	732482		30
	CHROMABOND® PCA adsorbent			
			730629	100 g

Glass columns, CHROMAFIX® cartridges and MULTI 96 on request.



PSA

propylsulfonic acid cation exchanger based on silica

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Propylsulfonic acid modified silica
- Very strong cation exchanger (capacity ~ 0.7 meq/g)
- Contrary to the SA phase no π–π interactions

- Recommended application: weak cations

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® PSA polypropylene columns				
		100 mg	500 mg	1 g	
	1 mL	730460			100
	3 mL	730462			50
	6 mL	730464			30
	CHROMABOND® PSA adsorbent				
				730630	100 g

Glass columns, LV columns, CHROMAFIX® cartridges and MULTI 96 on request.



Silica-based ion exchangers for SPE

SA

benzenesulfonic acid cation exchanger based on silica (SCX)

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Benzenesulfonic acid modified silica
- strongly acidic cation exchanger (capacity ~ 0.5 meq/g)
- adsorbent with hydrophobic and π–π interactions (benzene ring)
- Ion exchange of organic compounds from aqueous matrix
- elution of interesting compounds with solvent systems, which compensate the ionic and nonpolar interactions, e.g., methanolic HCl

- Recommended application:
- amino acids
- amines
- chlorophyll
- PCB

Sulfonamides in meat and kidney

B. Pacciarelli et al., Mitt. Gebiete Lebensm. Hyg. 82 (1991) 45–55
Compounds investigated: sulfaguanidine, sulfanilamide, sulfadiazine, sulfathiazole, sulfapyridine, sulfamerazine, sulfamethizole, sulfadimidine, sulfamethoxypyridazine, sulfachlorpyridazine, sulfadoxine, sulfadimethoxine

Column type:
 CHROMABOND® SA (= SCX), 3 mL, 500 mg
 REF 730077

Sample pretreatment: homogenize 10 g sample and 60 mL dichloromethane – acetone (1:1, v/v) for 30 s with a Polytron. Centrifuge the homogenate for 10 min at 2500 rpm. Filter the organic phase and wash the filter residue with a little dichloromethane – acetone. Add 5 mL glacial acetic acid to the filtered extract.

Column conditioning: apply 6 mL hexane and suck air until the column is dry (10 min). Then apply 6 mL dichloromethane – acetone – glacial acetic acid (10:10:1, v/v/v). Now the column must not run dry.

Sample application: 1/10 of the extract volume, flow rate about 2 mL/min; the column must not run dry

Column washing: 5 mL water, then 5 mL methanol; dry for 10 min under vacuum. Now suck NH₃ gas through the column until the acid is neutralized. To control the neutralization process, press air through the column: a wet pH paper should indicate a neutral or basic pH value.







Elution: 3 mL methanol (1–2 mL/min); carefully concentrate the eluate on a rotation evaporator (40 °C/100 mbar), dissolve the residue in 0.5 mL of 5.5 % acetonitrile in buffer (1.641 g sodium acetate in 1 L water, adjusted to pH 5 with glacial acetic acid) and centrifuge.

Further analysis: HPLC

MN Appl. No. 302710

Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight				Pack of
	CHROMABOND® SA polypropylene columns					
		100 mg	200 mg	500 mg	1 g	
	1 mL	730076				100
	3 mL		730275	730077		50
	6 mL			730425	730212	30
	CHROMABOND® SA polypropylene columns · BIGpack					
				500 mg		
	3 mL			730077.250		250
	CHROMABOND® LV-SA					
				500 mg		
	15 mL			732083		30
	CHROMAFIX® SA cartridges					
	Size	S	M	L		
	Adsorbent weight Ø	220 mg	450 mg	920 mg		
		731831	731832	731833		50
	CHROMABOND® MULTI 96 SA					
				96 x 100 mg		
				738141.100M		1
	CHROMABOND® SA adsorbent					
				730609		100 g

Glass columns on request



SB quaternary ammonium anion exchanger based on silica (SAX)

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Silica modified with quaternary amine
- strongly basic anion exchanger (capacity ~ 0.3 meq/g)
- Not suited for very strong anions such as sulfonic acids, because these are difficult to elute

- Recommended application:
 - organic acids
 - caffeine
 - saccharin

Vitamins: folic acid from food (e.g., wheat germs)

Column type:
CHROMABOND® SB (= SAX), 3 mL, 500 mg
REF 730079

Sample pretreatment:
homogenize 10 g food sample in 100 mL 0.01 mol/L phosphate buffer pH 7.4 and filter

Column conditioning: 2 column volumes *n*-hexane, then 2 column volumes methanol, finally 2 column volumes dist. water

Sample application: force or aspirate 10 mL of the filtrate through the column

Column washing: 2 column volumes dist. water

Elution: 5 mL 10 % sodium chloride in 0.1 mol/L sodium acetate buffer

MN Appl. No. 300650



Ordering information

Volume	Adsorbent weight				Pack of
	CHROMABOND® SB polypropylene columns				
		100 mg	200 mg	500 mg	1 g
	1 mL	730078			100
	3 mL		730322	730079	50
	6 mL			730426	730323
	CHROMABOND® SB polypropylene columns · BIGpack				
	3 mL			500 mg	250
	CHROMABOND® LV-SB				
	15 mL			500 mg	30
	CHROMAFIX® SB cartridges				
	Size	S	M	L	
	Adsorbent weight Ø	230 mg	460 mg	920 mg	
		731834	731835	731836	50
	CHROMABOND® MULTI 96 SB				
				96 x 100 mg	
	CHROMABOND® SB adsorbent				
				730610	100 g

Glass columns on request



Special SPE phases · pharmaceutical analyses

Drug

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Special bifunctional modification – C₈ / SA (strong cation exchanger – benzenesulfonic acid)

special silica phase for drug analysis

- Recommended application: enrichment of acidic, neutral and basic drugs from urine or plasma

Drugs from blood serum

W. Weinmann, M. Renz, C. Pelz, P. Brauchle, S. Vogt, S. Pollak, Blutalkohol 35 (1998), 1–9

Compounds investigated:

benzoylecgonine, amphetamine, codeine, morphine

Column type:

CHROMABOND® Drug, 3 mL, 200 mg
REF 730168

Sample pretreatment:

0.1 mL blood serum are mixed with 1.4 mL of a 0.1 mol/L KH₂PO₄ buffer (pH 6) and centrifuged

Column conditioning:

2 mL methanol, then 2 mL 0.1 mol/L KH₂PO₄ buffer (pH 6)

Sample application:

slowly force or aspirate the supernatant from the sample pretreatment through the column

Column washing:

2 mL 0.1 mol/L KH₂PO₄ buffer (pH 6), then 1 mL 0.1 mol/L acetic acid, then 2 mL methanol; finally dry the column first by centrifugation (2 min, 4000 U/min), then under vacuum for 10 min

Elution:

1.5 mL dichloromethane – 2-propanol – 25% ammonia solution (80:20:2, v/v/v)

Further analysis: HPLC with NUCLEOSIL® 100-5 C₁₈ AB (application 110240) or GC/MS after derivatization with perfluoropropanoic acid anhydride/pentafluoropropanol, e.g., with column OPTIMA® 5 MS, 0.25 mm film, 30 m x 0.25 mm ID, (REF 726220.30)





MN Appl. No. 302020



Poppy seeds as source of opiates

Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight			Pack of
	CHROMABOND® Drug polypropylene columns				
		100 mg	200 mg	500 mg	
	1 mL	730681			100
	3 mL	730168			50
	6 mL	730682			30
	CHROMABOND® Drug polypropylene columns · BIGpack				
		200 mg			
	3 mL	730168.250			250
	CHROMABOND® LV-Drug				
		200 mg			
	15 mL	732168			30
	CHROMABOND® MULTI 96 Drug				
		96 x 100 mg			
		738161.100M			1



Drug II

extraction of THC and derivatives, acidic analytes from biological fluids (urine, blood, etc.)

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Special bifunctional modification – C₈ / SB (strong anion exchanger – quaternary amine –NR₃⁺)
- Two primary retention mechanisms facilitate use of very strong interferant-eluting solvents, resulting in very pure extracts

- Recommended application: extraction of THC and derivatives from urine, blood, serum, plasma
- acidic analytes from biological fluids

11-nor-Δ⁹-THC-carboxylic acid from urine

Compounds investigated:

tetrahydrocannabinol, 11-nor-Δ⁹-THC-carboxylic acid

Column type:

CHROMABOND® Drug II, 3 mL, 200 mg
REF 730680

Sample pretreatment: add 300 µL 10 mol/L potassium hydroxide solution and internal standard (for GC/MS deuterium labeled 11-nor-9-THC-carboxylic acid) to 5 mL urine. Vortex the sample and then hydrolyze at 60 °C for 15 min. Cool sample and add 200 µL glacial acetic acid and 2 mL 50 mmol/L ammonium acetate solution. If necessary, adjust sample pH to 6–7.

Column conditioning: 2 mL methanol, 2 mL dist. water; equilibrate column with 2 mL 50 mmol/L ammonium acetate buffer
Sample application: slowly force or aspirate the sample through the column (1–2 mL/min)

Column washing: elute interferants with 10 mL methanol – water (1:1, v/v); dry the column for 10 min at high vacuum; further wash the column with 2 mL acetonitrile and dry for another 2 min
Elution: elute THC metabolites with 3 mL hexane – ethyl acetate – glacial acetic acid (75:25:1, v/v/v)

Further analysis: we recommend GC/MS on an OPTIMA® 5 MS column after derivatization with 50 µL Silyl-991 (REF 701480; BSTFA – TMCS 99:1) at 70 °C for 20 min; inject 1–2 µL onto the GC column.

Recovery rates: 70–80 %

MN Appl. No. 303880



Ordering information

Volume	Adsorbent weight			Pack of
	CHROMABOND® Drug II polypropylene columns			
		100 mg	200 mg	500 mg
	1 mL	730685		100
	3 mL		730680	50
	6 mL			30
	CHROMABOND® LV-Drug II			
	15 mL		200 mg	30
	CHROMABOND® MULTI 96 Drug II			
			96 x 100 mg	1
			738680.100M	



SPE phases for pharmaceutical applications

Crosslinks

special phase for enrichment of collagen crosslinks

- Special cellulose phase for enrichment of collagen crosslinks

- Recommended application: collagen crosslinks in urine

Pyridinoline and deoxypyridinoline are collagen crosslinks occurring in bones and cartilage. If these substances are released, they can be detected in the urine. In cases of increased bone catabolism (e.g., during osteoporosis) the urine concentrations of pyridinoline and deoxypyridinoline are increased.

Pyridinium crosslinks from urine

Compounds investigated: pyridinoline, deoxypyridinoline

Column type:

CHROMABOND® Crosslinks, 3 mL, 300 mg
REF 730458

Sample pretreatment: 250 µL urine and 50 µL of an internal standard (e.g., pyridoxine) are hydrolyzed in 250 µL conc. HCl at about 100–105 °C for 12–16 h. Then 2.5 mL wash solution (*n*-butanol – glacial acetic acid 80:20, v/v) are added to the hydrolyzate.

Column conditioning:

5 mL of the wash solution

Sample application:

force or aspirate the pre-treated sample through the column. Discard the flow-through. Wash with 15–25 mL of the wash solution.

Elution:

force or aspirate 3–5 mL dist. water through the column

MN Appl. No. 302070

Ordering information

Volume	Adsorbent weight	Pack of
CHROMABOND® Crosslinks polypropylene columns		
	300 mg	
3 mL	730458	50
Product for research purposes only (see page 325)		

Tetracycline

special phase for enrichment of tetracyclines

- Silica phase with special C₁₈ modification, tested for tetracyclines

Constant recovery rates for the title compounds (every batch individually tested)

- Recommended application: tetracyclines from biological samples

Tetracyclines from musculature

Private communication of Mr. Lippold, Chemisches Landesuntersuchungsamt (Chem. Research Agency) Freiburg, Germany

Compounds investigated:

tetracycline, oxytetracycline, chlorotetracycline (100–500 mg/kg)

Column type:

CHROMABOND® Tetracycline, 6 mL, 500 mg
REF 730315

Sample pretreatment:

see detailed description in appl. 302030 at

www.mn-net.com/apps

Column conditioning:

1 column volume methanol, 1 column volume dist. water, then 1 column volume EDTA – succinate buffer

CAUTION: DO NOT LET THE COLUMN RUN DRY!

Sample application:

force or aspirate 50 mL of the eluate from the sample pretreatment through the CHROMABOND® column

Column washing:

2 mL dist. water (removal of Cu ions), 2 mL *n*-hexane

Elution: 7.5 mL methanol into a 25-mL tapered flask. Add 1 mL of an ethylene glycol – methanol mixture (22 g ethylene glycol filled up to 100 mL with methanol) and evaporate to dryness with a rotation evaporator (max. 40 °C). Fill up the residue to 400 mL with 0.1 mol/L McIlvaine-EDTA buffer (52.5 g citric acid · H₂O, 44.5 g Na₂HPO₄ · H₂O and 93 g Titriplex III dissolved in 2.5 L dist. water, adjusted to pH 4 with NaOH).

Further analysis:


HPLC with column 250 x 4 mm NUCLEOSIL® 100-5 C₁₈ HD (application 110710)

Recovery rates: tetracycline, chlorotetracycline ~ 50–70 %, oxytetracycline ~ 60–80 %

MN Appl. No. 302030



Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® Tetracycline polypropylene columns		
		500 mg	
	6 mL	730315	30
	Product for research purposes only (see page 325)		

HR-P-AOX

AOX from waters with high salt loads (DIN 38409 – H22)

Special PS/DVB phase

Recommended application:
extraction of AOX (adsorbable organically bonded halogens) from waters containing high salt loads or organic pollutants in accordance with DIN 38409 – H22

AOX from water (DIN 38409 – H22)

Column type:
CHROMABOND® HP-P-AOX, 6 mL, 500 mg
REF 730111.AOX

Column conditioning:
5 mL methanol, 10 mL dist. water.
Do not let the column run dry!

Sample application:
force or aspirate 100 mL original or diluted sample (pH 1) through the column (3–5 mL/min).
Do not let the column run dry!


Column washing:
50 mL nitrate rinsing solution (dissolve 17 g NaNO₃ in 100 mL dist. water, add 1.4 mL HNO₃ 10 mol/L, fill up to 1000 mL; take 50 mL and fill to 1000 mL with dist. water). Discard the flow-through.

Elution:
slowly aspirate 1 x 1 mL, then 1 x 4 mL methanol and 10 mL dist. water through the column.
Collect eluates in 100 mL volumetric flask and fill to 100 mL with dist. water.

MN Appl. No. 302080



Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® HR-P-AOX polypropylene columns		
		200 mg	500 mg
	6 mL	730119.AOX	730111.AOX
			30



SPE phases for environmental analysis

C₁₈ PAH

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Special octadecyl modification for enrichment of PAH, not endcapped, carbon content 14%

octadecyl silica for PAH analysis

- Recommended application: PAHs from water

PAHs from water

Column type:
CHROMABOND® C₁₈ PAH, 6 mL, 2 g
REF 730166

Sample pretreatment:
mix 1000 mL water sample with 10 mL methanol

Column conditioning:
1 column volume methanol, then 1 column volume dist. water

Sample application: aspirate 1000 mL water sample through the column (~15–20 mL/min), then dry column (stream of nitrogen or 24 h in a desiccator over P₂O₅)

Elution: elute with 4 mL acetonitrile – toluene (3:1, v/v) and then evaporate or fill up to the volume required

Recovery rates (50 ng/L per component): naphthalene 87 %, acenaphthylene 89 %, acenaphthene 90 %, fluorene 82 %, phenanthrene 85 %, anthracene 90 %, fluoranthene 89 %, pyrene 89 %, benz[a]anthracene 87 %, chrysene 95 %, benzo[b]-fluoranthene 91 %, benzo[k]fluoranthene 89 %, benzo[a]pyrene 90 %, dibenz[ah]anthracene 97 %, benzo[ghi]perylene 91 %, indeno[1,2,3-cd]pyrene 96 %

MN Appl. No. 301250

Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® C₁₈ PAH polypropylene columns		
	6 mL	2 g 730166	30
	CHROMABOND® C₁₈ PAH glass columns		
	6 mL	730166G	30
	CHROMABOND® C₁₈ PAH adsorbent		
		730616	100 g

NH₂/C₁₈

combination phase for PAH analysis

- Special combination phase: aminopropyl phase for removal of interfering humic acids
- octadecyl phase for enrichment of PAH

- Recommended application: PAHs from water containing humic acids

PAHs from water containing humic acids

Column type:
CHROMABOND® NH₂/C₁₈, 6 mL, 500 mg/1 g glass column
REF 730620G

Sample pretreatment:
mix 500 mL water sample with 25 mL 2-propanol

Column conditioning: 10 mL dichloromethane, 10 mL methanol, then 10 mL dist. water – 2-propanol (9:1, v/v)

Sample application: aspirate 500 mL prepared water sample through the column (~5 mL/min)


Column washing: 2 mL dist. water – 2-propanol (9:1, v/v), then dry column (about 20 min, vacuum)

Elution: 4 x 0.5 mL CH₂Cl₂ (let percolate first 0.5 mL into the column packing without vacuum, then apply light vacuum), if necessary evaporate in a stream of N₂ and fill up with a suitable solvent

MN Appl. No. 301260



Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® NH₂/C₁₈ polypropylene columns		
		500/500 mg	500 mg/1 g
	6 mL	730618	730620 30
	CHROMABOND® NH₂/C₁₈ glass columns		
	6 mL	730618G	730620G 30

Na₂SO₄/Florisol®

hydrocarbons from water in accordance with DIN H-53 / ISO DIS 9377-4

Special combination phase of sodium sulfate and Florisol®

Recommended application: hydrocarbons from drinking, surface and waste waters

Hydrocarbons from water

Column type:
CHROMABOND® Na₂SO₄/Florisol®, 6 mL, 2 g/2 g, glass column, REF 730249G

Internal standard solution:
dissolve 20 mg *n*-tetracontane (C₄₀H₈₂) in petroleum ether, add 20 mL *n*-decane (C₁₀H₂₂) and fill up to one liter with petroleum ether. For preparation of the extraction solution dilute standard solution 1:10 with petroleum ether.

Sample pretreatment:
adjust 900 mL water (10 °C) with HCl (12 mol/L) to pH 2 and add 80 g MgSO₄. Add 50 mL of the extraction solution, close the bottle and stir the suspension intensely for 30 min. Add enough dist. water to separate the organic from the aqueous phase.

Column conditioning: 5 mL petroleum ether

Sample application:
slowly aspirate or force the sample through the column

Elution:
wash with 10 mL petroleum ether. Evaporate the combined solution from sample application and elution to 1 mL at about 75 °C. If necessary, fill up to 1 mL again. (If the hydrocarbon content is high, evaporation to 1 mL may not be necessary.)


Recovery rates: must be > 80 % for *n*-tetracontane.

MN Appl. No. 302090



Solid Phase Extraction

Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® Na₂SO₄/Florisol® polypropylene columns		
		2 g/2 g	
	6 mL	730249	30
	CHROMABOND® Na₂SO₄/Florisol® glass columns		
		2 g/2 g	
	6 mL	730249G	30
	CHROMABOND® Na₂SO₄/Florisol® glass columns · BIGpack		
		2 g/2 g	
	6 mL	730249G.250	250



SPE phases for environmental analysis

CN/SiOH

- Special combination phase
cyanopropyl phase for selective adsorption of polycyclic aromatics via π - π interactions
unmodified silica phase for removal of polar compounds

combination phase for PAH analysis

- Recommended application:
extraction of the 16 PAHs according to EPA from soil samples

PAHs from soil

Column type:
CHROMABOND® CN/SiOH, 6 mL, 500/1000 mg
REF 730135

Sample pretreatment:
dry 30 g soil with sodium sulfate and reflux 4 h with 250 mL petroleum ether in a Soxhlet extractor. For low PAH contents (colorless or weakly colored extracts) concentrate extract to 1/10 of its volume in a rotation evaporator.

Column conditioning:
4 mL petroleum ether

Sample application:
aspirate 20 mL of the extract through the column

Column washing: 2 mL petroleum ether

Elution:
2 x 2 mL acetonitrile – toluene (3:1, v/v), then evaporate or fill to the volume required

Further analysis: HPLC, e.g., with column 100 x 4 mm NUCLEODUR® C₁₈ PAH, 3 μ m, REF 760783.40 according to application 123820 (see page 168)

For recovery rates see application 301310 at www.mn-net.com
MN Appl. No. 301310



Ordering information

Volume	Adsorbent weight	Pack of
CHROMABOND® CN/SiOH polypropylene columns		
	500 mg/1 g	
3 mL	730112	50
6 mL	730135	30
CHROMABOND® CN/SiOH polypropylene columns · BIGpack		
	500 mg/1 g	
6 mL	730135.250	250
CHROMABOND® CN/SiOH glass columns		
	500 mg/1 g	
6 mL	730135G	30

NAN

- Special combination phase:
N: sodium sulfate for removal of trace water;
A: SiOH/AgNO₃ phase for removal of sulfur, sulfur-containing and polar compounds

special phase for PCB analysis

- Recommended application:
extraction of PCB from sludge



PCB from sludge

Compounds investigated: polychlorinated biphenyls (PCB)
This method can also be used for soil samples.

Column type:
CHROMABOND® NAN, 6 mL, 700/2000/700 mg
REF 730149

Sample pretreatment: extract 2 g lyophilized sludge with 70 mL *n*-hexane, evaporate extract and fill to 10 mL with *n*-hexane

Column conditioning: 10 mL *n*-hexane

Sample application: aspirate 2 mL extract into the column

Elution: slowly aspirate 40 mL *n*-hexane through the column with light vacuum, then evaporate and fill to 5 mL with *n*-hexane

Recovery rates:

PCB-28 104 %, PCB-52 100 %, PCB-101 99 %, PCB-138 98 %, PCB-153 101 %, PCB-180 98 %, PCB-209 104 %

MN Appl. No. 301400

Ordering information

Volume	Adsorbent weight	Pack of
CHROMABOND® NAN polypropylene columns		
	400/1400/400 mg	700/2000/700 mg
3 mL	730109	50
6 mL	730149	30
CHROMABOND® NAN polypropylene columns · BIGpack		
	700/2000/700 mg	
6 mL	730149.250	250
CHROMABOND® NAN glass columns		
	700/2000/700 mg	
6 mL	730149G	30
CHROMABOND® NAN adsorbent		
	730619*	100 g

* This product contains harmful substances which must be specially labeled as hazardous. For detailed information please see MSDS.

SA/SiOH

combination phase for PCB analysis

Special combination phase:

SA: strongly acidic cation exchanger based on silica with benzenesulfonic acid modification

SiOH: unmodified silica for removal of polar compounds

Recommended application:

extraction of PCBs from waste oil (hexane extract)

PCB from waste oil

Column type:
CHROMABOND® SA/SiOH, 3 mL, 500/500 mg
REF 730132

Column conditioning: 1 mL *n*-hexane

Sample application: apply 250 µL waste oil sample to the column and aspirate or force it into the adsorbent with 2 x 1 mL *n*-hexane

MN Appl. No. 301390

Elution: aspirate or force another 2 x 500 µL *n*-hexane through the column; collect all *n*-hexane fractions and if necessary adjust concentration for subsequent analysis by either evaporating *n*-hexane in a stream of nitrogen or by dilution with *n*-hexane

Recovery rates:

PCB-28 97 %, PCB-52 96 %, PCB-101 95 %, PCB-138 90 %, PCB-153 95 %, PCB-180 96 %, PCB-209 100 %

Ordering information

Volume	Adsorbent weight	Pack of
CHROMABOND® SA/SiOH polypropylene columns		
	500/500 mg	
3 mL	730132	50
6 mL	730235	50
CHROMABOND® SA/SiOH polypropylene columns · BIGpack		
	500/500 mg	
3 mL	730132.250	250



SPE phases for environmental analysis

SiOH-H₂SO₄/SA

Special combination phase

SiOH-H₂SO₄: H₂SO₄-impregnated silica phase for oxidation of accompanying compounds to ionic and/or polar compounds

SA: strongly acidic cation exchanger based on silica with benzenesulfonic acid modification for removal of ionic and sulfur-containing compounds

This combination column is used together with a SiOH column. Both columns together are available as Kombi-Kit PCB.

combination phase for PCB analysis

Recommended application:

extraction of PCBs from oil with reference to German industrial standard DIN 51527, part 1

PCB in oil samples

determination with reference to German industrial standard DIN 51527

Column type:

CHROMABOND® SiOH-H₂SO₄/SA, 3 mL, 500/500 mg and CHROMABOND® SiOH, 3 mL, 500 mg
REF 730085 and 730073
or Kombi-Kit PCB, REF 730125

Sample pretreatment:

extract oil-contaminated solids with *n*-hexane. Homogenize other oil samples and dissolve 1.5 to 2.0 g in 50 mL *n*-hexane. Water which may cause turbidity can be removed with sodium sulfate.

Column conditioning:

let 1 mL *n*-hexane flow through the CHROMABOND® SiOH-H₂SO₄/SA column

Sample application:

aspirate or force 500 µL sample through the CHROMABOND® SiOH-H₂SO₄/SA column. This phase offers better removal of interfering substances due to sulfonation. Place CHROMABOND® SiOH-H₂SO₄/SA column on top of the SiOH column with the aid of an adapter and after at least 30 s flush sample into the SiOH column with 2 x 1 mL *n*-hexane.

Elution:

elute SiOH column with 3 x 0.5 mL *n*-hexane; adjust to a suitable concentration for subsequent GC analysis by evaporation of *n*-hexane in a stream of nitrogen or by dilution with *n*-hexane

Recovery rates:

PCB-28 99 %, PCB-52 95 %, PCB-101 99 %, PCB-138 94 %, PCB-153 99 %, PCB-180 96 %, PCB-209 101 %

MN Appl. No. 301380



Solid Phase Extraction

Ordering information

Volume	Adsorbent weight	Pack of
CHROMABOND® SiOH-H₂SO₄/SA polypropylene columns		
3 mL	500/500 mg 730085	50
CHROMABOND® SiOH-H₂SO₄/SA polypropylene columns · BIGpack		
3 mL	500/500 mg 730085.250	250
CHROMABOND® SiOH-H₂SO₄/SA glass columns		
3 mL	500/500 mg 730085G	50
Kombi-Kit for extraction of PCB from oil with reference to DIN 51527, part 1		
	25 columns each of CHROMABOND® SiOH-H ₂ SO ₄ /SA and CHROMABOND® SiOH	730125 1 kit




Dry (Na₂SO₄)

special phase for drying of organic samples

- Orange diamond icon: Anhydrous high-purity sodium sulfate which forms Glauber's salt with traces of water
For removal of larger quantities of water several cartridges can be combined in series.

- Orange diamond icon: Recommended application:
removal of traces of water from organic solutions

Ordering information

Adsorbent weight					Pack of
	CHROMAFIX® Dry cartridges				
	Size	S	M	L	
	Adsorbent weight Ø	780 mg	1500 mg	2800 mg	
		731852	731853	731854	50

ABC18

special phase for analysis of acrylamide in food


- Orange diamond icon: Octadecyl silica phase with ion exchange functions for acrylamide analysis
- Orange diamond icon: Recommended application:
clean-up of acrylamide from ultra-heated starch-containing food, such as potato chips and other snacks, french fries, crispbread, cereals etc.

Important notes:

- Orange diamond icon: For "Determination of Acrylamide in Foods, SPE Clean-up Procedure for LC-MS-MS" please see application 303580 at www.mn-net.com/apps.
- Orange diamond icon: Acrylamide is created at temperatures above 100 °C from sugar and proteins, e.g., from potatoes or grain during the process of frying, baking, roasting or grilling. The formation depends on temperature, starting at 120 °C and increasing with more elevated temperatures. In cooked food, no acrylamide is found.
- Orange diamond icon: Minimum concentration of acrylamide should be 70 µg/kg.
- Orange diamond icon: The procedure includes no concentration step.
- Orange diamond icon: Acrylamide and the isotopically labeled form, is carcinogenic, mutagenic and neurotoxic.



Ordering information

	Volume	Adsorbent weight	Pack of
	CHROMABOND® ABC18 polypropylene columns		
		500 mg	
	6 mL	730533	30



Diamino

- Base material silica, pore size 60 Å, particle size 45 µm, specific surface 500 m²/g, pH stability 2–8
- Primary and Secondary Amine functions (PSA), 5 % C
- Removes polar compounds (e.g., organic acids, pigments, sugars) from matrices like fruit or vegetables
- Similar phases: Supelclean™ PSA, Bond Elut® PSA

special silica phase for determination of pesticides in food samples

- Recommended application:
special SPE phase for quick and cheap determination of pesticides in strongly matrix-contaminated samples by GC or HPLC (**QuEChERS** method = **Q**uick **E**asy **C**heap **E**ffective **R**ugged **S**afe)



QuEChERS method and pre-mixes

Within a few years after its development by Anastassiades et al. the QuEChERS method has gained a leading position for determination of pesticide residues in food samples by GC–MS or LC–MS, allowing rapid and cheap clean-up of strongly matrix-contaminated samples.

Standard clean-up of food samples

10 g sample are homogenized with 10 mL acetonitrile. After adding the internal standard the sample is shaken with 4 g MgSO₄ and 1 g NaCl and afterwards centrifuged.
1 mL of the supernatant is spiked with 25 mg CHROMABOND® Diamino and 150 mg MgSO₄ and shaken again. After centrifugation the supernatant is injected into GC/MS.

MN Appl. No. 303770

For optimizing the extraction of pH-dependent compounds, for minimizing decomposition of sensitive substances, and for broadening the matrix spectrum, different modifications of the QuEChERS method have been elaborated.

In addition to the required adsorbent CHROMABOND® Diamino MACHEREY–NAGEL offers a number of individually weighed and premixed extraction and buffer mixtures, specially composed for different sample matrices.

For extraction, the European standard EN 15662 recommends a citrate extraction mix (Mix I), while AOAC standard 2007.1 uses an acetate extraction mix (Mix II).

For clean-up, the Diamino phase (PSA) removes, e.g., sugars and organic acids. MgSO₄ removes water, C₁₈ removes nonpolar interferences such as fats and the Carbon phase removes pigments, sterols, and nonpolar interferences.

For selection of the proper clean-up mix see table on opposite page.

For detailed instructions please visit www.mn-net.com or the original references at www.quechers.com.



Ordering information

	Volume	Description	Composition	REF	Pack of
	CHROMABOND® QuEChERS extraction buffer mixes				
	15 mL*	Mix I	citrate extraction mix	4 g MgSO ₄ , 1 g NaCl, 0.5 g Na ₂ H citrate · 1.5 H ₂ O, 1 g Na ₃ citrate · 2 H ₂ O	730970 50
	15 mL*	Mix II	acetate extraction mix	6 g MgSO ₄ , 1.5 g Na acetate	730971 50
	CHROMABOND® QuEChERS clean-up mixes				
	15 mL*	Mix III	Diamino clean-up mix	0.15 g CHROMABOND® Diamino with 0.9 g MgSO ₄	730972 50
	15 mL*	Mix IV	Diamino/Carbon clean-up mix	0.15 g CHROMABOND® Diamino with 0.9 g MgSO ₄ and 15 mg Carbon	730973 50
	15 mL*	Mix V	Diamino/Carbon clean-up mix	0.15 g CHROMABOND® Diamino with 0.9 g MgSO ₄ and 45 mg Carbon	730975 50
	15 mL*	Mix VI	Diamino/C ₁₈ ec clean-up mix	0.15 g CHROMABOND® Diamino with 0.9 g MgSO ₄ and 150 mg C ₁₈ ec	730974 50
	CHROMABOND® Diamino polypropylene columns				
	3 mL	adsorbent weight 200 mg		730561	50
	6 mL	adsorbent weight 500 mg		730562	30
	CHROMABOND® Diamino adsorbent				
				730653.20	20 g
				730653	100 g
	CHROMABOND® QuEChERS accessories				
	50 mL	polypropylene centrifuge tube with screw cap		730223	50

* 15 mL centrifuge tubes with screw cap (2 mL or 50 mL centrifuge tubes on request)

A number of custom-made QuEChERS mixes is available on request.

QuEChERS mixes

Sample property			
Low fat content (e.g., apples, strawberries)	Moderate content of chlorophyll and carotinoids (e.g., carrots, lettuce)	Higher content of chlorophyll and carotinoids (e.g., bell peppers, spinach)	Higher fat content (e.g., avocado)
CHROMABOND® QuEChERS extraction mixes			
Citrate or acetate extraction	Citrate or acetate extraction	Citrate extraction	Citrate extraction
Mix I or Mix II	Mix I or Mix II	Mix I	Mix I
CHROMABOND® QuEChERS clean-up mixes			
Diamino clean-up	Diamino/Carbon clean-up	Diamino/Carbon clean-up (higher Carbon content)	Diamino/C ₁₈ ec clean-up
Mix III	Mix IV	Mix V	Mix VI



Accessories for SPE

CHROMABOND® vacuum manifolds

- For simultaneous preparation of up to 12, 16 or 24 samples
- Replacement parts and accessories for special applications



Vacuum manifold for 12 columns

- 1 Rectangular glass cabinet; 2 sizes available: small for up to 12 CHROMABOND® columns or CHROMAFIX® cartridges; large for up to 16 CHROMABOND® LV columns or up to 24 CHROMABOND® columns or CHROMAFIX® cartridges (depending on lid)
- 2 Polypropylene lid
- 3 Vacuum gauge for pressure reading
- 4 Control valve for adjustment of vacuum
- 5 Replaceable valves for vacuum control of individual SPE columns
- 6 Variable rack with exchangeable partitions, which accept a wide variety of vessels like test tubes, measuring flasks, scintillation vials, autosampler vials, plastic vials etc.
- 7 CHROMABOND® LV columns with 15 mL sample reservoir for medium size samples
- 8 Polypropylene sample reservoirs (30 or 70 mL)
- 9 Adapter for sample reservoirs
- 10 CHROMABOND® tubing adapters

Full description and manual can be downloaded from www.mn-net.com

Ordering information

Description	Pack of	REF
Vacuum manifold complete		
consists of glass cabinet with lid and lid gasket, removable needles on lower side of lid, vacuum gauge, control valve, valves and caps, variable rack:		
for up to 12 columns or cartridges (including PP tank)	1	730150
for up to 16 LV columns	1	730360
for up to 24 columns or cartridges	1	730151
Glass cabinets without accessories (1)		
for 12 columns	1	730173
for 16 LV or 24 columns	1	730174
Lids with gaskets (2)		
for 12 columns (including Luer fittings and valves (5))	1	730175
for 16 LV columns (including Luer fittings and valves (5))	1	730365
for 24 columns (including Luer fittings and valves (5))	1	730176
Gaskets for lid, for 12 columns	2	730177
Gaskets for lid, for 24 columns	2	730178




Ordering information

Description	Pack of	REF
General accessories for vacuum manifolds		
Luer stoppers for vacuum manifold, blue	12	730194
Luer fittings for lid, female	12	730183.12
Luer fittings for lid, male	12	730184.12
Valves, plastic	12	730185
Stainless steel needles	12	730152
Polypropylene needles	12	730154
PP tanks for vacuum manifold for 12 columns (not available for 16- or 24-position manifold)	2	730233
Vacuum gauge, complete with accessories (3+4)	1	730179
Drying attachment and collecting racks for evaporation of eluates		
Drying attachment, with 12 positions (11)	1	730187
Drying attachment, with 16 positions	1	730990
Drying attachment, for 24 columns	1	730188
Collecting rack for 12 columns (6)	1	730157
Collecting rack for 16 LV columns	1	730366
Collecting rack for 24 columns	1	730153
Products for protection from cross contamination		
Valve, brass, tarnished	1	730189.1
Valves, as above	12	730189.12
Stainless steel connectors	12	730106
PTFE connectors	12	730564
Tubing adapters for application of large sample volumes (10)		
for 3 and 6 mL glass columns	4	730387
for 1, 3 and 6 mL polypropylene columns	4	730243
for 15, 45 and 70 mL polypropylene columns (PTFE tube length approx. 1 m)	4	730386

Protection from cross contamination

For special applications, which require maximum protection from cross contamination we supply chrome-plated brass valves and stainless steel or PTFE connectors, the application of which is shown below. These special connectors are fitted through the lid; thus the sample only has contact with the inert connector and can flow directly into the receptacle.



Drying attachment

If the eluate has to be evaporated, this can be performed with the so-called drying attachment (11, see below). This special lid has a gas connector on one side (12), from which the gas is fed simultaneously to the 12, 16, or 24 stations (13). Thus 12, 16, or 24 eluates can be evaporated simultaneously by just changing the lid and applying a stream of inert gas, e.g., nitrogen.





Accessories for SPE

CHROMABOND® empty columns and accessories

For individual packing of SPE columns with CHROMABOND® adsorbents

Ordering information

Description	Pack of	REF
Empty polypropylene columns with PE frits, 1 mL	100	730159
Empty polypropylene columns with PE frits, 3 mL	50	730160
Empty polypropylene columns with PE frits, 6 mL	30	730161
Empty polypropylene columns with PE frits, 15 mL	20	730230
Empty polypropylene columns with PE frits, 30 mL	20	730380
Empty polypropylene columns with PE frits, 45 mL	20	730355
Empty polypropylene columns with PE frits, 70 mL	20	730158
Empty polypropylene columns with PE frits, 150 mL	20	730474
PE frits for polypropylene columns 1 mL	250	730164
PE frits for polypropylene columns 3 mL	250	730162
PE frits for polypropylene columns 6 mL	250	730163
PE frits for polypropylene columns 15 mL	250	730351
PE frits for polypropylene columns 30 mL	250	730034
PE frits for polypropylene columns 45 mL	250	730356
PE frits for polypropylene columns 70 mL	250	730026
PE frits for polypropylene columns 150 mL	250	730475
Empty glass columns with glass fiber frits, 3 mL	50	730171
Empty glass columns with glass fiber frits, 6 mL	30	730172
Glass fiber frits for glass columns 3 mL	250	730191
Glass fiber frits for glass columns 6 mL	250	730192
Empty LV polypropylene columns with PE frits, 15 mL, for 100 mg adsorbent weight	50	732500
Empty LV polypropylene columns with PE frits, 15 mL, for 200/500 mg adsorbent weight	50	732501
PE frits for LV polypropylene columns 15 mL for 100 mg adsorbent weight	250	732019
PE frits for LV polypropylene columns 15 mL for 200/500 mg adsorbent weight	250	732020
Adapters (PVDF) for glass columns (3 and 6 mL)	4	730104.4
Adapters as above	10	730105
Adapters (PP) for polypropylene columns (1, 3 and 6 mL)	4	730100.4
Adapters as above	10	730101
Adapters (PE) for polypropylene columns (15, 45, 70 mL)	4	730350.4
Adapters as above	10	730385
Adapter (PE) for polypropylene columns (30 and 70 mL)	1	730566
Reservoir columns for application of medium-size samples		
Reservoir column 30 mL, polypropylene, with one adapter for 1, 3, 6 mL CHROMABOND® polypropylene columns	1	730102
10 Reservoir columns 30 mL, polypropylene with one adapter for 1, 3, 6 mL CHROMABOND® polypropylene columns	1 kit	730103
Reservoir column 70 mL, polypropylene, with one adapter for 1, 3, 6 mL CHROMABOND® polypropylene columns	1	730381
10 Reservoir columns 70 mL, polypropylene with one adapter for 1, 3, 6 mL CHROMABOND® polypropylene columns	1 kit	730382
Reservoir column 70 mL, polypropylene, with one adapter for 15, 45, 70 mL CHROMABOND® polypropylene columns	1	730388
10 Reservoir columns 70 mL, polypropylene with one adapter for 15, 45, 70 mL CHROMABOND® polypropylene columns	1 kit	730389



Automated and on-line SPE

Performing Solid Phase Extraction (SPE) manually can be time consuming and nerve-racking, especially when recovery and reproducibility are lacking due to sample variability. If SPE can be reliably automated, it becomes a much more efficient and reproducible process.

On-line SPE is a powerful method in automated sample preparation where the SPE hardware is technically integrated into a HPLC system. Crude samples are placed in an autosampler and processed fully automatic prior to injection into a GC (MS) or LC (MS) system. MN offers different on-line column configurations designed to fit your on-line SPE needs and filled with a choice of different adsorbents, modifications and particle sizes:

- Ready-to-use EC columns or ChromCart® cartridges for on-line SPE (standard dimensions 20 x 2 mm or 20 x 4 mm, resp.), filled with CHROMABOND® HR-Xpert phases (15 µm particles) or with NUCLEODUR® C₁₈ ec, C₈ ec, CN (20 µm particles)



EC columns



CC cartridges

- Columns for Gilson ASPEC™ systems are ready-to-use assembled with caps. In addition to the columns and phases listed below, all 1, 3 and 6 mL CHROMABOND® polypropylene columns from our program can be supplied assembled with ASP caps.



Columns for the Gilson ASPEC™

Ordering information

Gilson ASPEC™ columns

Column size	Weight [mg]	Pack of [columns]	REF
CHROMABOND® SiOH			
1 mL	100	100	730071ASP
3 mL	500	100	730073ASP
6 mL	1000	100	730075ASP

Other dimensions and adsorbents on request

- Special SPE columns equipped with caps and needles to be used in the SPE unit of the Gerstel MultiPurposeSampler (MPS), available in 1, 3, 6 mL.



SPE cartridges for Gerstel MPS system



Gerstel MPS system

Ordering information

Gerstel MPS columns

Column size	Weight [mg]	Pack of [columns]	REF
CHROMABOND® SiOH			
3 mL	200	50	730214MPS
3 mL	500	50	730073MPS
6 mL	1000	30	730075MPS
CHROMABOND® C₁₈ ec			
1 mL	100	100	730011MPS
3 mL	200	50	730012MPS
3 mL	500	50	730013MPS
CHROMABOND® HR-X			
1 mL	100	30	730935MPS
3 mL	200	30	730931MPS
6 mL	500	30	730939MPS



High-throughput SPE

CHROMABOND® MULTI 96 for robot systems

Alternatively CHROMABOND® MULTI 96 plates provide a means of high throughput sample preparation by processing 96 samples in a standard 8x12 microcolumn plate format compatible with standard 96-well plate liquid handling technologies and injection systems. MULTI 96 plates are available for solid phase extraction (SPE) and for filtration (see page 76).

CHROMABOND® MULTI 96 - SPE in microtiter format

- 96-well PP microtiter plates with PE filter elements
- Cavity volume 1.5 mL
- Adsorbent weights 10, 25, 50, 100 mg per micro-column
- Supplied with any CHROMABOND® SPE adsorbents
- For simultaneous preparation of 96 samples
- Easy method transfer from CHROMABOND® columns or CHROMAFIX® cartridges to CHROMABOND® MULTI 96


Advantages of this high-throughput system:

- Simultaneous preparation of 96 samples; this means a 4-fold increase over traditional 24-position SPE processors
- Economical by saving time and solvent
- Use of multi-channel pipettors facilitates liquid transfer steps
- Readily adaptable to all common automated and robotic handling systems
- Minimized dead volume ($\leq 40 \mu\text{L}$)


Instrument compatibility

CHROMABOND® MULTI 96 SPE microtiter plates as well as CHROMAFIL® MULTI 96 filtration plates are compatible with, e.g., the following liquid handling and SPE automation systems:

- Perkin Elmer MultiProbe® II
- Tomtec Quadra 3® and Quadra 3® SPE
- Hamilton Microlab® SPE Workstation
- Beckman Coulter Biomek® 2000
- Caliper Life Science RapidTrace®
- Gilson ASPEC™ XL4 and ASPEC™ XL
- Gilson 215 SPE Liquid Handler
- Tecan Genesis™ FE500
- Eppendorf epMotion®



Multiprobe® II (Perkin-Elmer)



Biomek® 2000 (Beckman Coulter)



CHROMABOND® MULTI 96 vacuum manifold

For handling of CHROMABOND® MULTI 96 SPE plates for up to 96 samples

CHROMABOND® MULTI 96 is designed for use in common robotic workstations or commercially available liquid handling systems. Alternatively, use of multi-channel pipetters facilitates a manual liquid transfer. Extraction is carried out using the CHROMABOND® MULTI 96 vacuum manifold. With the help of the control valve the vacuum of the manifold can be adjusted leading to an optimum flow rate through the CHROMABOND® MULTI 96 SPE plate.




A reservoir tank and 96-well collection plates (96 x 0.5 or 96 x 2 mL) made of polypropylene can be supplied as accessories. An interesting alternative for collection of the eluates is a collection rack, which can be fitted with twelve 8-well strips of polypropylene tubes (each 1 mL). If you have to work on less than 96 samples, you can seal individual rows of the 96-well plate with a PTFE-covered rubber pad.



Ordering information

Description	Pack of	REF
CHROMABOND® MULTI 96 vacuum manifold with reservoir tank, vacuum gauge, and control valve	1	738630.M
96-well microtiter plates (polypropylene) 96 x 0.25 mL	10	738651
96-deep-well collecting plate (polypropylene) 96 x 2 mL	5	738650.5
Collection racks with polypropylene tube strips (twelve 8-well strips) 96 x 1.0 mL	5	738637
Polypropylene tube strips (twelve 8-well strips) 96 x 1.0 mL	10	738652
8-well strip sealing caps for PP tube strips (REF 738652)	30	738638
Reservoir tanks (polypropylene)	2	738639.M
Butyl rubber pad, PTFE covered for sealing of individual rows of the 96-well plate, 125 x 85 mm	1	738645



For CHROMAFIL® MULTI 96 filter plates see page 76. The ordering information of 96-well plates packed with individual CHROMABOND® adsorbents is listed with the respective phases.



Kieselguhr phase for liquid-liquid extraction

CHROMABOND® XTR

for liquid-liquid extraction

- ✦ Base material coarse-grained kieselguhr (also known as diatomaceous earth, hydromatrix, celite)
Large pore size, high pore volume, constantly high batch-to-batch quality
pH working range 1-13
- ✦ **Application:**
liquid-liquid extraction of highly viscous aqueous solutions such as physiological fluids (blood, plasma, and serum) in clinical chemistry, dyes in textiles, environmental and food analysis without use of a separation funnel
High water loadability without breakthrough of water during elution with organic solvents
also suited for removing small amounts of water from solvents which are not miscible with water
- ✦ **Advantages:**
fast, reproducible and economical
simultaneous preparation of several samples
no problems with phase separation · no formation of emulsions
high recovery rates
saving of time and solvents
organic solutions need not to be dried after separation

Liquid-liquid extraction

Solvents applicable for elution

- ✓ diethyl ether
- ✓ tert-butyl methyl ether
- ✓ ethyl acetate
- ✓ *n*-hexane
- ✓ cyclohexane
- ✓ toluene
- ✓ dichloromethane (methylene chloride)
- ✓ trichloromethane (chloroform)
- ✓ trichloromethane – methanol (90:10, v/v)
- ✓ trichloromethane – methanol (85:15, v/v)
- ✓ diethyl ether – ethanol (90:10, v/v)
- ✓ diethyl ether – ethanol (80:20, v/v)
- ✓ dichloromethane – 2-propanol (90:10, v/v)
- ✓ dichloromethane – 2-propanol (85:15, v/v)

Eluents with too high alcohol contents cause an increase in volume of the aqueous phase on the CHROMABOND® XTR. Here the column could be overloaded and the aqueous phase displaced from the column. In this case, a greater capacity column should be used.

Depending on the concentration of the analytes eluates can be analyzed immediately, or the organic solvent is evaporated. The pH value of the aqueous solution can be altered on the column, which enables elution of different compounds of a sample under optimized conditions. Under certain circumstances, acidic, neutral, and basic compounds can be fractionated in this way.

General column parameters

CHROMABOND® XTR Volume	Amount of adsorbent	Max. volume capacity of aq. solution	Waiting period before elution	Elution volume
1 mL	250 mg	0.25 mL	5 min	3 mL
3 mL	500 mg	0.5 mL	5 min	6 mL
6 mL	1 g	1 mL	5-10 min	8 mL
15 mL	3 g	3 mL	5-10 min	12 mL
30 mL	4.5 g	5 mL	5-10 min	16 mL
45 mL	8.3 g	10 mL	10-15 min	24 mL
70 mL	14.5 g	20 mL	10-15 min	40 mL
150 mL	37.5 g	50 mL	10-15 min	90 mL



Sample application



Adsorption of the sample



Sample elution



Determination of azo dyes and aromatic amines in colored textile materials with reference to § 64 LFGB (formerly § 35 LMBG)

Sample pretreatment:

Weigh about 1 g cut-up textile sample (colored textiles about 0.1 g) in a 100 mL threaded vial. (Degrease leather samples before processing: cover sample with technical purity *n*-hexane and put the vial in an ultrasonic bath for 20 min. After decanting the *n*-hexane rinse with little *n*-hexane and dry sample by gentle heating and blowing with air or N₂.)

Add 250 µL internal standard (IS: 1.2 mg/mL tetramethylbenzidine in methanol – ethyl acetate (1:1, v/v)), 17.0 mL citrate buffer (pH 6) (25.05 g citric acid and 12.64 g NaOH, fill up with deionized water to 2 L) and heat 30 min at 70 °C. Then add 3 mL of a freshly prepared solution of 0.2 g/mL sodium dithionite in water and heat for exactly 30 min to 70 °C while shaking occasionally.

Sample application:

Cool the solution immediately (put vial in water – stopping of reductive cleavage). After 5–10 min pour it onto the CHROMABOND® XTR column (squeeze textile remains).






Elution:

Allow solution to be soaked up by the adsorbent for 15 min. Then elute four times with 20 mL each of diethyl ether or diethyl ether – ethanol (90:10, v/v) (depending on recovery rates), using the first 40 mL to rinse the sample remains. Evaporate eluates to 3 mL with a rotation evaporator and transfer the solution into a 10 mL measuring flask using a pasteur pipette and rinsing with methanol. Fill up to the marking with methanol, shake, and pipette about 1 mL into a vial.

Further analysis: Fast GC on OPTIMA® δ-3, 10 m, 0.1 mm ID, 0.1 µm film, REF 726410.10 (application 210820) or HPLC on NUCLEOSIL® 100-5 C₁₈ HD (application 110500 at www.mn-net.com/apps)

MN Appl. No. 302100

Ordering information

	Column volume	1 mL	3 mL	6 mL	15 mL	30 mL	45 mL	70 mL	150 mL
	Adsorbent weight	250 mg	500 mg	1 g	3 g	4.5 g	8.3 g	14.5 g	37.5 g
	Max. volume capacity of aqueous solution	0.25 mL	0.5 mL	1 mL	3 mL	5 mL	10 mL	20 mL	50 mL
	Pack of	100	50	30	30	30	30	30	10
	CHROMABOND® XTR polypropylene columns								
	730501	730502	730487	730489	730505	730506	730507	730509	
	CHROMABOND® XTR polypropylene columns · BIGpacks								
		730487.250 (250 col.)					730507.100 (100 col.)		
	CHROMABOND® MULTI 96 XTR								
	96–well plates 96 x 150 mg, packs of 1 plate, for max. 96 x 0.2 mL aqueous solution								
		738131.150M							
	CHROMABOND® XTR adsorbent								
	50 bags of 14.5 g, for max. 20 mL aqueous solution each)								
	for 70 mL PP columns with 100 PE filter elements		for NT20 with 50 PE filter elements (10 mm dia.)		500 g		1 kg		5 kg
	730585	730586	730595.500		730595.1000		730595.5000		
	Accessories for liquid–liquid extraction with CHROMABOND® XTR								
	variable polypropylene rack for 24 positions, incl. 24 PP stopcocks and 24 PP needles								730508

For parallel processing of up to 24 CHROMABOND® XTR columns 1–150 mL we recommend the polypropylene rack REF 730508 consisting of:

two side walls (1), middle part including stopcocks and needles (2), bottom part (3), top part for stabilizing 45 mL and 70 mL CHROMABOND® XTR columns (4).

This rack can be adjusted to various heights depending on the CHROMABOND® XTR columns and the collection vials used. Each position of the middle part is equipped with a polypropylene stopcock on the top (REF 730185) and a polypropylene needle on the bottom (REF 730154).

For collection of the sample, vessels such as vials, test tubes, round bottom or tapered flasks, can be used. For our program of sample vials, please see the chapter "Vials and accessories" from page 77.





Columns for gravity flow phase separation

CHROMABOND® PTS and PTL

columns for phase separation

- ❖ Automatic separation of a two-phase mixture without separation funnel
 Two-phase mixtures are completely applied to the column and the phase boundary is determined without further work. The special membrane automatically stops the flow when the lower phase has passed. The upper phase remains in the column, thus both phases are available for further analysis.
 Columns **must not** be run with vacuum or pressure
- ❖ **PTS**
 for solvents **heavier** than water, e.g., trichloromethane, dichloromethane
 maximum size 150 mL
- ❖ **PTL**
 for solvents **lighter** than water, e.g., diethyl ether, hexane
 maximum size 70 mL

Ordering information

Column volume [mL]	Pack of [columns]	REF
CHROMABOND® PTS		
for solvents heavier than water		
1	100	730710
3	100	730712
6	100	730714
15	100	730716
30	100	730718
45	50	730720
70	50	730722
150	20	730724
CHROMABOND® PTL		
for solvents lighter than water		
1	100	730730
3	100	730732
6	100	730734
15	100	730736
30	100	730738
45	50	730740
70	50	730742



Ideal tool for breaking emulsions



CHROMABOND® PTL in action: organic upper phase (colorless), aqueous lower phase (red)



Glass columns and accessories for Flash chromatography

- Economic low-tech method for the synthesis laboratory
Suited for the separation of compounds up to gram levels
No expensive equipment required

- MN flash chromatography kits include a glass column, eluent reservoir, silica 60 and accessories.
Glass columns of different sizes and accessories can be ordered separately.

These columns are normally filled to a height of about 15 cm, working pressures are 1.5 to 2 bar.

The most used adsorbent is silica 60 with particle size 40–63 μm (see page 204), however, you may also use our ranges of other LC adsorbents and of POLYGOPREP silica phases (see page 203). Particle sizes < 25 μm should only be used with very low-viscosity mobile phases, because otherwise flow rates will be very low.

These columns are to be packed by the user.



Ordering information

Designation	Pack of	REF
Flash chromatography kits		
Flash chromatography kit I, consists of 1 glass column 20 mm ID x 400 mm length, one 1-L eluent reservoir, 100 g silica 60 (40–63 μm), sea sand, silanized glass fiber wadding, 1 m PTFE tubing	1 kit	727450
Flash chromatography kit II, consists of 1 glass column 40 mm ID x 450 mm length, one 2-L eluent reservoir, 100 g silica 60 (40–63 μm), sea sand, silanized glass fiber wadding, 1 m PTFE tubing	1 kit	727451
Flash chromatography columns		
complete with adapter and PTFE tap, fitted with a polyethylene net to protect against bursting		
20 mm ID x 200 mm length	1 column	727400
20 mm ID x 400 mm length	1 column	727401
25 mm ID x 200 mm length	1 column	727402
25 mm ID x 400 mm length	1 column	727403
30 mm ID x 300 mm length	1 column	727404
30 mm ID x 400 mm length	1 column	727405
40 mm ID x 300 mm length	1 column	727406
40 mm ID x 450 mm length	1 column	727407
Accessories for flash chromatography glass columns		
1-L eluent reservoir with adapter, covered with a protective plastic sleeve for burst protection; this also prevents build-up of UV-induced radicals in the eluent	1	727420
2-L eluent reservoir as above	1	727421
Pressure gauge for controlling flow rates	1	727422
PTFE tubing, 3 mm OD, 2 mm ID, length 1 m	1 m	727424
Sea sand, acid washed and calcined	1000 g	727423
Glass fiber wadding, silanized	25 g	718002



CHROMABOND® Flash RS cartridges

CHROMABOND® Flash RS cartridges


ideal for Flash separations from 10 mg up to 160 g

- ❖ **For convenient operation and reliable upscaling**
 the complete program of ready-to-use Flash cartridges for the ISCO® Companion® and other Teledyne Isco CombiFlash® systems, or as stand-alone version for all pump/detector combinations, e.g., from Biotage®, Büchi®
 Adsorbent weights of 4 g to 1600 g from one of the leading companies in silica and TLC business
- ❖ **Increases flexibility**
 considerable program of different phases and modifications
- ❖ **Saves time and money**
 convenient prices, short delivery times
- ❖ **Increases analytical safety**
 high pressure stability of 15 bar/220 psi (12 bar for cartridges > 200 g), excellent separation efficiency, good reproducibility



Technical features

- ❖ **Distribution of eluent stream** via highly porous frits
- ❖ **Cartridge material and geometry:**
 organic solvent resistant, low bleed polypropylene, thick column walls, one piece body, sophisticated length to diameter ratio for high plate numbers and excellent separation efficiencies
- ❖ **Column connections**
 CHROMABOND® RS cartridges are 100% compatible with the ISCO® Companion®, no additional hardware is needed for this type of purification systems.
 CHROMABOND® RS cartridges (except RS 800 and RS 1600 with Maxi Luers) can also be used as stand alone system with any pump/detector/fraction collector combination using the CHROMABOND® Flash Starter Kit (REF 730798) or the CHROMABOND® Flash Stand Alone Kit (REF 732903).



Accessories for CHROMABOND® Flash columns · Ordering information

Description	Pack of	REF
CHROMABOND® Flash Starter Kit		
consists of 1/8" PTFE tubing, 1.5 mm ID, 3 m long; 5 x 1/4"-28 PP nuts; 5 x 1/8" ETFE ferrules; 5 x 1/4"-28 nylon unions; 2 x 1/4"-28 PP Luer lock, female; 1 x 1/4"-28 PP Luer lock, male; 1 x 1/4"-28 PP Luer tip, male	1 kit	730798



Description	Pack of	REF
CHROMABOND® Flash Stand Alone Kit, Luer		
consist of 1 x 1/4"-28 PP Luer lock, female; 1 x 1/4"-28 PP Luer lock, male; 2 x 1/8" ETFE ferrules; 2 x 1/4"-28 nylon unions; 2 x 1/4"-28 PP nuts	1 kit	732903

CHROMABOND® Flash solutions for Flash instruments

- Product range designed for use in the Teledyne Isco CombiFlash® systems (Companion®, Rf etc.) and Flash systems of Biotage® AB (FlashMaster™) without additional connectors or capillaries

On request all column types listed below can be packed with any adsorbent as described on page 8 (please note that other packings often result in differing adsorbent weights).

Ordering information

Designation	Column length [cm]	ID [mm]	Adsorbent weight [g]	Pack of [columns]	REF
CHROMABOND® Flash RS columns for Teledyne Isco® systems					
All CHROMABOND® Flash RS cartridges can be directly used in the Teledyne Isco Companion®, Rf, etc.					
CHROMABOND® Flash RS 4 SiOH	9.8	12.4	4	20	732800
CHROMABOND® Flash RS 15 SiOH	11.6	21.2	15	20	732801
CHROMABOND® Flash RS 25 SiOH	16.5	21.2	25	15	732802
CHROMABOND® Flash RS 40 SiOH	17.1	26.4	40	15	732803
CHROMABOND® Flash RS 80 SiOH	24.0	30.8	80	12	732804
CHROMABOND® Flash RS 120 SiOH	25.5	36.0	120	10	732805
CHROMABOND® Flash RS 200 SiOH	20.0	60.0	200	6	732806
CHROMABOND® Flash RS 330 SiOH	27.0	60.0	330	4	732807
CHROMABOND® Flash RS 800 SiOH	38.5	82.0	800	2	732808
CHROMABOND® Flash RS 1600 SiOH	43.0	104.0	1600	2	732809
Corresponding TLC plates: silica, see page 213					
CHROMABOND® Flash RS 4 C ₁₈ ec	9.8	12.4	4.3	2	732810
CHROMABOND® Flash RS 15 C ₁₈ ec	11.6	21.2	16.4	1	732811
CHROMABOND® Flash RS 25 C ₁₈ ec	16.5	21.2	26	1	732812
CHROMABOND® Flash RS 40 C ₁₈ ec	17.1	26.4	43	1	732813
CHROMABOND® Flash RS 80 C ₁₈ ec	24.0	30.8	86	1	732814
CHROMABOND® Flash RS 120 C ₁₈ ec	25.5	36.0	130	1	732815
CHROMABOND® Flash RS 200 C ₁₈ ec	20.0	60.0	220	1	732816
CHROMABOND® Flash RS 330 C ₁₈ ec	27.0	60.0	360	1	732817
CHROMABOND® Flash RS 800 C ₁₈ ec	38.5	82.0	880	1	732818
CHROMABOND® Flash RS 1600 C ₁₈ ec	43.0	104.0	1760	1	732819
Corresponding TLC plates: RP-18 W/UV ₂₅₄ , see page 220					
CHROMABOND® Flash RS 4 NH ₂	9.8	12.4	4.3	2	732820
CHROMABOND® Flash RS 15 NH ₂	11.6	21.2	16.4	1	732821
CHROMABOND® Flash RS 25 NH ₂	16.5	21.2	26	1	732822
CHROMABOND® Flash RS 40 NH ₂	17.1	26.4	43	1	732823
CHROMABOND® Flash RS 80 NH ₂	24.0	30.8	86	1	732824
CHROMABOND® Flash RS 120 NH ₂	25.5	36.0	130	1	732825
CHROMABOND® Flash RS 200 NH ₂	20.0	60.0	220	1	732826
CHROMABOND® Flash RS 330 NH ₂	27.0	60.0	360	1	732827
Corresponding TLC plates: Nano-SIL NH ₂ , see page 222					



CHROMABOND® Flash RS cartridges

Designation	Column length [cm]	ID [mm]	Adsorbent weight [g]	Pack of [columns]	REF
CHROMABOND® Flash RS 4 OH (Diol)	9.8	12.4	4.3	2	732830
CHROMABOND® Flash RS 15 OH (Diol)	11.6	21.2	16.4	1	732831
CHROMABOND® Flash RS 25 OH (Diol)	16.5	21.2	26	1	732832
CHROMABOND® Flash RS 40 OH (Diol)	17.1	26.4	43	1	732833
CHROMABOND® Flash RS 80 OH (Diol)	24.0	30.8	86	1	732834
CHROMABOND® Flash RS 120 OH (Diol)	25.5	36.0	130	1	732835
CHROMABOND® Flash RS 200 OH (Diol)	20.0	60.0	220	1	732836
CHROMABOND® Flash RS 330 OH (Diol)	27.0	60.0	360	1	732837
Corresponding TLC plates: Nano-SIL DIOL, see page 223					
CHROMABOND® Flash RS 4 CN	9.8	12.4	4.3	2	732840
CHROMABOND® Flash RS 15 CN	11.6	21.2	16.4	1	732841
CHROMABOND® Flash RS 25 CN	16.5	21.2	26	1	732842
CHROMABOND® Flash RS 40 CN	17.1	26.4	43	1	732843
CHROMABOND® Flash RS 80 CN	24.0	30.8	86	1	732844
CHROMABOND® Flash RS 120 CN	25.5	36.0	130	1	732845
CHROMABOND® Flash RS 200 CN	20.0	60.0	220	1	732846
CHROMABOND® Flash RS 330 CN	27.0	60.0	360	1	732847
Corresponding TLC plates: Nano-SIL CN, see page 221					
CHROMABOND® Flash RS 4 Alox A	9.8	12.4	8	20	732870
CHROMABOND® Flash RS 4 Alox N	9.8	12.4	8	20	732871
CHROMABOND® Flash RS 4 Alox B	9.8	12.4	8	20	732872
CHROMABOND® Flash RS 15 Alox A	11.6	21.2	30	20	732874
CHROMABOND® Flash RS 15 Alox N	11.6	21.2	30	20	732873
CHROMABOND® Flash RS 15 Alox B	11.6	21.2	30	20	732875
CHROMABOND® Flash RS 25 Alox A	16.5	21.2	50	15	732876
CHROMABOND® Flash RS 25 Alox N	16.5	21.2	50	15	732877
CHROMABOND® Flash RS 25 Alox B	16.5	21.2	50	15	732878
CHROMABOND® Flash RS 40 Alox A	17.1	26.4	80	15	732879
CHROMABOND® Flash RS 40 Alox N	17.1	26.4	80	15	732881
CHROMABOND® Flash RS 40 Alox B	17.1	26.4	80	15	732880
Corresponding TLC plates: Alox, see page 224					

CHROMABOND® Flash columns for Biotage® FlashMaster™ systems

CHROMABOND® Flash FM 15/2 SiOH	9.0	15.8	2.0	50	730881
CHROMABOND® Flash FM 25/5 SiOH	10.0	20.5	5.0	50	730891
CHROMABOND® Flash FM 25/10 SiOH	10.0	20.5	10.0	50	730666
CHROMABOND® Flash FM 70/10 SiOH	15.4	26.8	10.0	30	730885
CHROMABOND® Flash FM 70/20 SiOH	15.4	26.8	20.0	30	730915
CHROMABOND® Flash FM 70/25 SiOH	15.4	26.8	25.0	30	730892
CHROMABOND® Flash FM 150/25 SiOH	17.0	38.2	25.0	20	730667
CHROMABOND® Flash FM 150/50 SiOH	17.0	38.2	50.0	20	730887
CHROMABOND® Flash FM 150/70 SiOH	17.0	38.2	70.0	10	730880
CHROMABOND® Flash FM 15/2 C ₁₈ ec	9.0	15.8	2.0	50	730890
CHROMABOND® Flash FM 25/5 C ₁₈ ec	10.0	20.5	5.0	20	730884
CHROMABOND® Flash FM 70/10 C ₁₈ ec	15.4	26.8	10.0	20	730886
CHROMABOND® Flash FM 150/50 C ₁₈ ec	17.0	38.2	50.0	10	730888
CHROMABOND® Flash FM 70/10 NH ₂	15.4	26.8	10.0	20	730768
CHROMABOND® Flash FM 70/20 NH ₂	15.4	26.8	20.0	20	730767



Technical support

Loadability

Due to the narrow particle size distribution, the excellent packing quality and the optimized stationary phases (acid washed silica, reduced particulate matter) our cartridges can realize highest loadability at best possible separation efficiency. Additionally, the large range of different cartridge lengths and diameters eases to find the optimum in loadability for a given sample amount.

Rule of thumb for the loadability

Separation	Loadability	g sample / g adsorbent
difficult	low	≤ 1%
easy	high	≥ 10%

Loadability table CHROMABOND® Flash RS

SiOH cartridge	Average loadability per cartridge [g]	
	difficult separation	easy separation
RS 4	0.04	0.4
RS 15	0.15	1.5
RS 25	0.25	2.5
RS 40	0.4	4
RS 80	0.8	8
RS 120	1.2	12
RS 200	2	20
RS 330	3.3	33
RS 800	8	80
RS 1600	16	160

Upscaling of the optimum flow rate


This depends on the eluent and the separation problem.

For RS cartridges the upscaling relation is simple:

The silica weight in g is proportional to the flow rate (for equal eluent polarity), e.g.,

4 g silica → optimum flow: ~ 6–12 mL/min

40 g silica → optimum flow: ~ 60–120 mL/min



Back pressure and pressure stability

The back pressure always depends on flow rate and viscosity of the eluent mixture, column length and diameter and the particle size. The high performance CHROMABOND® Flash RS cartridges up to 200 g silica

are stable up to 15 bar (220 psi, > 200 g: 12 bar). We recommend using a pressure guard, because short time pressure peaks (viscosity of eluent or gradient changes) can exceed the pressure limit.

Back pressure of CHROMABOND® Flash RS SiOH cartridges (eluent hexane – ethyl acetate 9:1 or 8:2)

Cartridge	Flow rate						
	20 mL/min	40 mL/min	80 mL/min	120 mL/min	160 mL/min	200 mL/min	240 mL/min
RS 4	0.75 bar	1.5 bar					
RS 15	0.25 bar	0.75 bar	1.5 bar	2.0 bar			
RS 25	0.5 bar	1.0 bar	1.75 bar	3.0 bar	4.0 bar	5.0 bar	
RS 40		0.75 bar	1.5 bar		3.0 bar		3.5 bar
RS 80			1.5 bar	2.5 bar	3.0 bar	3.5 bar	4.0 bar
RS 120			1.0 bar	1.5 bar	2.0 bar	2.5 bar	3.0 bar
RS 200			1.0 bar		2.0 bar		3.0 bar
RS 330			1.5 bar		3.0 bar		4.0 bar

Conditioning volumes for CHROMABOND® Flash RS cartridges (normally 1.5 column volumes of the eluent)

Cartridge	Volume of eluent for conditioning	Cartridge	Volume of eluent for conditioning
RS 4	20 mL	RS 120	440 mL
RS 15	60 mL	RS 200	750 mL
RS 25	90 mL	RS 330	1100 mL
RS 40	140 mL	RS 800	2900 mL
RS 80	280 mL	RS 1600	5000 mL



CHROMABOND® Flash RS cartridges

TLC upscaling

TLC is often used for the development of a selective and reproducible method in Flash chromatography, because it is often necessary to test a large number of eluent and/or adsorbent combinations. MN TLC plates and sheets are coated with the same base silica, which is used in our CHROMABOND® Flash cartridges. This is an important prerequisite for the reproducible transfer

of a TLC separation to the Flash column, because the parameters are identical in both systems.



MN TLC and Flash product ensure:

- Same selectivity and easy upscaling from TLC to Flash separations
- Saving time and money, because expensive optimizations are not required

Examples for transfer of a TLC separation to a Flash column:

R_f values of the TLC separation should be in the range of 0.1–0.4 (low height).

ΔR_f values on the TLC plate should be as high as possible.



During TLC optimization always use solvents, which are well suited for the following Flash chromatography!

MN adsorbents

a unique variety of phases


- As with our SPE products, all Flash columns and cartridges from MN are available with our whole range of CHROMABOND® phases (more than 40, e.g., C₁₈, C₈, OH, Alox as listed on page 8)
- Additionally you can choose from our range of POLYGOPREP silica packings in particle sizes from 20 to 130 μm and pore sizes from 60 to 4000 Å (see page 203).
- For high performance Flash separations you can order columns packed with spherical NUCLEODUR® featuring very high separation efficiency and extremely increased column lifetime (particle size > 20 μm as listed on page 198)

Technical silica information

Specification of modified and plain silica:
acid-washed irregular silica, pore size 60 Å, particle size 45 μm , specific surface 500 m²/g, pH stability 2–8

Additionally available silicas and particle sizes:

- Irregular POLYGOPREP silica with particle sizes of 20 to 130 μm and pore sizes of 60 to 4000 Å
- Spherical high performance silica (NUCLEODUR®, 110 Å) with particle sizes of 20 or 30 μm for high separation efficiency and very long column life



Comparison of separation efficiency and price of irregular versus spherical silica




Separation efficiency and reproducibility

Our optimized automatic packing process leads to an excellent packing quality, irrespective of the phase or particle size distribution (normal phase or reversed phase, spherical or irregular particles).

MN, as a manufacturer of silicas, has decades of experience in the production of first class separation phases and columns. This leads to highest separation efficiencies of the columns, a constant back pressure (via controlled narrow particle size distribution) and good reproducibilities from cartridge to cartridge.

The separation efficiency is in the first step not influenced by the dimension or the geometry of the Flash RS cartridges. The chromatograms below show an identical resolution and peak form for different column dimensions, when flow and sample amount is adjusted correctly. This is advantageous for optimization and upscaling experiments.





Syringe filters CHROMAFIL®



Sample clarification

Syringe filters are used for filtration of suspended matter from liquid samples or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Special manipulations are not required. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment.

Advantages:

✦ Polypropylene housing

Considerably better solvent stability compared to acrylate and polystyrene filters, featuring a low content of extractable substances

✦ Lowest content of extractable substances

The housing of every CHROMAFIL® filter is **ultrasonically sealed (welded), not glued**, because glue may have extractable ingredients. Welding leads to a tight connection between both parts, thus the filter can be used in both directions. The special **thick rim** of the housing is ideal for use in laboratory robots (e.g., SOTAX®, Benchmate™).

✦ Luer lock on the side of entry

For a safe connection on the high-pressure side every filter provides a Luer lock on the side of entry.

✦ Luer exit

For 25 and 3 mm filters: standard Luer exit
For 15 mm filters: minispikes · This Luer configuration offers a low hold-up volume and easy filtration into autosampler vials and NMR tubes.

With the aid of a special adapter, filter inlet and filter exit can be fitted to all CHROMABOND® columns and accessories for selective sample preparation.

✦ No rupture of membrane due to the impact plate

The input solvent stream is broken and distributed by the impact plate, and does not directly hit the membrane: this prevents rupture of the membrane. The high pressure stream is diverted into four lanes.

✦ Optimum flow geometry because of the star-shaped distribution device

The stream of liquid is broken into 4 lanes by the impact plate and then further distributed to 8 slots in the form of a star connected with 5 or 8 circular channels (for 15 mm and 25 mm filters, respectively). Thus, the fluid is able to penetrate the membrane on the whole surface, not only on a small region; the filter is not plugged up rapidly, which results in a high flow efficiency.

✦ Color coded filters

Filters with 0.2 µm pores have a yellow upper shell, that of filters with 0.45 µm pores is colorless; the different membrane types are distinguished by different colors of the lower shell.

✦ Different pore sizes for versatile filtration

Standard pore sizes 0.2 and 0.45 µm (additionally: PET filters with 1.2 µm, glass fiber filters with 1 µm, PES filters with 5 µm). Filters with 0.45 µm pore size efficiently remove fine particles that can plug chromatography columns. Filters with 0.2 µm pore size are excellent for filtration of UHPLC samples or other techniques requiring high purity samples.

✦ Filter sizes

25, 15 and 3 mm diameter: the small diameter filters are especially recommended for very small samples, which require extremely low dead volumes: 5 µL for 3 mm Ø, 35 µL for 15 mm Ø, 80 µL for 25 mm Ø

Recommended filter size depending on sample volume

Sample volume	Recommended filter diameter
≤ 1 mL	3 mm
1–5 mL	15 mm
5–100 mL	25 mm

Filters can be **autoclaved** at 121 °C, 1.1 bar for 30 min. All 25 mm CHROMAFIL® filters are designed to be 100% compatible and reliable for use with the SOTAX® AT70 smart fully automated dissolution testing systems.




Depending on your filtration task you can choose filter membranes made from different materials:

Material	Page
Combi Filters with glass fiber prefilters	
Polyester (GF/PET)	68
Regenerated cellulose (GF/RC)	68
Polyvinylidene difluoride (GF/PVDF)	68
Syringe filters without prefilters	
Polyester (PET)	69
Regenerated cellulose (RC)	69
Polytetrafluoroethylene (PTFE)	70
Cellulose mixed esters (MV)	70
Cellulose acetate (CA) · sterile and non-sterile	71
Polyamide / Nylon (PA)	72
Polyethersulfone (PES)	71
Polyvinylidene difluoride (PVDF)	72
Glass fiber (GF)	73

CHROMAFIL® BIG-BOX

- 400 (25 mm) or 800 (15 mm) color-coded quality syringe filters · 400 labeled Xtra syringe filters
- Food safe PE box with screw cap
- Economical prices

CHROMAFIL® Combi filters



Combi syringe filters with a coarse glass fiber prefilter and a small-pore membrane as main filter

User benefits:

- For solutions with a high load of particulate matter: lower back pressure, easy filtration
- For high yields of filtrate: more mL of pure filtrate per filter

The technology:

The glass fiber membrane (1.0 µm) removes coarse particles, before they can block the fine main membrane. This results in a better filtration efficiency, especially for highly contaminated samples.

Housing:	Solvent-resistant, ultra low bleed polypropylene
Inlet:	Luer lock
Exit:	Luer
Pore diameter:	1.0 / 0.20 µm or 1.0 / 0.45 µm
Filter diameter:	25 mm
Void volume:	< 80 µL
Packing unit:	100 filters; BIG-BOX with 400 filters

CHROMAFIL® Xtra

labeled for method validation and certification

- Xtra:** imprint for direct identification of the membrane type, diameter and pore size
- Xtra:** low bleeding PP housing
- Xtra:** color-free plain polypropylene





CHROMAFIL® Combi filters

Polyester with glass fiber prefilter (GF/PET)

- Hydrophilic multipurpose membrane for polar as well as nonpolar solvents
The HPLC filter with glass fiber prefilter, especially suited for mixtures of water and organic solvents
- Recommended for solutions with a high load of particulate matter or for highly viscous solutions



Ordering information

Type	Pore size [μm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
			top	bottom	filters/pack	REF	filters/pack	REF
GF/PET-20/25	1.0/0.20	25	blue	orange	100	729032	400	729032.400
GF/PET-45/25	1.0/0.45	25	black	orange	100	729033	400	729033.400

Regenerated cellulose with glass fiber prefilter (GF/RC)

- Hydrophilic membrane for aqueous and organic-aqueous liquids, i.e. polar and medium polar sample solutions
- Recommended for solutions with a high load of particulate matter or for highly viscous aqueous solutions



Ordering information

Type	Pore size [μm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
			top	bottom	filters/pack	REF	filters/pack	REF
GF/RC-20/25	1.0/0.20	25	blue	blue	100	729050	400	729050.400
GF/RC-45/25	1.0/0.45	25	black	blue	100	729051	400	729051.400

Polyvinylidene difluoride with glass fiber prefilter (GF/PVDF)

- Hydrophilic membrane
- Recommended for filtration of biological samples with high particle loads. This filter features a high binding capacity for proteins.
- Also suited for filtration of polar and nonpolar solutions



Ordering information

Type	Pore size [μm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
			top	bottom	filters/pack	REF	filters/pack	REF
GF/P-45/25	1.0/0.45	25	black	white	100	729039	400	729039.400



Polyester (PET)

- Hydrophilic multipurpose membrane for polar as well as nonpolar solvents
- The HPLC filter**, especially suited for mixtures of water and organic solvents
- For TOC/DOC determination
- Not cytotoxic, does not inhibit the growth of microorganisms and higher cells



Ordering information - CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
				filters/pack	REF	filters/pack	REF
PET-20/25	0.20	25	labeled	100	729221	400	729221.400
PET-45/25	0.45	25	labeled	100	729220	400	729220.400
PET-120/25	1.2	25	labeled	100	729229	400	729229.400

Ordering information - CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
			top	bottom	filters/pack	REF	filters/pack	REF
PET-20/15 MS	0.20	15	yellow	orange	100	729022	800	729022.800
PET-45/15 MS	0.45	15	colorless	orange	100	729023	800	729023.800
PET-20/25	0.20	25	yellow	orange	100	729021	400	729021.400
PET-45/25	0.45	25	colorless	orange	100	729020	400	729020.400

MS = minispikes on filter exit

Regenerated cellulose (RC)

- Hydrophilic membrane with very low adsorption for aqueous and organic-aqueous liquids, i.e. polar and medium polar sample solutions
- Binding capacity for proteins 84 µg per 25 mm filter



Ordering information - CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
				filters/pack	REF	filters/pack	REF
RC-20/25	0.20	25	labeled	100	729230	400	729230.400
RC-45/25	0.45	25	labeled	100	729231	400	729231.400

Ordering information - CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
			top	bottom	filters/pack	REF	filters/pack	REF
RC-20/15 MS	0.20	15	yellow	blue	100	729036	800	729036.800
RC-45/15 MS	0.45	15	colorless	blue	100	729037	800	729037.800
RC-20/25	0.20	25	yellow	blue	100	729030	400	729030.400
RC-45/25	0.45	25	colorless	blue	100	729031	400	729031.400

MS = minispikes on filter exit



CHROMAFIL® syringe filters

Polytetrafluoroethylene (PTFE)

- Hydrophobic membrane for nonpolar liquids and gases
- Very resistant towards all kinds of solvents as well as acids and bases
Flushing with alcohol, followed by water, makes the originally hydrophobic membrane more hydrophilic.



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
PTFE-20/25	0.20	25	labeled	100	729207	400	729207.400
PTFE-45/25	0.45	25	labeled	100	729205	400	729205.400
PTFE-100/25	1.0	25	labeled	100	729247	400	729247.400

Ordering information · CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code top	Color code bottom	Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
O-20/3	0.20	3	colorless	colorless	100	729014		
O-45/3	0.45	3	colorless	colorless	100	729015		
O-20/15 MS	0.20	15	yellow	colorless	100	729008	800	729008.800
O-45/15 MS	0.45	15	colorless	colorless	100	729009	800	729009.800
O-20/25	0.20	25	yellow	colorless	100	729007	400	729007.400

MS = minispikes on filter exit

Cellulose mixed esters (MV)

- Hydrophilic membrane with very low adsorption for aqueous or polar solutions



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
MV-20/25	0.20	25	labeled	100	729206	400	729206.400
MV-45/25	0.45	25	labeled	100	729204	400	729204.400

Ordering information · CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code top	Color code bottom	Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
A-20/25	0.20	25	yellow	yellow	100	729006	400	729006.400
A-45/25	0.45	25	colorless	yellow	100	729004	400	729004.400



Cellulose acetate (CA)

- Hydrophilic membrane for filtration of water-soluble oligomers and polymers, especially suited for biological macromolecules
- Very high shape stability in aqueous solutions
- Extremely low binding capacity for proteins (21 µg/filter)
- Also available in a sterile package (S) for filtration under sterile conditions (each filter individually sealed)



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	REF	BIG-BOX filters/pack	REF
CA-20/25	0.20	25	labeled	100	729226	400	729226.400
CA-45/25	0.45	25	labeled	100	729227	400	729227.400

Ordering information · CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code top	Color code bottom	Standard pack filters/pack	REF	BIG-BOX filters/pack	REF
CA-20/15 MS	0.20	15	yellow	red	100	729054	800	729054.800
CA-45/15 MS	0.45	15	colorless	red	100	729055	800	729055.800
CA-20/25	0.20	25	yellow	red	100	729026	400	729026.400
CA-45/25	0.45	25	colorless	red	100	729027	400	729027.400

Sterile filters

CA-20/15 MS (S)	0.20	15	yellow	red	50	729052		
CA-45/15 MS (S)	0.45	15	colorless	red	50	729053		
CA-20/25 (S)	0.20	25	yellow	red	50	729024		
CA-45/25 (S)	0.45	25	colorless	red	50	729025		

MS = minispikes on filter exit

Polyethersulfone (PES)

- Hydrophilic membrane for aqueous liquids and aqueous liquids with low organic contents
- Very low adsorption for pharmaceuticals and proteins good stability against acids and bases
- Binding capacity for proteins 29 µg per 25 mm filter



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	REF	BIG-BOX filters/pack	REF
PES-20/25	0.20	25	labeled	100	729240	400	729240.400
PES-45/25	0.45	25	labeled	100	729241	400	729241.400
PES-500/25	5.0	25	labeled	100	729242	400	729242.400



Syringe filters CHROMAFIL®

Polyamide (PA) = Nylon

- ⚡ Rather hydrophilic membrane for aqueous and organic-aqueous medium polar liquids



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
PA-20/25	0.20	25	labeled	100	729212	400	729212.400
PA-45/25	0.45	25	labeled	100	729213	400	729213.400

Ordering information · CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code top	Color code bottom	Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
AO-20/3	0.20	3	light beige	light beige	100	729010		
AO-45/3	0.45	3	light beige	light beige	100	729011		
AO-20/15 MS *	0.20	15	yellow	green	100	729048	800	729048.800
AO-45/15 MS *	0.45	15	colorless	green	100	729049	800	729049.800
AO-20/25	0.20	25	yellow	green	100	729012	400	729012.400
AO-45/25	0.45	25	colorless	green	100	729013	400	729013.400

Polyvinylidene difluoride (PVDF)

- ⚡ Hydrophilic membrane for polar and nonpolar solutions, water-soluble oligomers and polymers like proteins
- ⚡ Binding capacity for proteins 82 µg per 25 mm filter



Ordering information · CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
PVDF-20/25	0.20	25	labeled	100	729218	400	729218.400
PVDF-45/25	0.45	25	labeled	100	729219	400	729219.400

Ordering information · CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code top	Color code bottom	Standard pack filters/pack	Standard pack REF	BIG-BOX filters/pack	BIG-BOX REF
PVDF-20/15 MS	0.20	15	yellow	white	100	729043	800	729043.800
PVDF-45/15 MS	0.45	15	colorless	white	100	729044	800	729044.800

MS = minispikes on filter exit



Glass fiber (GF)

- Orange Inert filter, nominal pore size 1 µm, allows higher flow rates than small pore filters
- Orange For solutions with high loads of particulate matter or for highly viscous solutions (e.g., soil samples, fermentation broths)
- Orange As prefilters for other CHROMAFIL® filters, they prevent plugging of the membrane.



Ordering information - CHROMAFIL® Xtra

Type	Pore size [µm]	Membrane diameter [mm]		Standard pack filters/pack	REF	BIG-BOX filters/pack	REF
GF-100/25	nom. 1.0	25	labeled	100	729228	400	729228.400

Ordering information - CHROMAFIL®

Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack filters/pack	REF	BIG-BOX filters/pack	REF
			top	bottom				
GF-100/15 MS	nom. 1.0	15	blue	colorless	100	729034		
GF-100/25	nom. 1.0	25	yellow	black	100	729028	400	729028.400

MS = minispikes on filter exit



Sample clarification



CHROMAFIL® materials · compatibility

Chemical compatibility of filter materials

The following table lists the chemical compatibility of our CHROMAFIL® materials. The chemical compatibility depends on several parameters such as time, pressure, temperature and concentration. In most cases, CHROMAFIL® filters will have only short contact with a

solvent. In these cases they may be used despite of limited compatibility.

For example, a PTFE filter with PP housing does not liberate any UV-detectable substances during filtration of 5 mL THF, although PP shows only limited resistance towards THF.

Solvent	Material										PP
	MV	CA	RC	PA	PTFE	PVDF	PES	PET	GF		
Acetaldehyde	⊖	⊖	⊕	⦿	⊕	⊕		⊕	⊕	⦿	
Acetic acid, 100 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	
Acetone	⊖	⊖	⊕	⊕	⊕	⊖	⊖	⊕	⊕	⊕	
Acetonitrile	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Ammonia, 25 %	⊖	⊖	⦿	⊖	⊕	⊕	⊕	⦿	⊕	⊕	
Benzene	⊕	⊕	⊕	⊕	⊕	⦿		⊕	⊕	⦿	
n-Butanol	⊕	⊕	⊕	⦿	⊕	⊕	⊕	⊕	⊕	⊕	
Cyclohexane	⊕	⊕	⊕	⦿	⊕	⊕	⊕	⊕	⊕	⊕	
Dichloromethane	⊕	⊖	⊕	⊖	⊕	⊕	⊖	⊕	⊕	⊖	
Diethyl ether	⦿	⦿	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⦿	
Dimethylformamide	⊖	⊖	⦿	⊕	⊕	⊖	⊖	⊕	⊕	⊕	
1,4-Dioxane	⊖	⊖	⊕	⊕	⊕	⦿	⊖	⊕	⊕	⦿	
Ethanol	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Ethyl acetate	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⦿	
Ethylene glycol	⦿	⦿	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Formic acid, 100 %	⊕	⊖	⦿	⊖	⊕	⊕	⊕	⦿	⊕	⊕	
Hydrochloric acid, 30 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊖	⊕	⊕	
Methanol	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Nitric acid, 65 %	⊖	⊖	⊖	⊖	⦿	⦿		⦿	⊕	⊖	
Oxalic acid, 10 % aqueous	⊕	⊖	⊕	⊖	⊕	⊕		⊕	⊕	⊕	
Petroleum ether	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Phosphoric acid, 80 %	⊖	⊖	⦿	⊖	⊕	⦿		⊕	⊕	⊕	
Potassium hydroxide, 1 mol/L	⊖	⊖	⦿	⊕	⊕	⦿	⊕	⦿	⊕	⊕	
2-Propanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Sodium hydroxide, 1 mol/L	⊖	⊖	⦿	⊕	⊕	⦿	⦿	⦿	⦿	⊕	
Tetrachloromethane	⊕	⊖	⊕	⊕	⊕	⦿		⊕	⊕	⦿	
Tetrahydrofuran	⊖	⊖	⊕	⦿	⊕	⊕	⊖	⊕	⊕	⦿	
Toluene	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⦿	
Trichloroethene	⊕	⊕	⊕	⦿	⊕	⊕		⊕	⊕	⦿	
Trichloromethane (chloroform)	⊕	⊖	⊕	⊖	⊕	⊕	⊖	⊕	⊕	⊖	
Urea	⊕	⊕	⊕	⊕	⊕	⊕		⊕	⊕	⊕	
Water	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	
Xylene	⊕	⊕	⊕	⊕	⊕	⦿		⊕	⊕	⦿	

Data not guaranteed.

⊕ resistant, ⊖ not resistant, ⊙ limited resistance

MV = cellulose mixed esters, CA = cellulose acetate, RC = regenerated cellulose, PA = polyamide, PTFE = polytetrafluoroethylene, PVDF = polyvinylidene difluoride, PES = polyethersulfone, PET = polyester, GF = glass fiber, PP = polypropylene (housing material)



Hints for using CHROMAFIL® syringe filters

For optimum filtration results we recommend to keep the following in mind:

- ✦ Either discard the first mL or rinse the filter unit with 1 mL of the solvent prior to filtration
- ✦ Before filling the syringe, draw about 1 mL air into the syringe in order to minimize the liquid remaining in the filter
- ✦ Start filtration with a slight pressure; this will optimize the throughput of the filter. As soon as particles accumulate on the filter, filtration will become more difficult and the pressure on the filter will increase.
- ✦ Change the filter, whenever the resistance becomes too large in order to prevent rupture of the housing.
- ✦ Do not apply CHROMAFIL® syringe filters on humans; they are only intended for lab use!
- ✦ Always use syringes ≥ 10 mL; smaller syringes can easily cause pressures above the 6 bar limit of the filters.
- ✦ The temperature should not exceed 55 °C.
- ✦ Do not re-use the filters.

CHROMAFIL® filtration cartridges

- ✦ Filtration cartridges for sample clarification under vacuum (e.g., using the CHROMABOND® vacuum manifold or SPE automation systems like Gilson Aspec™, Rapidtrace®) or by gravity
- ✦ Cartridge sizes 3 mL and 6 mL
- ✦ Different membranes (PET, RC, PTFE, PVDF, GF) and pore sizes (0.2, 0.45 and 1.0 μm). Membrane materials correspond to the respective CHROMAFIL® syringe filters.



Sample clarification

Ordering information

Description	Pore size [μm]	Pack of [cartridges]	Column volume	
			3 mL	6 mL
Filtration cartridges PET (polyester)	0.20	100	730578.320	730578.620
Filtration cartridges PET (polyester)	0.45	100	730578.345	730578.645
Filtration cartridges RC (regenerated cellulose)	0.20	100	730068.320	730068.620
Filtration cartridges RC (regenerated cellulose)	0.45	100	730068.345	730068.645
Filtration cartridges PTFE (polytetrafluoroethylene)	0.20	100	730570.320	730570.620
Filtration cartridges PTFE (polytetrafluoroethylene)	0.45	100	730570.345	730570.645
Filtration cartridges PVDF (polyvinylidene difluoride)	0.20	100	730579.320	730579.620
Filtration cartridges PVDF (polyvinylidene difluoride)	0.45	100	730579.345	730579.645
Filtration cartridges GF (glass fiber)	nom. 1.0	100	730517.3100	730517.6100



96-well filter plates CHROMAFIL® MULTI 96

CHROMAFIL® MULTI 96 filter plates

- 96-well polypropylene plates for simultaneous filtration of 96 samples
- Advantages of this high-throughput system are:
Economical by saving time and solvent
Use of multi-channel pipetters facilitates liquid transfer steps
Readily adaptable to all common automated and robotic handling systems
Minimized dead volume ($\leq 40 \mu\text{L}$)
- Membrane materials correspond to the respective CHROMAFIL® syringe filters.



Sample clarification

Ordering information

Description	Pack of	REF
Filter plates with cellulose mixed ester filter elements (0.20 μm)	1	738770.M
Filter plates with cellulose mixed ester filter elements (0.45 μm)	1	738771.M
Filter plates with RC filter elements (regenerated cellulose, 0.2 μm)	1	738656.M
Filter plates with RC filter elements (regenerated cellulose, 0.45 μm)	1	738657.M
Filter plates with PTFE filter elements (0.2 μm)	1	738660.M
Filter plates with PTFE filter elements (0.45 μm)	1	738661.M
Filter plates with PTFE filter elements (1.0 μm)	1	738662.M
Filter plates with PTFE filter elements (3.0 μm)	1	738663.M
Filter plates with PE filter elements (20 μm)	1	738655.M
Filter plates with PE filter elements (50 μm)	1	738659.M
Filter plates with glass fiber filter elements (nominal 1 μm)	1	738655.2M
Filter plates with glass fiber filter elements (nominal 3 μm)	1	738658.M
CHROMABOND® MULTI 96 vacuum manifold for monoblocks, with reservoir tank, vacuum gauge, and control valve, for filtration with 96-well filter plates	1	738630.M

Disposable syringes with Luer tip (non-sterile, body and piston made from polypropylene)

Volume	Pack of	REF	
2 mL	100	729100	
5 mL	100	729101	
10 mL	100	729102	



Glass vials and inserts

According to the high requirements of chemical analyses, especially with regard to reproducibility combined with high detection sensitivity, the container material for the respective samples is of great importance. In chromatography generally vials made from 1st hydrolytic class glass are being used. This type includes borosilicate glasses like Duran®, Pyrex®, Fiolax®, and others. Glass of this class, often called neutral glass, has a very good chemical resistance towards acidic and neutral solutions. The relatively low alkali content permits good values for the resistance towards alkaline solutions, too. Except for the snap cap vials for storage of powdery samples the vials of our program are made from glass of 1st hydrolytic class (manufactured in accordance with Eu.Ph., U.S.P., DAB, Ph. Jap.).

The dimensions stated in this catalog with respect to vial diameter and height are exact values. Please note that other suppliers often list rounded values (e.g., 12 x 32 mm instead of 11.6 x 32 mm), the actual dimensions are, however, identical due to the required fit in the instrument. Our data concerning the volume are defined realistically usable volumes, not calculated values. For reasons of safety we state rather low values. Here, too, deviations from data of other suppliers may occur, which either use the calculated volume (e.g., 2 mL instead of 1.5 mL) or a defined, realistically usable volume in the upper range (e.g., 1.8 mL instead of 1.5 mL).

Septa guide

	Temperature resistance from / to	Analytical purity	Fragmentation due to hardness and molecular structure (coring)	Hardness (needle penetration)	Resealability (in case of multiple injections)
PTFE virginal	-200 °C / 260 °C	very high		very hard (but very thin material)	no resealability
Natural rubber / PTFE	-40 °C / 120 °C	low	high, big particles	very hard	high
Red Rubber / TEF (FEP)	-40 °C / 110 °C	medium	medium	medium hard	medium
Butyl	-40 °C / 120 °C	medium	medium	medium hard	medium
Butyl / PTFE	-40 °C / 120 °C	medium	medium	medium hard	medium
Silicone / PTFE	-60 °C / 200 °C	high	low to medium	soft	low to medium
PTFE / Silicone / PTFE	-60 °C / 200 °C	high	very low	soft	very low

In HPLC the septa often need to have a (cross-) slit as a penetration aid for thick and blunt needles. Furthermore, the slit avoids a vacuum in the vial during sample removal and thus guarantees constant volumes in case of multiple injections. The PTFE lamination protects the elastomeric carrier material of the septa from decomposition by aggressive samples or solvents, but also – in the other direction – the sample from possible contaminations through substances of the carrier material.

Autosampler compatibility charts (see page 103)

The autosampler compatibility charts generally show the most typical vials and closures for use on the instruments of a given manufacturer. In addition to the products listed in those charts, our catalog may contain other technically and functionally suitable products for use on a given autosampler which are not marketed actively as accessories by the respective manufacturer. We look forward to recommend any suitable product.

Compatibility charts have been compiled for the following instrument manufacturers: Agilent, CTC, Dionex, PerkinElmer, Shimadzu, Thermo Scientific, Varian (Agilent), VWR (Merck® / Hitachi), Waters®. Where applicable, each chart is divided into fields of use (GC, HPLC, Headspace).

We generally recommend that you ask for cost-free samples for testing purposes, as even technically comparable products may differ in their optical appearance.

We kindly ask for your understanding that we do not take over any guarantee for the correctness and completeness of the data indicated here.

Miscellaneous

Should you need more information concerning this product range, you can ask for our separate brochure "Vials and caps", which – among others – features 1:1 drawings of all glass products.

Except where explicitly mentioned, septa are assembled ready to use. Septa beneath or beside a cap are shown for illustration purposes only, and they are pictured upside down.

All drawings in this chapter are scale 1:2.

Should you wish to translate article numbers of other manufacturers into MACHERY-NAGEL REFs, we recommend to use our VialFinder on the internet under www.mn-net.com. Detailed information for its use are available as download.



Crimp neck vials and closures N 8



Crimp neck vials and closures N 8

Micro-vials with 8 mm crimp neck

0.2–1.2 mL usable volume

















Adapter required for use in an autosampler







Available with flat, round or conical bottom, made of clear and amber glass

Aluminium crimp closures N 8

Economic versions: three-layer septum Natural rubber / Butyl / TEF or two-layer septum Red Rubber / FEP

For more demanding analyses: high purity Silicone / PTFE septa

Description	Dimensions (scale 1:2)			Pack of	REF		
Crimp neck vials N 8							
							
							
Type of vial	Usable volume		OD x height				
Clear, conical	0.2 mL		5.5 x 31.5 mm		100	70286	
Clear, round bottom	0.3 mL		5.5 x 31.5 mm		100	70282	
Clear, flat bottom	0.8 mL		8.2 x 30 mm		100	70251	
Clear, conical	0.6 mL		7 x 40 mm		100	70212	
Amber, conical	0.6 mL		7 x 40 mm		100	70212.1	
Clear, flat bottom	0.7 mL		7 x 40 mm		100	702002	
Amber, flat bottom	0.7 mL		7 x 40 mm		100	702003	
Clear, flat bottom	1.2 mL		8.2 x 40 mm		100	702880	











Ready assembled crimp closures N 8 and plain crimp caps N 8						
						
Cap description	Septum description		Hardness	Thickness		
N 8 aluminium crimp cap, silver, center hole	PTFE virginal, white		53° shore D	0.25 mm	100	70283
N 8 aluminium crimp cap as above	Natural rubber / Butyl red-orange / TEF colorless		45° shore A	1.0 mm	100	70252.1
N 8 aluminium crimp cap as above	Red Rubber / FEP colorless		40° shore A	1.0 mm	100	702025
N 8 aluminium crimp cap as above	Silicone white / PTFE red		40° shore A	1.0 mm	100	70289
N 8 aluminium crimp cap as above	PTFE red / Silicone white / PTFE red		40° shore A	1.0 mm	100	702878
N 8 aluminium crimp cap as above	no liner		–	–	100	702800

Crimping tools N 8		
Manual crimper for 8 mm aluminium crimp caps	1	735126
Manual decapper for 8 mm aluminium crimp caps	1	735408











Screw neck vials and closures N 8

- Are among the oldest vials for HPLC and GC (besides crimp neck vials N 11)
- More and more replaced by screw neck vials N 9, which are easier to fill due to the wide opening compared to screw neck vials N 8 with small opening
- Due to the cap design not universally usable on all autosamplers in GC and HPLC – however, often used on instruments of VWR (Merck®) / Hitachi, Varian, Knauer, Gilson, Shimadzu and others
- In combination with closed top screw closures also used for sample storage (see page 92)

Description	Dimensions (scale 1:2)				Pack of	REF
Screw neck vials N 8, small opening (8–425 thread), and compatible inserts						
 70213  70213.2  702004  702893	 702968  702968.1  702974.1  702824  702005  702860					
Type of vial	Usable volume		OD x height			
Clear, flat bottom	1.5 mL		11.6 x 32 mm		100	70213
Amber, flat bottom	1.5 mL		11.6 x 32 mm		100	70213.2
Clear, flat bottom, label and scale	1.5 mL		11.6 x 32 mm		100	702004
Amber, flat bottom, label and scale	1.5 mL		11.6 x 32 mm		100	702893
Insert for small opening vials, clear, conical, 15 mm tip	0.1 mL		5 x 31 mm		100	702968*
Insert for small opening vials, clear, conical, 9 mm tip	0.15 mL		5 x 31 mm		100	702968.1*
Metal spring for conical inserts 5 x 31 mm	–		–		100	702974.1
Insert for small opening vials, clear, with plastic spring	0.1 mL		5 x 29 mm		100	702824
Insert for small opening vials, clear, flat bottom	0.25 mL		5 x 31 mm		100	702005
Micro-vial, clear, conical	1.1 mL		11.6 x 32 mm		100	702860
* Optionally you may use metal springs 702974.1 in combination with these products to push them up in the vial.						

Ready assembled screw closures N 8 and plain screw caps N 8

 702067  702068  70245  702066  702437  702069  70249  70250						
Cap description	Septa description		Hardness	Thickness		
N 8 PP screw cap, black, center hole as above, but with closed top	Red Rubber / FEP colorless		40° shore A	1.0 mm	100	702067
N 8 PP screw cap, black, center hole as above, but with closed top	Red Rubber / FEP colorless		40° shore A	1.0 mm	100	702068
N 8 PP screw cap, black, center hole as above, but with closed top	Silicone white / PTFE red		40° shore A	1.3 mm	100	70245
N 8 PP screw cap, black, center hole as above, but with closed top	Silicone white / PTFE red		40° shore A	1.3 mm	100	702066
N 8 PP screw cap, black, center hole as above, but with closed top	Silicone white / PTFE blue, slit		40° shore A	1.0 mm	100	702437
N 8 PP screw cap, black, center hole as above, but with closed top	PTFE red / Silicone white / PTFE red		40° shore A	1.0 mm	100	702069
N 8 PP screw cap, black, center hole as above, but with closed top	no liner		–	–	100	70249
N 8 PP screw cap, black, center hole as above, but with closed top	no liner		–	–	100	70250



Screw neck vials, inserts, and closures N 9

Description	Dimensions (scale 1:2)	Pack of	REF
N 8 Septa for screw caps N 8			
Material	Illustration	Hardness	Thickness
PTFE virginal, white		53° shore D	0.25 mm
Red Rubber / FEP colorless		40° shore A	1.0 mm
Silicone white / PTFE red		40° shore A	1.3 mm
Silicone white / PTFE blue, slit		40° shore A	1.0 mm






Screw neck vials and closures N 9


























- Can be used on almost all HPLC and GC autosamplers
- Large range of vials and closures
- Also available as bonded closures (advantage: thick (blunt) HPLC needles cannot push the septum into the vial)
- Also available as convenient Vial Kits with 100 vials and 100 caps and as pre-sealed vial-closure combinations

Description	Dimensions (scale 1:2)	Pack of	REF
Screw neck vials N 9, wide opening (short thread), and compatible inserts			
 702282 702293 702283 702284 silanized 702078 silanized 702079 silanized 702077 702813 702716 702818 702825 702006 702007 702008 702135 702088 702009 PP / glass PP	 32 11.6 32 11.6 32 11.6 32 11.6 32 11.6 32 11.6 32 11.6 31 6 31 6 29 5.7 31 6 32 11.6 32 11.6 32 11.6 32 11.6 32 11.6 32 11.6		
Type of vial	Usable volume	OD x height	
Clear, flat bottom	1.5 mL	11.6 x 32 mm	100 702282
Amber, flat bottom	1.5 mL	11.6 x 32 mm	100 702293
Clear, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702283
as above, silanized	1.5 mL	11.6 x 32 mm	100 702078
Amber, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702284
as above, silanized	1.5 mL	11.6 x 32 mm	100 702079
Insert for wide opening vials, clear, conical, 15 mm tip	0.2 mL	6 x 31 mm	100 702813
as above, silanized	0.2 mL	6 x 31 mm	100 702077
Insert for wide opening vials, clear, conical, 12 mm tip	0.25 mL	6 x 31 mm	100 702716
Insert for wide opening vials, clear, with plastic spring	0.1 mL	5.7 x 29 mm	100 702818
Insert for wide opening vials, clear, flat bottom	0.3 mL	6 x 31 mm	100 702825
Micro-vial, clear, 15 µL funnel in solid glass bottom	1.1 mL	11.6 x 32 mm	100 702006
Micro-vial, clear, with integrated 0.2 mL insert	0.2 mL	11.6 x 32 mm	100 702007
Micro-vial, amber, with integrated 0.2 mL insert	0.2 mL	11.6 x 32 mm	100 702008
Micro-vial, polypropylene, transparent, with integrated 0.15 mL glass insert, conical	0.15 mL	11.6 x 32 mm	100 702135
Micro-vial, clear, conical, with a round pedestal glass plate	1.1 mL	11.6 x 32 mm	100 702088
Micro-vial, polypropylene, transparent, with inner cone	0.3 mL	11.6 x 32 mm	100 702009

Screw neck vials, inserts, and closures N 9






Description	Dimensions (scale 1:2)		Pack of	REF
Bonded screw closures N 9 (septa firmly connected with the cap, cannot be removed)				
 702028	 702026	 702027		
Cap description	Septa description	Hardness	Thickness	
N 9 PP bonded screw cap, blue, center hole	Red Rubber / TEF colorless	65° shore A	1.0 mm	100 702028
N 9 PP bonded screw cap, blue, center hole	Silicone beige / PTFE white	45° shore A	1.3 mm	100 702026
N 9 PP bonded screw cap, blue, center hole	Silicone beige / PTFE white, slit	45° shore A	1.3 mm	100 702027

Ready assembled screw closures N 9						
 702029	 702031	 702032				
Cap description	Septa description	Hardness	Thickness			
N 9 PP screw cap, transparent, center hole	PTFE virginal, white	53° shore D	0.25 mm	100		702029
N 9 PP screw cap, blue, center hole	PTFE virginal, white	53° shore D	0.25 mm	100		702031
N 9 PP screw cap blue, closed top	PTFE virginal, white	53° shore D	0.25 mm	100		702032
 702030	 702732	 702080	 702081	 702082	 702033	
N 9 PP screw cap, transparent, center hole	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702030
N 9 PP screw cap, blue, center hole	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702732
N 9 PP screw cap, black, center hole	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702080
N 9 PP screw cap, red, center hole	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702081
N 9 PP screw cap, green, center hole	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702082
N 9 PP screw cap blue, closed top	Red Rubber / FEP colorless	40° shore A	1.0 mm	100		702033
 702287	 702287.1	 702036	 702037	 702038	 702034	
 702107						
N 9 PP screw cap, transparent, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702287
N 9 PP screw cap, blue, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702287.1
N 9 PP screw cap, black, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702036
N 9 PP screw cap, red, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702037
N 9 PP screw cap, green, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702038
N 9 PP screw cap, yellow, center hole	Silicone white / PTFE red	40° shore A	1.0 mm	100		702107
N 9 PP screw cap blue, closed top	Silicone white / PTFE red	40° shore A	1.0 mm	100		702034
 702288	 702288.1	 702039	 702040	 702083	 702109	
N 9 PP screw cap, transparent, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702288
N 9 PP screw cap, blue, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702288.1
N 9 PP screw cap, black, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702039
N 9 PP screw cap, red, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702040
N 9 PP screw cap, green, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702083
N 9 PP screw cap, yellow, center hole	Silicone white / PTFE blue, slit	40° shore A	1.0 mm	100		702109
 702286	 702035	 702084	 702085			
N 9 PP screw cap, transparent, center hole	PTFE red / Silicone white / PTFE red	40° shore A	1.0 mm	100		702286
N 9 PP screw cap, blue, center hole	as above	40° shore A	1.0 mm	100		702035
N 9 PP screw cap, red, center hole	as above	40° shore A	1.0 mm	100		702084
N 9 PP screw cap, green, center hole	as above	40° shore A	1.0 mm	100		702085



Screw neck vials, inserts, and closures N 9

Description	Dimensions (scale 1:2)			Pack of	REF
N 9 septa for screw caps N 9					
Material	Illustration	Hardness	Thickness		
PTFE virginal, white		53° shore D	0.25 mm	100	702043
Red Rubber / FEP colorless		40° shore A	1.0 mm	100	702041
Silicone white / PTFE red		40° shore A	1.0 mm	100	702042



Vial Kits screw neck N 9

Packs of 100 vials and 100 closures, each

Closure →					
	702287.1	702288.1	702732	702026	702027
Vial ↓					
702282: 1.5 mL, clear, flat bottom	702201	702204	702207		
702283: 1.5 mL, clear, flat bottom, label and scale	702202	702205	702208	702211	702213
702284: 1.5 mL, amber, flat bottom, label and scale	702203	702206	702209	702212	702214
702009: 0.3 mL, PP, transparent, with inner cone		702226			

Other Vial Kits on request

Pre-sealed vial-closure combinations with screw neck N 9















Vial description	Closure description	Pack of	REF
Pre-sealed vials 702282: 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening	pre-screwed with 702732: N 9 PP screw cap, blue, center hole, Red Rubber / FEP colorless, 40° shore A, 1.0 mm	100	702857
Pre-sealed vials 702283: 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening, label and scale	pre-screwed with 702732: N 9 PP screw cap, blue, center hole, Red Rubber / FEP colorless, 40° shore A, 1.0 mm	100	702858
Pre-sealed vials 702282: 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening	pre-screwed with 702287.1: N 9 PP screw cap, blue, center hole, Silicone white / PTFE red, 40° shore A, 1.0 mm	100	702874
Pre-sealed vials 702283: 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening, label and scale	pre-screwed with 702288.1: N 9 PP screw cap, blue, center hole, Silicone white / PTFE blue, slit, 40° shore A, 1.0 mm	100	702863
Pre-sealed vials 702284: 1.5 mL screw neck vial N 9, 11.6 x 32 mm, amber, flat bottom, wide opening, label and scale	pre-screwed with 702288.1: N 9 PP screw cap, blue, center hole, Silicone white / PTFE blue, slit, 40° shore A, 1.0 mm	100	702873

Other pre-sealed vial-closure combinations on request









Screw neck vials and closures N 10

- Wide opening for easy filling
- Due to the cap height not universally suitable for all instruments
- Often used on Waters®, Shimadzu, PerkinElmer and Jasco instruments

Description	Dimensions (scale 1:2)			Pack of	REF	
Screw neck vials N 10, wide opening (10–425 thread), and compatible inserts						
						
702011	702012	702013	702813 silanized 702077	702716	702818	702825
						
Type of vial	Usable volume		OD x height			
Clear, flat bottom	1.5 mL		11.6 x 32 mm			
Clear, flat bottom, label and scale	1.5 mL		11.6 x 32 mm			
Amber, flat bottom, label and scale	1.5 mL		11.6 x 32 mm			
Insert for wide opening vials, clear, conical, 15 mm tip as above,	0.2 mL		6 x 31 mm			
	0.2 mL		6 x 31 mm			
Insert for wide opening vials, clear, conical, 12 mm tip	0.25 mL		6 x 31 mm			
Insert for wide opening vials, clear, with plastic spring	0.1 mL		5.7 x 29 mm			
Insert for wide opening vials, clear, flat bottom	0.3 mL		6 x 31 mm			

Screw closures N 10 and plain screw caps N 10

					
702044	702045	702046	702047	702048	702049
Cap description		Septa description		Hardness	Thickness
N 10 PP bonded screw cap*, black, center hole		Red Rubber / TEF colorless		65° shore A	1.0 mm
N 10 PP bonded screw cap* as above		Silicone white / PTFE beige		45° shore A	1.5 mm
N 10 PP bonded screw cap* as above		Silicone white / PTFE red		45° shore A	1.5 mm
N 10 PP bonded screw cap* as above		Silicone white / PTFE blue, slit		45° shore A	1.5 mm
* Septum firmly connected with the cap, cannot be removed					
N 10 PP screw cap, black, center hole		PTFE red / Silicone white / PTFE red		45° shore A	1.0 mm
N 10 PP screw cap, black, center hole		no liner		–	–
















Crimp neck vials, inserts, and closures N 11



Crimp neck vials and closures N 11








- Broad variety of standard crimp neck vials (with small or wide opening), as well as crimp neck micro-vials for smaller sample volumes
- Economic closures: Natural rubber / TEF (2 layers), Natural rubber / Butyl / TEF (3 layers) and Red Rubber / FEP (2 layers)
- For more demanding analyses: analytically pure Silicone / PTFE septa with lower fragmentation
- Magnetic closure: REF 702879 for use on CTC GC PAL
















Description	Dimensions (scale 1:2)		Pack of	REF	
Crimp neck vials N 11, small opening, and compatible inserts					
 70201CG	 70214CG	 702968	 702968.1	 702824	 702005
Type of vial	Usable volume		OD x height		
Clear, flat bottom, small opening	1.5 mL		11.6 x 32 mm		
Amber, flat bottom, small opening	1.5 mL		11.6 x 32 mm		
Insert for small opening vials, clear, conical, 15 mm tip	0.1 mL		5 x 31 mm		
Insert for small opening vials, clear, conical, 9 mm tip	0.15 mL		5 x 31 mm		
Insert for small opening vials, clear, with plastic spring	0.1 mL		5 x 29 mm		
Insert for small opening vials, clear, flat bottom	0.25 mL		5 x 31 mm		
* Optionally you may use metal springs 702974.1 in combination with these products to push them up in the vial.					





Crimp neck vials N 11, wide opening, and compatible inserts						
 70201HP	 702885 silanized 702075	 702892 silanized 702076	 702813 silanized 702077	 702716	 702818	 702825
Type of vial	Usable volume		OD x height			
Clear, flat bottom, wide opening	1.5 mL		11.6 x 32 mm		100	70201HP
Clear, flat bottom, wide opening, label and scale	1.5 mL		11.6 x 32 mm		100	702885
as above, silanized	1.5 mL		11.6 x 32 mm		100	702075
Amber, flat bottom, wide opening, label and scale	1.5 mL		11.6 x 32 mm		100	702892
as above, silanized	1.5 mL		11.6 x 32 mm		100	702076
Insert for wide opening vials, clear, conical, 15 mm tip	0.2 mL		6 x 31 mm		100	702813
as above, silanized	0.2 mL		6 x 31 mm		100	702077
Insert for wide opening vials, clear, conical, 12 mm tip	0.25 mL		6 x 31 mm		100	702716
Insert for wide opening vials, clear, with plastic spring	0.1 mL		5.7 x 29 mm		100	702818
Insert for wide opening vials, clear, flat bottom	0.3 mL		6 x 31 mm		100	702825

Crimp neck vials, inserts, and closures N 11



Description	Dimensions (scale 1:2)				Pack of	REF
Crimp neck micro-vials N 11						
						
702888	702891	702014	702134 PP / glass	702015	702016	702141
Type of vial	Usable volume		OD x height			
Micro-vial, clear, flat bottom 15 µL funnel in solid glass bottom	1.1 mL		11.6 x 32 mm		100	702888
Micro-vial, clear, with integrated 0.2 mL insert	0.2 mL		11.6 x 32 mm		100	702891
Micro-vial, amber, with integrated 0.2 mL insert	0.2 mL		11.6 x 32 mm		100	702014
Micro-vial, polypropylene, transparent, with integrated 0.15 mL glass insert, conical	0.15 mL		11.6 x 32 mm		100	702134
Micro-vial, clear, conical with round pedestal glass plate	1.1 mL		11.6 x 32 mm		100	702015
Micro-vial, amber, conical with round pedestal glass plate	1.1 mL		11.6 x 32 mm		100	702016
Micro-vial, clear, conical	1.1 mL		11.6 x 32 mm		100	702141

Ready assembled aluminium crimp closures N 11														
														
70284	70231	702001	702730	70256	70231.1	70231.2	70231.3	70231.4	70288	702823	702995	702879	702801	702401
Cap description				Septa description				Hardness		Thickness				
N 11 aluminium crimp cap, silver, center hole				PTFE virginal, white				53° shore D		0.25 mm		100	70284	
N 11 crimp cap as above, silver				Natural rubber / Butyl red-orange / TEF colorless				45° shore A		1.3 mm		100	70231	
N 11 crimp cap as above, silver				Natural rubber red-orange / TEF colorless (corresponds to Agilent quality)				60° shore A		1.0 mm		100	702001	
N 11 crimp cap as above, silver				Red Rubber / FEP colorless				40° shore A		1.0 mm		100	702730	
N 11 crimp cap as above, silver				Natural rubber / Butyl red-orange / TEF colorless				45° shore A		1.0 mm		100	70256	
N 11 crimp cap as above, green				as above				45° shore A		1.0 mm		100	70231.1	
N 11 crimp cap as above, red				as above				45° shore A		1.0 mm		100	70231.2	
N 11 crimp cap as above, blue				as above				45° shore A		1.0 mm		100	70231.3	
N 11 crimp cap as above, gold				as above				45° shore A		1.0 mm		100	70231.4	
N 11 crimp cap as above, silver				Silicone white / PTFE red				40° shore A		1.3 mm		100	70288	
N 11 crimp cap as above, silver				Silicone white / PTFE blue, cross-slit				40° shore A		1.5 mm		100	702823	
N 11 crimp cap as above, silver				PTFE red / Silicone white / PTFE red				40° shore A		1.0 mm		100	702995	
N 11 magnetic crimp cap, gold, center hole				Silicone white / PTFE red				55° shore A		1.0 mm		100	702879	
N 11 aluminium crimp cap, silver, center hole				no liner				-		-		100	702801	
N 11 PE cap, transparent, closed top, with thin piercing area												100	702401	





N 11 Septa for crimp caps N 11					
Material	Illustration	Hardness	Thickness		
PTFE virginal, white		53° shore D	0.25 mm	100	70262
Red Rubber / FEP colorless		40° shore A	1.0 mm	100	702065
Silicone white / PTFE red		40° shore A	1.3 mm	100	70263
PTFE red / Silicone white / PTFE red		40° shore A	1.0 mm	100	70264



Crimp neck vials, inserts, and closures N 11

Vial Kits crimp neck N 11

Packs of 100 vials and 100 closures, each

Closure →				
Vial ↓	70288	702995	70256	
70201HP: 1.5 mL, clear, flat bottom	702215	702218	702222	
702885: 1.5 mL, clear, flat bottom, label and scale	702216	702219	702223	
702892: 1.5 mL, amber, flat bottom, label and scale	702217	702221	702224	

Other Vial Kits on request

Pre-sealed vial-closure combinations with crimp neck N 11

Vial description	Closure description	Pack of	REF
Pre-sealed vials 70201CG: 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, small opening	crimped with 70256: N 11 aluminium crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.0 mm	100	702881
Pre-sealed vials 70201HP: 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, wide opening	crimped with 70256: N 11 aluminium crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.0 mm	100	702101HP
Pre-sealed vials 702892: 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, amber, flat bottom, wide opening, label and scale	crimped with 70256: N 11 aluminium crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.0 mm	100	702859

Other pre-sealed vial-closure combinations on request

Crimping tools N 11


















Description	Pack of	REF
Manual ergonomic crimper for 11 mm aluminium crimp caps	1	735211
Manual ergonomic decapper for 11 mm aluminium crimp caps	1	735311
Manual crimper, height adjustable, for 11 mm aluminium crimp caps	1	735111
Manual decapper for 11 mm aluminium crimp caps	1	735911
Pneumatic crimping tool for 11 mm aluminium crimp caps (complete with hand switch)	1	735114
Pneumatic crimping tool for 11 mm aluminium crimp caps (complete with foot switch)	1	735117
Crimping head N 11 (without pneumatic basic tool)	1	735121
Decapping head N 11 (without pneumatic basic tool)	1	735434.1
Pneumatic basic tool with hand switch	1	735124
Pneumatic basic tool with foot switch	1	735125

For more crimping tools N 11 see pages 100 and 101












Snap ring vials and closures N 11

- Quick, convenient sealing method which, however, should only be used in HPLC
- Can be used on all common HPLC autosamplers
- Alternatively crimp closures N 11 can be used (see preceding pages).
- 0.3 mL PP snap ring vial for special applications, e.g., in ion chromatography
- Most frequent closure: with cross-slit Silicone / PTFE septum, which supports easy penetration with the relatively thick, blunt HPLC needle

Description	Dimensions (scale 1:2)				Pack of	REF
Snap ring vials N 11, wide opening, and compatible inserts						
 702714	 702713	 702712	 702813 silanized 702077	 702716	 702818	 702825
 702709	 702134 PP / glass	 702809 PP				
 32 11.6	 32 11.6	 32 11.6	 31 6	 31 6	 29 5.7* 6	 31 6
Type of vial			Usable volume		OD x height	
Clear, flat bottom			1.5 mL		11.6 x 32 mm	
Clear, flat bottom, label and scale			1.5 mL		11.6 x 32 mm	
Amber, flat bottom, label and scale			1.5 mL		11.6 x 32 mm	
Insert for wide opening vials, clear, conical, 15 mm tip as above, silanized			0.2 mL		6 x 31 mm	
			0.2 mL		6 x 31 mm	
Insert for wide opening vials, clear, conical, 12 mm tip			0.25 mL		6 x 31 mm	
Insert for wide opening vials, clear, with plastic spring			0.1 mL		5.7 x 29 mm	
Insert for wide opening vials, clear, flat bottom			0.3 mL		6 x 31 mm	
Micro-vial, clear, with integrated 0.2 mL insert			0.2 mL		11.6 x 32 mm	
Micro-vial, polypropylene, transparent, with integrated 0.15 mL glass insert, conical			0.15 mL		11.6 x 32 mm	
Micro-vial, polypropylene, transparent, with conical insert			0.3 mL		11.6 x 32 mm	







Snap ring vials, inserts, and closures N 11

Description		Dimensions (scale 1:2)		Pack of		REF		
Ready assembled snap ring closures N 11								
 702731	 702063	 702710	 702710.1	 702064	 702717.2	 702718	 702718.1	 702401
Cap description		Septa description		Hardness		Thickness		
N 11 PE snap ring cap, blue, center hole		Red Rubber / TEF colorless		65° shore A		1.0 mm		
as above, transparent cap		Red Rubber / TEF colorless		65° shore A		1.0 mm		
N 11 PE snap ring cap, blue, center hole		Silicone white / PTFE red		55° shore A		1.0 mm		
as above, transparent cap		Silicone white / PTFE red		55° shore A		1.0 mm		
N 11 PE snap ring cap, blue, center hole		Silicone white / PTFE blue, cross-slit		55° shore A		1.0 mm		
as above, transparent cap		Silicone white / PTFE blue, cross-slit		55° shore A		1.0 mm		
N 11 PE snap ring cap, blue, center hole		PTFE red / Silicone white / PTFE red		45° shore A		1.0 mm		
as above, transparent cap		PTFE red / Silicone white / PTFE red		45° shore A		1.0 mm		
N 11 PE cap, transparent, closed top, with thin piercing area						100		
						702401		



Vial Kits snap ring N 11

Packs of 100 vials and 100 closures, each

Closure →					
		 702710	 702064	 702731	 702718
Vial ↓					
702714: 1.5 mL, clear, flat bottom		702225	702228	702232	702235
702713: 1.5 mL, clear, flat bottom, label and scale		702719	702229	702233	702236
702712: 1.5 mL, amber, flat bottom, label and scale		702227	702231	702234	702237

Other Vial Kits on request

Containers for screw neck vials N 8 / N 9 / N 10 as well as crimp neck and snap ring vials N 11











Description	Pack of	REF
81 position container blue for vials 11.6 x 32 mm, outer length 130 mm, outer width 130 mm, outer height 45 mm, coded, with transparent lid (suitable for freezers)	1	702514

Crimp neck vials and closures N 13



Crimp neck vials and closures N 13

Description	Dimensions (scale 1:2)		Pack of	REF
Crimp neck vials N 13				
 70255	 70203	 70258		
Type of vial	Usable volume	OD x height		
Clear, conical	1 mL	11 x 40 mm	100	70255
Clear, flat bottom	2 mL	13.75 x 35 mm	100	70203
Clear, flat bottom	2 mL	11 x 43 mm	100	70258
Ready assembled crimp closures N 13 and plain crimp caps N 13				
 70257	 70232	 702802	 702803	 702820
Cap description	Septa description	Hardness	Thickness	
N 13 aluminium crimp cap, silver, center hole	Butyl dark gray / PTFE gray (only centrally laminated, typically called Pharma-Fix)	50° shore A	2 mm	100 70257
N 13 aluminium center tear off cap, gold	Butyl dark gray / PTFE gray (only centrally laminated, typically called Pharma-Fix)	50° shore A	2 mm	100 70232
N 13 aluminium crimp cap, silver, center hole	no liner	-	-	100 702802
N 13 aluminium center tear off cap, coppery	no liner	-	-	100 702803
Stoppers N 13				
	N 13 Bromobutyl stopper, gray	45° shore A	-	100 702820

Crimping tools N 13

Manual crimper, height adjustable, for 13 mm aluminium crimp caps	1	735113
Manual crimper, height adjustable, for 13 mm Flip Top/Flip Off caps	1	735133
Manual decapper for 13 mm aluminium crimp caps	1	735913

Containers for crimp and screw neck vials N 13



Description	Pack of	REF
49 position container blue for crimp and screw neck vials N 13, outer length 130 mm, outer width 130 mm, outer height 50 mm, with transparent lid (suitable for freezers)	1	702515















Screw neck vials, inserts, and closures N 13






Screw neck vials and closures N 13

- Generally used for large volume samples in HPLC
- In combination with closed top screw closures suitable for sample storage (see page 92)
- Compatible insert requires metal spring for central alignment
- Range of ready assembled closures and plain caps with center hole or with closed top as well as separate septa (PTFE virginal, Red Rubber / FEP and Silicone / PTFE) are available.

Description	Dimensions (scale 1:2)		Pack of	REF
Screw neck vials N 13 (13–425 thread) and compatible insert				
 702962	 702973	 702089	 702972	 702974
Type of vial	Usable volume		OD x height	
Clear, flat bottom	4 mL		14.75 x 45 mm	
Amber, flat bottom	4 mL		14.75 x 45 mm	
Amber, flat bottom, label and scale	4 mL		14.75 x 45 mm	
Insert, clear, conical, metal spring required	0.3 mL		6 x 40 mm	
Metal spring for 702972	–		–	
			100	702974

Ready assembled screw closures and plain screw caps N 13						
 702103	 702050	 702051	 702926	 702052	 702963	 702966
Cap description	Septa description		Hardness	Thickness		
N 13 screw cap (13-425), green, closed top	F217 white / PTFE beige (firmly fixed)			1.5 mm	100	702103
N 13 PP screw cap, black, center hole	Red Rubber / FEP colorless		40° shore A	1.5 mm	100	702050
as above, but closed top	Red Rubber / FEP colorless		40° shore A	1.5 mm	100	702051
N 13 PP screw cap, black, center hole	Silicone white / PTFE red		40° shore A	1.3 mm	100	702926
as above, but closed top	Silicone white / PTFE red		40° shore A	1.3 mm	100	702052
N 13 PP screw cap, black, center hole	no liner		-	-	100	702963
as above, but closed top	no liner		-	-	100	702966

N 12 septa for screw caps N 13						
Material	Illustration	Hardness	Thickness			
PTFE virginal, white		53° shore D	0.25 mm	100	70260	
Red Rubber / FEP colorless		40° shore A	1.5 mm	100	702053	
Silicone white / PTFE red		40° shore A	1.3 mm	100	702292	



Micro reaction vials and closures

Description	Dimensions (scale 1:2)	Pack of	REF
Micro reaction vials with screw neck N 13 / N 20 and closures			
N 13		N 20	
<p>702210</p> <p>702220</p>		<p>702230</p> <p>702240</p>	

Micro reaction vials complete with screw caps and septa

Type of vial	Usable volume	OD x height		
Screw neck vials N 13 (13–425 thread), clear, with inner conical funnel in solid, flat glass bottom; complete with screw closure	0.25 mL	14 x 33 mm	1	702210
Screw neck vials N 13 (13–425 thread), clear, with inner conical funnel in solid, flat glass bottom; complete with screw closure	0.75 mL	14 x 46 mm	1	702220
Screw neck vials N 20 (20–400 thread), clear, with inner conical funnel in solid, flat glass bottom; complete with screw closure	3 mL	20 x 46 mm	1	702230
Screw neck vials N 20 (20–400 thread), clear, with inner conical funnel in solid, flat glass bottom; complete with screw closure	4.5 mL	20 x 60 mm	1	702240

Replacement screw caps with septa for micro reaction vials N 13 and N 20

Cap description	Septa description	Hardness	Thickness		
N 13 PP screw cap, black, center hole, assembled	Silicone white / PTFE red	40° shore A	1.3 mm	100	702926
N 20 phenolic screw cap, black, center hole (unassembled)	Butyl red / PTFE gray	55° shore A	1.4 mm	48	702280

Replacement septa N 12 and N 18 for screw caps N 13 and N 20, respectively

Material	Illustration	Hardness	Thickness		
Silicone white / PTFE red		40° shore A	1.3 mm	100	702292
Butyl red / PTFE gray		55° shore A	1.4 mm	48	702300



Special vials and closures

Screw neck vials for storage of liquid samples

- Usable volumes of 1.5 up to 24 mL
- Available neck sizes N 8, N 9, N 13, N 15, N 18 and N 20
- Corresponding closed top screw closures with different septa materials

N 8				N 9				N 13		
70213	70213.2	702004	702893	702282	702293	702283	702284	702962	702973	702089
can be used with				can be used with				can be used with		
702068		702066		702032	702033	702034		702103	702051	702052

Description	Dimensions	Pack of	REF
Screw neck vials N 8, small opening (8–425 thread)			
Type of vial	Usable volume	OD x height	
Clear, flat bottom	1.5 mL	11.6 x 32 mm	100 70213
Amber, flat bottom	1.5 mL	11.6 x 32 mm	100 70213.2
Clear, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702004
Amber, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702893
Closed top screw closures N 8			
Cap description	Septa description	Hardness Thickness	
N 8 PP screw cap, black, closed top	Red Rubber / FEP colorless	40° shore A 1.0 mm	100 702068
N 8 PP screw cap, black, closed top	Silicone white / PTFE red	40° shore A 1.3 mm	100 702066
For drawings of vials see page 79			

Description	Dimensions	Pack of	REF
Screw neck vials N 9, wide opening (short thread)			
Type of vial	Usable volume	OD x height	
Clear, flat bottom	1.5 mL	11.6 x 32 mm	100 702282
Amber, flat bottom	1.5 mL	11.6 x 32 mm	100 702293
Clear, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702283
as above, silanized	1.5 mL	11.6 x 32 mm	100 702078
Amber, flat bottom, label and scale	1.5 mL	11.6 x 32 mm	100 702284
as above, silanized	1.5 mL	11.6 x 32 mm	100 702079
Closed top screw closures N 9			
Cap description	Septa description	Hardness Thickness	
N 9 PP screw cap blue, closed top	PTFE virginal, white	53° shore D 0.25 mm	100 702032
N 9 PP screw cap blue, closed top	Red Rubber / FEP colorless	40° shore A 1.0 mm	100 702033
N 9 PP screw cap blue, closed top	Silicone white / PTFE red	40° shore A 1.0 mm	100 702034
For drawings of vials see page 80			



Description	Dimensions	Pack of	REF
Screw neck vials N 13 (13–425 thread)			
Type of vial	Usable volume	OD x height	
Clear, flat bottom	4 mL	14.75 x 45 mm	100 702962
Amber, flat bottom	4 mL	14.75 x 45 mm	100 702973
Amber, flat bottom, label and scale	4 mL	14.75 x 45 mm	100 702089
Closed top screw closures N 13			
Cap description	Septa description	Hardness	Thickness
N 13 screw cap (13–425), green, closed top	F217 white / PTFE beige (firmly fixed)		1.5 mm 100 702103
N 13 PP screw cap, black, closed top	Red Rubber / FEP colorless	40° shore A	1.5 mm 100 702051
N 13 PP screw cap, black, closed top	Silicone white / PTFE red	40° shore A	1.3 mm 100 702052
For drawings of vials see page 90			

Description	Dimensions (scale 1:2)	Pack of	REF
Screw neck vials N 15, N 18, and N 20 for storage of liquid samples			
<div> <div> <p>N 15</p> <p>702096 / 702311</p> <p>70285 / 702097</p> <p>702104 702104</p> </div> <div> <p>N 18</p> <p>702098</p> <p>702105</p> </div> <div> <p>N 20</p> <p>702099</p> <p>702106</p> </div> </div>			
Type of vial	Usable volume	OD x height	
Screw neck vial N 15 (15–425 thread), clear, flat bottom	8 mL	16.6 x 61 mm	100 702096
Screw neck vial N 15 (15–425 thread), amber, flat bottom	8 mL	16.6 x 61 mm	100 702311
Screw neck vial N 15 (15–425 thread), clear, flat bottom	12 mL	18.5 x 66 mm	100 70285
Screw neck vial N 15 (15–425 thread), amber, flat bottom	12 mL	18.5 x 66 mm	100 702097
Screw neck vial N 18 (18–400 thread), clear, flat bottom	16 mL	20.6 x 71 mm	100 702098
Screw neck vial N 20 (20–400 thread), clear, flat bottom	24 mL	22.7 x 86 mm	100 702099
Closed top screw closures N 15, N 18, and N 20			
Cap description	Septa description	Thickness	
N 15 screw cap (15–425), green, closed top	F217 white / PTFE beige (firmly fixed)	1.5 mm	100 702104
N 18 screw cap (18–400), green, closed top	F217 white / PTFE beige (firmly fixed)	1.5 mm	100 702105
N 20 screw cap (20–400), green, closed top	F217 white / PTFE beige (firmly fixed)	1.5 mm	100 702106

For screw neck vials with even larger volumes please see page 102.



Special vials and closures

Snap cap vials for storage of powdery samples

- Available sizes N 18 and N 22
- Usable volumes from 5 up to 25 mL
- Glass of 3rd hydrolytic class

Description	Dimensions (scale 1:2)	Pack of	REF
Snap cap vials N 18 and N 22 for storage of powdery samples			
 70271	 70272	 70274	for 70271 and 70272:
 702019	 70273	 70275	for 702019 and 70273:
Snap cap vials			
Type of vial	Usable volume	OD x height	
N 18, clear, flat bottom	5 mL	20 x 40 mm	100 70271
N 18, clear, flat bottom	10 mL	22 x 50 mm	100 70272
N 22, clear, flat bottom	15 mL	26 x 48 mm	100 702019
N 22, clear, flat bottom	25 mL	26 x 65 mm	100 70273
PE snap caps			
N 18 PE snap cap, transparent, for 70271 and 70272			100 70274
N 22 PE snap cap, transparent, for 702019 and 70273			100 70275

Shell vials N 8 and N 12

- Economic combination of vials and closures for uncritical HPLC applications






Description	Dimensions (scale 1:2)	Pack of	REF
Shell vials N 8 and N 12 with PE plug			
 70202.1	 702017	 702807	for 70202.1 and 702017:
 702018	 702054	 702054	for 702018:
Shell vials			
Type of vial	Usable volume	OD x height	
N 8, clear, flat bottom	1 mL	8.2 x 40 mm	100 70202.1
N 8, amber, flat bottom	1 mL	8.2 x 40 mm	100 702017
N 12, clear, flat bottom	2 mL	11.6 x 31.5 mm	100 702018
PE plugs			
N 8 PE plug, transparent, for 70202.1 and 702017			100 702807
N 12 PE plug, transparent, for 702018			100 702054



Screw neck vials and magnetic screw closures N 18



- Headspace vials for convenient, safe and consistent handling
- High tightness and better reproducibility of the sealing process (as compared to crimping)
- Thinner septum (1.5 mm instead of 3 mm septum thickness in crimp caps), thus safe penetration of the needle and less fragmentation (especially important for SPME applications)
- Improved run in autosamplers with magnets (CTC Combi PAL and equivalent instruments), since a flat surface for the magnet is ensured, thus avoiding that the filled vial can drop from the magnet

Description	Dimensions (scale 1:2)		Pack of	REF
Headspace screw neck vials N 18				
 702866	 702826	 702827	 702055	 702072
Type of vial	Usable volume	OD x height		
Clear, rounded bottom	10 mL	22.5 x 46 mm	100	702866
Clear, rounded bottom	20 mL	22.5 x 75.5 mm	100	702826

Ready assembled, magnetic screw closures N 18					
Cap description	Septa description	Hardness	Thickness		
N 18 magnetic screw cap, silver, center hole	Silicone blue transparent / PTFE white	45° shore A	1.5 mm	100	702827
N 18 magnetic screw cap, silver, center hole	Silicone white / PTFE blue	55° shore A	1.5 mm	100	702055
N 18 magnetic screw cap, silver, center hole	Red Rubber / TEF colorless	65° shore A	1.5 mm	100	702072

N 17 septa for magnetic screw caps N 18					
Material	Illustration	Hardness	Thickness		
Silicone blue transparent / PTFE white		45° shore A	1.5 mm	100	702981
Silicone white / PTFE blue		55° shore A	1.5 mm	100	702110

Containers for screw neck vials N 18 and crimp neck vials N 20



Description	Pack of	REF
25 position container blue, with removable divider for headspace screw neck vials N 18 and crimp neck vials N 20; outer length 130 mm, outer width 130 mm, outer height 80 mm, with transparent lid (suitable for freezers)	1	702516



















Crimp neck vials and closures N 20

Crimp neck vials and closures N 20



- Large range of Headspace crimp neck vials with different volumes and diameters
 - Flat DIN crimp neck with stable bearing surface for the septum (especially suited for high vial pressures) as well as beveled HS crimp neck for instruments of certain manufacturers (PerkinElmer).
 - Assignment to respective instrument manufacturers in parentheses
 - Different types of crimp closures depending on instrument and application
- Please consider our various crimping tools on page 100 and 101.

Description	Dimensions (scale 1:2)		Pack of	REF
Headspace crimp neck vials N 20 (volume 5–10 mL)				
 70204.36		 70215.36	 702917	 702020
 70205.36		 70216.36	 702918	 702924
 38 20.5		 38 20.5	 38.25 22	 38.25 21.7
 54.5 20.5		 54.5 20.5	 46 22.5	 46 22.5
Type of vial	Usable volume	OD x height		
Clear, flat bottom, flat DIN crimp neck (Varian)	5 mL	20.5 x 38 mm	100	70204.36
Amber, flat bottom, flat DIN crimp neck (Varian)	5 mL	20.5 x 38 mm	100	70215.36
Clear, rounded bottom, beveled HS crimp neck (PerkinElmer)	6 mL	22.0 x 38.25 mm	100	702917
Clear, flat bottom, beveled HS crimp neck (Metrohm, Karl–Fischer titration)	5 mL	21.7 x 38.25 mm	100	702020
Clear, flat bottom, flat DIN crimp neck (Varian)	10 mL	20.5 x 54.5 mm	100	70205.36
Amber, flat bottom, flat DIN crimp neck (Varian)	10 mL	20.5 x 54.5 mm	100	70216.36
Clear, flat bottom, flat DIN crimp neck (Dani, Agilent)	10 mL	22.5 x 46 mm	100	702918
Clear, rounded bottom, beveled HS crimp neck (CTC)	10 mL	22.5 x 46 mm	100	702924

Crimping tools N 20

Manual ergonomic crimper for 20 mm aluminium crimp caps	1	735220
Manual ergonomic decapper for 20 mm aluminium crimp caps	1	735320
Manual crimper, height adjustable, for 20 mm aluminium crimp caps	1	735120
Manual crimper, height adjustable, for 20 mm Flip Top / Flip Off caps	1	735132
Manual decapper for 20 mm aluminium crimp caps	1	735920

For electronic crimping tools N 20 see page 101; pneumatic crimping tools are available on request.












Crimp neck vials and closures N 20













Description	Dimensions (scale 1:2)	Pack of	REF
Crimp neck vials N 20 (volume 20 mL)			
Type of vial	Usable volume	OD x height	
Clear, flat bottom, flat DIN crimp neck	20 mL	23.25 x 75.5 mm	100 70206.36
Amber, flat bottom, flat DIN crimp neck	20 mL	23.25 x 75.5 mm	100 70217.36
Clear, rounded bottom, beveled HS crimp neck (PerkinElmer)	20 mL	23 x 75.5 mm	100 70254
Amber, rounded bottom, beveled HS crimp neck (PerkinElmer)	20 mL	23 x 75.5 mm	100 702541
Clear, rounded bottom, beveled HS crimp neck, label (PerkinElmer)	20 mL	23 x 75.5 mm	100 702540
Clear, flat bottom, flat DIN crimp neck (Dani, Agilent)	20 mL	22.5 x 75.5 mm	100 702261
Clear, rounded bottom, flat DIN crimp neck (CTC)	20 mL	22.5 x 75.5 mm	100 702263
Clear, rounded bottom, beveled HS crimp neck, graduation at 15 mL	20 mL	22 x 75 mm	100 702086
Crimp neck vials N 20 (volume > 20 mL)			
Type of vial	Usable volume	OD x height	
Clear, flat bottom, flat DIN crimp neck	25 mL	30 x 65 mm	100 70210.36
Clear, flat bottom, flat DIN crimp neck	50 mL	31 x 101 mm	100 70208.36
Amber, flat bottom, flat DIN crimp neck	50 mL	31 x 101 mm	100 702100
Clear, flat bottom, flat DIN crimp neck (3 rd hydrolytic class)	100 mL	51.6 x 94.5 mm	88 70209.1








Crimp neck vials and closures N 20

Cap description		Septa description					Hardness	Thickness	Pack of	REF
Center hole caps										
	with assembled septum	 702773	 702775	 70234.9	 70234	 702056	 70237	 702093	 702094	no liner 702804
N 20 aluminium crimp cap, silver, center hole				Butyl red / PTFE gray			50° shore A	3 mm	100	702773
as above				Butyl light gray / PTFE dark gray			50° shore A	3 mm	100	702775
as above				Molded septum Butyl / PTFE gray			50° shore A	3 mm	100	70234.9
as above				Butyl dark gray / PTFE gray*			50° shore A	3 mm	100	70234
N 20 aluminium crimp cap, gold , center hole				 Butyl dark gray / PTFE gray*			50° shore A	3 mm	100	702056
N 20 aluminium crimp cap, silver, center hole				Butyl stopper gray, unassembled (separate parts)			37° shore A	–	100 each	70237
as above				Silicone blue transparent / PTFE colorless			40° shore A	3 mm	100	702093
as above				Silicone white / PTFE beige			40° shore A	3 mm	100	702094
as above				no liner			–	–	100	702804
N 20 aluminium crimp cap, gold , center hole				 no liner			–	–	100	702112






Pressure release caps									
	with assembled septum	 702836	 702829	 70234.8	 702071	 702927	 702835	no liner 702799	
N 20 aluminium pressure release cap, silver, center hole		Butyl red / PTFE gray				50° shore A	3 mm	100	702836
as above		Butyl light gray / PTFE dark gray				50° shore A	3 mm	100	702829
as above		Molded septum Butyl / PTFE gray				50° shore A	3 mm	100	70234.8
as above		Butyl dark gray / PTFE gray*				50° shore A	3 mm	100	702071
as above		Silicone blue transparent / PTFE colorless				40° shore A	3 mm	100	702927
as above		Silicone white / PTFE beige				40° shore A	3 mm	100	702835
as above		no liner				-	-	100	702799







Center tear off caps								
 70233	 70236	 70236.1	no liner					
N 20 aluminium center tear off cap, gold	Butyl dark gray / PTFE gray*		50° shore A	3 mm	100	70233		
N 20 aluminium center tear off cap, silver	Butyl stopper gray, unassembled (separate parts)		37° shore A	–	100	70236		
N 20 aluminium center tear off cap, silver	no liner		–	–	100	70236.1		








Complete tear off caps								
 70235		 70238		 702805	no liner			
N 20 aluminium complete tear off cap, silver	Butyl dark gray / PTFE gray*			50° shore A	3 mm	100	70235	
N 20 aluminium complete tear off cap, silver	Butyl stopper gray, unassembled (separate parts)			37° shore A	–	100 each	70238	
N 20 aluminium complete tear off cap, silver	no liner			–	–	100	702805	

* only centrically laminated, typically called Pharma-Fix



Cap description	Septa description	Hardness	Thickness	Pack of	REF
Bi-metal crimp caps					
 with assembled septum	 702838	 702834	 702837	 no liner 702833	
N 20 bi-metal crimp cap, blue / silver, center hole	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100	702838
as above	Silicone blue transp. / PTFE colorless	40° shore A	3 mm	100	702834
as above	Silicone white / PTFE beige	40° shore A	3 mm	100	702837
as above	no liner	-	-	100	702833





Magnetic crimp caps					
 with assembled septum	 702774	 702928	 702928.9	 702929	 no liner 702808
N 20 magnetic crimp cap, silver, 8 mm center hole	Butyl red / PTFE gray	50° shore A	3 mm	100	702774
as above	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100	702928
as above	Butyl dark gray / PTFE gray*	50° shore A	3 mm	100	702928.9
as above	Silicone blue transp. / PTFE colorless	40° shore A	3 mm	100	702929
as above	no liner	-	-	100	702808



N 20 septa for crimp caps N 20					
Material	Illustration	Hardness	Thickness	Pack of	REF
Butyl red / PTFE gray		50° shore A	3 mm	100	70277
Butyl light gray / PTFE dark gray		50° shore A	3 mm	100	702057
Molded septum Butyl / PTFE gray		50° shore A	3 mm	100	702101
Butyl dark gray / PTFE gray*		50° shore A	3 mm	100	702D20TB
Silicone blue transparent / PTFE colorless		40° shore A	3 mm	100	702780
Silicone white / PTFE beige		40° shore A	3 mm	100	70278
Silicone white / Aluminium foil silver		50° shore A	3 mm	100	70279

Stopper N 20					
Butyl gray		37° shore A	-	100	702931
Bromobutyl red		45° shore A	-	100	702931.1

* only centrally laminated, typically called Pharma-Fix

PE caps N 20

Description	Illustration	Hardness	Thickness	Pack of	REF				
PE caps N 20									
height 8.4 mm		70266		702128	height 9.1 mm		70267		702129
N 20 PE cap, transparent, for beveled HS crimp neck N 20, 4.3 mm center hole (no liner)				100	70266				
as above, but with assembled septum natural rubber red-orange / TEF colorless, 45° shore A, 1.3 mm				100	702128				
N 20 PE cap, transparent, for flat DIN crimp neck N 20, 4.3 mm center hole (no liner)				100	70267				
as above, but with assembled septum natural rubber red-orange / TEF colorless, 45° shore A, 1.3 mm				100	702129				

N 19 septa for PE caps N 20					
Butyl beige / PTFE gray		55° shore A	1.3 mm	100	70269
Natural rubber red-orange / TEF colorless		45° shore A	1.3 mm	100	702904



Crimping tools

Manual crimping tools

Advanced ergonomic version



- ✦ **Available for 11 mm and 20 mm crimp caps**
 - More lightweight than complete steel crimpers
 - Ergonomically designed handles
 - Adjustment by a knob on the crimping head that is easily accessible and visible
 - Activated by bottom handle motion only which allows a steadier and safer hold of the tool during crimping
 - Due to design and alignment of the crimping head better vertical clearance over the vial
- ✦ Advanced ergonomic decappers allow safe removal of caps; no adjustment required

Standard version



- ✦ **Available for 8, 11, 13, and 20 mm crimp caps**
 - Adjustable crimping height via hexagon key, which allows to move the inner part of the crimping head up and down (not possible for manual crimpers N 8)
 - Crimping pressure adjustable via screw in the handle
 - Manual crimpers for N 13 and N 20 Flip Top / Flip Off caps (pharmaceutical closures) available
 - Long life time and convenient handling
- ✦ Manual decappers (standard version) allow safe removal of caps; no adjustment required

Description	Pack of	REF
Manual crimpers (ergonomic)		
(crimping pressure adjustable by knob on the crimping head)		
Manual ergonomic crimper for 11 mm aluminium crimp caps	1	735211
Manual ergonomic crimper for 20 mm aluminium crimp caps	1	735220
Manual decappers (ergonomic)		
Manual ergonomic decapper for 11 mm aluminium crimp caps	1	735311
Manual ergonomic decapper for 20 mm aluminium crimp caps	1	735320
Manual crimpers (standard)		
Manual crimper for 8 mm aluminium crimp caps	1	735126
Manual crimper, height adjustable, for 11 mm aluminium crimp caps	1	735111
Manual crimper, height adjustable, for 13 mm aluminium crimp caps	1	735113
Manual crimper, height adjustable, for 13 mm Flip Top / Flip Off caps	1	735133
Manual crimper, height adjustable, for 20 mm aluminium crimp caps	1	735120
Manual crimper, height adjustable, for 20 mm Flip Top / Flip Off caps	1	735132
Manual decappers (standard)		
Manual decapper for 8 mm aluminium crimp caps	1	735408
Manual decapper for 11 mm aluminium crimp caps	1	735911
Manual decapper for 13 mm aluminium crimp caps	1	735913
Manual decapper for 20 mm aluminium crimp caps	1	735920



Electronic crimping tools

Battery-powered electronic crimping tools

for 11 mm and 20 mm aluminium crimp caps (not suitable for 20 mm magnetic / bi-metal crimp caps)



- Mobile tools for consistent and reproducible crimping results
- Crimping pressure adjustable by pushing +/- buttons of the control unit on top of the tool
- Long lasting lithium ion cell batteries (full battery charge for several hundred vials, life time of battery > 1500 charges)
- CE certificate of conformity along with one year warranty
- One tool each necessary for crimping and for decapping

Electronic high power crimping tool

for 11 mm and 20 mm crimp caps (also suitable for magnetic / bi-metal crimp caps)



- Due to a more powerful motor also suitable for magnetic and bi-metal crimp caps
- Fixed power supply
- Exchangeable crimping / decapping heads
- Digital LED display of crimp settings; different jaw settings can be stored in separate programs
- CE certificate of conformity along with one year warranty
- For more convenient handling a stand is optionally available

Description	Pack of	REF
Electronic crimpers (battery-powered)		
Electronic crimper for 11 mm aluminium crimp caps	1	735511
Electronic crimper for 20 mm aluminium crimp caps (not suitable for magnetic / bi-metal crimp caps)	1	735520
Electronic decappers (battery-powered)		
Electronic decapper for 11 mm aluminium crimp caps	1	735611
Electronic decapper for 20 mm aluminium crimp caps (not suitable for magnetic / bi-metal crimp caps)	1	735620
Accessories for battery-powered electronic crimping/decapping tools		
Replacement battery 6.4 Volt, 8.6 Wh		735500

Electronic high power crimping tool		
Electronic high power crimping tool with power supply (exchangeable crimping / decapping heads separately available)	1	735700
Accessories for 735700		
Crimping head for 11 mm crimp caps (for electronic high power crimping tool 735700)	1	735711
Crimping head for 20 mm crimp caps (aluminium, magnetic, bi-metal) (for electronic high power crimping tool 735700)	1	735720
Decapping head for 11 mm crimp caps (for electronic high power crimping tool 735700)	1	735811
Decapping head for all 20 mm crimp caps (for electronic high power crimping tool 735700)	1	735820
Stand for electronic crimping tools	1	735501







Screw neck vials and closures N 24

Screw neck vials and closures N 24 (EPA)












- Recommended for VOC and TOC analyses
- Closed top screw closures for sample storage
- Most frequently used: 40 mL clear glass
- Often called EPA vials, since they are defined in the regulations of the US Environmental Protection Agency
- Due to their size mainly used as bonded closure for a firm fit of the septum
- Recommended for environmental analysis: screw closure with center hole and Silicone / PTFE septum

Description	Dimensions (scale 1:2)		Pack of	REF
Screw neck vials N 24 (EPA)				
 702021 / 702022	 702132 / 702133	 702023 / 702024	 702074 / 702131	
Type of vial	Usable volume	OD x height	Pack of	REF
Clear, flat bottom	20 mL	27.5 x 57 mm	100	702021
Amber, flat bottom	20 mL	27.5 x 57 mm	100	702022
Clear, flat bottom	30 mL	27.5 x 72.5 mm	100	702132
Amber, flat bottom	30 mL	27.5 x 72.5 mm	100	702133
Clear, flat bottom	40 mL	27.5 x 95 mm	100	702023
Amber, flat bottom	40 mL	27.5 x 95 mm	100	702024
Clear, flat bottom	60 mL	27.5 x 140 mm	100	702074
Amber, flat bottom	60 mL	27.5 x 140 mm	100	702131

Screw neck vials and closures N 24



Description	Dimensions (scale 1:2)				Pack of	REF
Screw closures N 24, plain screw caps N 24 and single septa N 22						
 702058	 702059	 702073	 702130	 702102	 702060	 702061
Cap description		Septum description	Hardness	Thickness		
N 24 PP bonded* screw cap, white, center hole		Silicone white / PTFE beige	45° shore A	3.2 mm	100	702058
as above, but closed top		Silicone white / PTFE beige	45° shore A	3.2 mm	100	702059
N 24 PP bonded* screw cap, white, center hole		Red Rubber / TEF colorless	65° shore A	2.5 mm	100	702073
* Septum firmly connected with the cap, cannot be removed						
N 24 PP screw cap, white, center hole		Butyl red / PTFE gray	50° shore A	2.4 mm	100	702130
as above, but closed top		Butyl red / PTFE gray	50° shore A	2.4 mm	100	702102
N 24 PP screw cap, white, center hole		no liner			100	702060
as above, but closed top		no liner			100	702061
N 22 septum, Silicone natural / PTFE beige			45° shore A	3.2 mm	100	702062
N 22 septum, Butyl red / PTFE gray			50° shore A	2.4 mm	100	702791

Autosampler compatibility

Agilent

Application / Type of vial	Most popular MN products fur use on Agilent instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70282, 70286		70289	78
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079, 702007, 702008, 702088	702813, 702818, 702825, 702077	702732, 702082, 702081, 702080, 702287.1, 702038, 702037, 702035, 702085, 702084, 702028, 702026	80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076, 702891, 702014, 702015, 702016	702813, 702818, 702825, 702077	70256, 70231.3, 70231.1, 70231.2, 702001, 702730, 70288, 702995, 702879 (for GC PAL)	84
HPLC:				
N 9 screw (standard samples)	As indicated under GC, but additionally the following closures with slit septum: 702288.1, 702083, 702040, 702027			80
N 11 crimp (standard samples)	As indicated under GC, but additionally the following closure with slit septum: 702823			84
N 11 snap ring (standard samples)	702714, 702713, 702712, 702709	702813, 702818, 702825, 702077	702063, 702731, 702710.1, 702064, 702718, 702401	87
Headspace:	Vials:		Closures:	
N 18 screw (Combi PAL + G 1888A)	702866, 702826		702055	95
N 20 crimp	702918, 702261, 702263		70234, 702094, 702093, 702071, 702835, 702927	96



Autosampler compatibility

CTC

Application / Type of vial	Most popular MN products for use on CTC instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70282, 70286, 70212, 70212.1, 702002, 702003		70289, 702878	78
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079	702813, 702818, 702825, 702077	702287.1, 702035, 702026, 702027	80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076	702813, 702818, 702825, 702077	702879 (for GC PAL), 70288, 702995	84
HPLC:				
N 9 screw (standard samples)	As indicated under GC, but additionally the following closures with slit septum: 702288.1, 702027			80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076	702813, 702818, 702825, 702077	70288, 702995, 702823	84
N 11 snap ring (standard samples)	702714, 702713, 702712	702813, 702818, 702825, 702077	702710.1, 702717.2, 702718.1	87
Headspace:				
N 18 screw (Combi PAL)	702866, 702826		702055, 702827	95
N 20 crimp	702924, 702263		702929, 702834	96

Dionex

Application / Type of vial	Most popular MN products for use on Dionex instruments			Page
HPLC:	Vials:	Inserts:	Closures:	
N 8 screw (microsampling)	702880, 70286, 70282		702025, 70289	78
N 8 screw (standard samples)	70213, 70213.2, 702004, 702893, 702860	702968, 702824, 702005	70245, 702437	79
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079, 702007, 702008	702813, 702818, 702825, 702077	702287.1, 702288.1, 702026, 702027	80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076, 702891, 702014	702813, 702818, 702825, 702077	70288, 702823, 70256	84
N 11 snap ring (standard samples)	702714, 702713, 702712, 702709	702813, 702818, 702825, 702077	702710.1, 702717.2	87
IC:				
N 9 screw	702009		702288.1, 702027	80



PerkinElmer

Application / Type of vial	Most popular MN products for use on PerkinElmer instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70251, 70286		70252.1, 70283, 702025	78
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702078, 702079	702818, 702825	702732, 702287.1, 702026, 702028	80
N 10 screw (standard samples)	702012, 702013	702818, 702825	702045, 702046	83
N 11 crimp (standard samples)	70201CG*, 70214CG*	702824*, 702005*	70256, 702730, 70231.1, 70231.2, 70231.3, 70288, 702995	84
(* small opening; ** wide opening)	70201HP**, 702885**, 702892**, 702888**, 702075**, 702076**	702818**, 702825**		84
HPLC:				
N 8 crimp (microsampling)	70286		70252.1, 702025	78
N 9 screw (standard samples)	As indicated under GC, but additionally the following closures with slit septum: 702288.1, 702027			80
N 10 screw (standard samples)	As indicated under GC, but additionally the following closure with slit septum: 702047			83
N 11 crimp (standard samples)	As indicated under GC, but additionally the following closure with slit septum: 702823			84
N 11 snap ring (standard samples)	702714, 702713, 702712	702818, 702825	702064, 702718, 702710, 702401	87
Headspace:				
N 18 screw (CTC Combi PAL + TurboMatrix™ HS 16 + 40)	702866, 702826		702055, 702827, 702072	95
N 20 crimp (CTC Combi PAL)	702924, 702263		702929, 702834, 702928.9, 702928, 702774	96
N 20 crimp (TurboMatrix™ HS 16, 40 + 110) *** not suited for TurboMatrix™ 110	702917***, 70254, 702540, 702541		702829, 702836, 702071, 70234.8, 702835, 702927, 702775, 702773, 70234, 70234.9, 702093, 702094, 70237, 702931	96



Autosampler compatibility

Vials and accessories

Shimadzu

Application / Type of vial	Most popular MN products for use on Shimadzu instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70282, 70286, 70212, 70212.1, 702002, 702003		70289, 702878	78
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079, 702007, 702008	702813, 702818, 702825, 702077	702287.1, 702035, 702026	80
N 10 screw (standard samples)	702011, 702012, 702013	702813, 702818, 702825, 702077	702045, 702046, 702048	83
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076, 702891, 702014	702813, 702818, 702825, 702077	702879 (for AOC5000), 70288, 702995	84
N 13 screw (large volume samples)	702962, 702973, 702089	702972 + spring 702974	702926	90
HPLC:				
N 8 crimp (microsampling)	70282, 70286		702025, 70289, 702878	78
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079, 702009, 702007, 702008	702813, 702818, 702825, 702077	702287.1, 702037, 702038, 702036, 702026, 702288.1, 702040, 702083, 702039, 702027, 702031	80
N 10 screw (standard samples)	702011, 702012, 702013	702813, 702818, 702825, 702077	702045, 702046, 702047	83
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702141, 702075, 702076, 702891, 702014	702813, 702818, 702825, 702077	702730, 70288, 702823	84
N 11 snap ring (standard samples)	702714, 702713, 702712, 702809	702813, 702818, 702825, 702077	702064, 702717.2, 702710, 702710.1	87
N 8 + N 11 shell vials (standard samples)		Vials + closures:	70202.1 + 702807, 702017 + 702807, 702018 + 702054	94
Headspace:				
N 18 screw (AOC 5000)	702866, 702826		702055, 702827	95
N 20 crimp (AOC 5000)	702924, 702263		702929, 702834, 702928, 702774	96
N 20 crimp (HT200H)	702918, 702263		702094, 702093	96

Thermo Scientific

Application / Type of vial	Most popular MN products for use on Thermo Scientific instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70251, 70282, 70286, 702880, 70212, 70212.1, 702002, 702003		70252.1, 702025, 70289, 702878	78



Application / Type of vial	Most popular MN products for use on Thermo Scientific instruments			Page
N 8 screw (standard samples)	70213, 70213.2, 702004, 702893, 702860	702968, 702968.1, 702824, 702005	702067, 70245, 702069	79
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079, 702007, 702008	702813, 702818, 702716, 702825, 702077	702732, 702287.1, 702035, 702026	80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702141, 702075, 702076, 702891, 702014	702813, 702818, 702716, 702825, 702077	702879 (GC PAL), 702001, 70256, 702730, 70288, 702995	84
HPLC:				
N 8 crimp (microsampling)	70282, 70286, 702880, 70212, 70212.1, 702002, 702003		70252.1, 702025, 70289, 702878	78
N 8 screw (standard samples)	As indicated under GC			79
N 9 screw (standard samples)	As indicated under GC			80
N 11 crimp (standard samples)	As indicated under GC, however, not closure 702879			84
N 11 snap ring (standard samples)	702714, 702713, 702712, 702809	702813, 702818, 702716, 702825, 702077	702064, 702717.2, 702710, 702710.1, 702731, 702063	87
Headspace:				
N 18 screw (Combi PAL)	702866, 702826		702055, 702827	95
N 20 crimp (Combi PAL)	702924, 702263		702929, 702834	96
N 20 crimp (HS850/HS200)	702924, 702263		702094, 702093, 702773 / 702775 / 70234.9	96

Varian (now Agilent)

Application / Type of vial	Most popular MN products for use on Varian instruments			Page
GC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70282, 70286, 702880, 70212, 70212.1, 702002, 702003		70289, 702878	78
N 8 screw (standard samples)	70213, 70213.2, 702004, 702893, 702860	702968, 702824	702067, 70245, 702069	79
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702078, 702079	702813, 702818, 702716, 702825, 702077	702732, 702287.1, 702037, 702035, 702026	80
N 11 crimp (standard samples)	70201HP, 702885, 702892, 702888, 702075, 702076	702813, 702818, 702716, 702825, 702077	702879 (GC PAL), 70256, 70288, 702995	84
HPLC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	As indicated under GC, but additionally closures 70252.1, 702025			78
N 8 screw (standard samples)	As indicated under GC, but additionally closure 702437			79
N 9 screw (standard samples)	As indicated under GC, but additionally closure 702040			80



Autosampler compatibility

Varian (now Agilent)

Application / Type of vial	Most popular MN products for use on Varian instruments			Page
N 10 screw (standard samples)	702011, 702012, 702013	702813, 702818, 702716, 702825, 702077	702045 / 702046, 702048, 702044	83
N 11 crimp (standard samples)	As indicated under GC, however, not closure 702879			84
N 11 snap ring (standard samples)	702714, 702713, 702712, 702809	702813, 702818, 702716, 702825, 702077	702064, 702717.2, 702710, 702710.1, 702731, 702063, 702718	87
Headspace:				
N 18 screw (Combi PAL)	702866, 702826		702055, 702827, 702073	95
N 20 crimp (Combi PAL)	702924, 702263		702929, 702834	96
N 20 crimp (CP-9020/9025, CP-9060, Genesis)	702924, 702918, 702261		702093, 70234, 702773 / 702775	96

VWR (Merck® / Hitachi)

Application / Type of vial	Most popular MN products for use on VWR instruments			Page
HPLC:	Vials:	Inserts:	Closures:	
N 8 crimp (microsampling)	70286, 70282		70289, 702878	78
N 8 screw (standard samples)	70213, 70213.2, 702004, 702893, 702860	702968, 702968.1, 702824, 702005	702067, 70245, 702437	79
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702078, 702079	702813, 702716, 702818, 702077	702287.1, 702288.1, 702026, 702027, 702031	80
N 11 snap ring (standard samples)	702714, 702713, 702712	702813, 702716, 702818, 702077	702063, 702710.1, 702717.2	87
N 13 screw (large volume samples)	702962, 702973	702972 + spring 702974	702926, 702963 + 70260	90

Waters®

Application / Type of vial	Most popular MN products for use on Waters® instruments			Page
HPLC:	Vials:	Inserts:	Closures:	
N 9 screw (standard samples)	702282, 702293, 702283, 702284, 702006, 702007, 702008, 702078, 702079, 702009	702818	702026, 702027, 702287.1, 702288.1, 702037, 702040, 702038, 702083, 702287, 702288	80
N 10 screw (standard samples)	702011, 702012, 702013	702818	702045, 702046, 702047	83
N 11 snap ring (standard samples)	702714, 702713, 702712, 702809, 702709	702818	702710.1, 702717.2	87
N 8 shell vials (standard samples)		Vials + closures:	70202.1 + 702807, 702017 + 702807	94
N 13 screw (large volume samples)	702962, 702973, 702089	702972 + spring 702974	702926, 702963 + 70260	90