



MASS SPEC



Compatible w/MODELS:
3200, 4000, 4500, 5000,
5500, 6500

33-1034
OEM 025392

TIS Capillary Electrode



Compatible w/MODELS:
2000, QTRAP™

CTS-21212
OEM 018782

TIS Capillary Electrode



Compatible w/MODELS:
4500, 5500, 6500

CTS-21309
OEM 4465678

**Rotor Seal, Valco
Cheminert**



Compatible w/MODELS:
3200, 3500, 4000, 4500,
5000, 5500, 6500

CTS-21134
OEM 027471

**1/16" Short Hex
PEEK™ Fitting**



Compatible w/MODELS:
3200, 3500, 4000, 4500,
5000, 5500, 6500

CTS-21173
OEM 003403

**O-Ring
(3/16" ID X 1/16", Viton™)**



Compatible w/MODELS:
3200, 3500, 4000, 4500,
5000, 5500, 6500

CTS-21137
OEM 025348, 5045984

**1/16" Union for TIS Probe/
APCI for Turbo V™ Source**



Compatible w/MODELS:
4000, 4500, 5000, 5500,
6500

CTS-21523
OEM 017809

**O-Ring, Viton™
(5.3 ID X 1.77 W)**



Compatible w/MODELS:
4000, 3200, 3500, 4500,
5500, 6500

CTS-21500 **KIT**
OEM 5058491

Electrode Turbo Kit, MS


- Peek™ Fitting-Natural (CTS-21446) 1
- Spring, Turbo V™ Source (CTS-21189) 1
- O-Ring (3/16" ID X 1/16", Viton™) (CTS-21173) 1
- 1/16" Union for TIS Probe and APCI (CTS-21137) 1
- TIS Capillary Electrode (33-1034) 5



Compatible w/MODELS:
4000, 4500, 5000, 5500,
6500

CTS-21189
OEM 026626

Spring, Turbo V™ Source



Compatible w/MODELS:
4000, 4500, 5000, 5500,
6500

CTS-21102
OEM 025388

APCI Capillary Electrode

MASS SPEC



Compatible w/MODELS: 3200, 3500, 4000, 4500, 5500, 6500

CTS-21286 KIT

OEM 5058495,
(Old# 5043951)

PM KIT, MS

- 1/16" Short Hex PEEK™ Fitting (CTS-21134) 1
- Spring, Turbo V™ Source (CTS-21189) 1
- O-Ring (3/16" ID X 1/16", Viton®) (CTS-21173) 1
- TIS Capillary Electrode (33-1034) 2
- 1/16" Union for TIS Probe and APCI (CTS-21137) 1
- Fingertight One Piece 10-32 PEEK™ Male Nut/Ferrule for 1/16" OD Tubing (CTS-1818) 2
- Tubing PEEK™ Red, 1/16" OD X 0.005" ID X 31.5 in (CTS-21172) 1



Compatible w/MODELS: 4600, 5600

CTS-21573

OEM 1032854

Air Intake Filter, 4600/5600 Systems




Compatible w/MODELS: 4500, 5500, 6500

CTS-21459

OEM 1027708

Air Intake Filter, 4/pk



Compatible w/MODELS: 3000, 4000, 5000

CTS-21444

OEM 1003712

Air Intake Filter-Yellow



Compatible w/MODEL: LC-20AD

CTS-21771

OEM 4443034

PM Kit, LC-20AD



Compatible w/MODEL: 3200

CTS-21572

OEM 1017396

Air Intake Filter



Compatible w/MODEL: LC-20AD/AB

CTS-21822

OEM 4444114

PM Kit, LC-20AD/AB



Compatible w/MODELS: 2000, 3000, 3200, 3200 QTRAP™, 4000, 4000 QTRAP™, QTRAP™ 4500, QTRAP, 4600, 5000, 5500+, 5600, 6600, QTRAP™ 6500

CTS-21936

OEM 1017396

Foam Swab, Anti-Static (50/pk)



Compatible w/MODEL: LC-20ADXR

CTS-21823

OEM 4448440

PM Kit, LC-20ADXR

MASS SPEC




Compatible w/MODELS: 4500, 5500, 6500

CTS-21568 KIT

OEM 5054105

PM Kit (4500, 5500, 6500)

- | | | | |
|---|---|--|----|
| • Air Intake Filter, 4/pk (CTS-21459) | 1 | • Spring, Turbo V™ Source (CTS-21189) | 1 |
| • PEEK™ Tubing-Red 1/6" OD X .005" ID X 19.685" (CTS-21592) | 1 | • O-Ring Viton® Brown .208 ID X .070" Thick (CTS-21523) | 2 |
| • Label, Field Service-Sciencix (CTS-21447) | 1 | • O-Ring .750" ID X .093" Thickness (CTS-21569) | 2 |
| • Fingertight One Piece 10-32 PEEK™ Male Nut/Ferrule for 1/16" OD Tubing (CTS-1818) | 2 | • Anti-Static Foam Swab-White Large Square Tip 6" (CTS-21595) | 20 |
| • 1/16" Short Hex PEEK™ Fitting (CTS-21134) | 1 | • TIS Capillary Electrode (33-1034) | 1 |
| • 1/16" Union for TIS Probe and APCI (CTS-21137) | 1 | • Anti-Static Foam Swab-Blue Small Round Tip 6.25" (CTS-21594) | 20 |
| • O-Ring (3/16" ID X 1/16", Viton®) (CTS-21173) | 1 | | |



Compatible w/MODELS: 4000, 5000, 4000 QTRAP™

CTS-21445 KIT

OEM 028447

PM Kit, 4000, 5000

- | | | | |
|--|----|---|---|
| • O-Ring Viton® Brown .208 ID X .070" (CTS-21523) | 2 | • Spring, Turbo V™ Source (CTS-21189) | 1 |
| • PEEK™ Tubing-Red 1/6" OD X .005" ID X 19.685" (CTS-21592) | 1 | • O-Ring (3/16" ID X 1/16", Viton®) (CTS-21173) | 1 |
| • Anti-Static Foam Swab-Blue Small Round Tip 6.25" (CTS-21594) | 20 | • 1/16" Union for TIS Probe and APCI (CTS-21137) | 1 |
| • Anti-Static Foam Swab-White Large Square Tip 6" (CTS-21595) | 10 | • 1/16" Short Hex PEEK™ Fitting (CTS-21134) | 1 |
| • Label, Field Service-Sciencix (CTS-21447) | 1 | • Fingertight One Piece 10-32 PEEK™ Male Nut/Ferrule for 1/16" OD Tubing (CTS-1818) | 2 |
| • Air Intake Filter-Yellow (CTS-21444) | 1 | • TIS Capillary Electrode (33-1034) | 1 |

Distribútor pre SR:

AZ CHROM s.r.o., Robotnícka 10, 831 03 Bratislava
 Tel. 0907 244526 azetchrom@hplc.sk, www.azetchrom.sk