

PUMP PARTS



MODEL: LC-10AS, LC-7A

CTS-10319
OEM 228-17019-93

Sapphire Plunger



MODEL: LC-10ATvp

CTS-10584
OEM 228-35009-92

Ceramic Plunger Assembly



MODEL: LC-600, LC-9A, LC-10AD

CTS-10320
OEM 228-18523-91

Sapphire Plunger



MODEL: LC-20AT

CTS-11425
OEM 228-35009-93

Sapphire Plunger



MODEL: LC-20AD/AB, LC-20ADXR, LC-10ADvp, LC-2010, LC-30ADSF

CTS-11297
OEM 228-35601-93

Sapphire Plunger



MODEL: LC-6A

CTS-0422
OEM 228-11999-00

Plunger Seal



MODEL: LC-20AD/AB, LC-20ADXR, LC-10ADvp, LC-2010, LC-30ADSF

11-3795
OEM 228-35281-95

Sapphire Plunger with Holder



MODEL: LC-7A, LC-10Ai, LC-10AS, LC-10AT, LC-10ATvp

CTS-10340
OEM 228-21975-00

Plunger Seal



MODEL: LC-10ADvp

CTS-A138
OEM 228-35601-91

Ceramic Plunger Assembly



MODEL: LC-600, LC-9A, LC-10AD, LC-10ADvp, LC-20AD/AB, LC-2010 A/C HT

CTS-10693
OEM 228-35146-00
Old# 228-18745-00

Plunger Seal



MODEL: LC-10ADvp, LC-20AD, LC-20AB, LC-20ADXR, LC-30AD

11-3897
OEM 228-35602-91

Plunger Holder



MODEL: LC-10ATvp, LC-20AT, LC-2010 A/C HT, SIL-20A/AC, SIL-2010 A/C HT, SIL-10ADvp, LC-HT, SIL-30AC, SIL-30ACMP

CTS-10694
OEM 228-35145-00

Plunger Seal

PUMP PARTS



MODEL: LC-10Ai, LC-10AS,
LC-10AT, LC-10ATvp

CTS-10583
OEM 228-28499-00
Plunger Rinse Seal



MODEL: LC-600,
LC-9A, LC-10AD

CTS-A102
OEM 228-33492-91
Old# 228-18522-91
Inlet Check Valve



MODEL: LC-10ADvp,
LC-20AD/AB

CTS-10963
OEM 228-32628-00
Gold Plunger Seal



MODEL: LC-600, LC-9A,
LC-10AD, LC-10AT

CTS-A103
OEM 228-32531-92
Old# 228-18522-92
Outlet Check Valve



MODEL: LC-10AD, LC-10ADvp,
LC-2010, LC-2010 HT,
LC-20AD/AB, LC-30ADSF

11-2891
OEM 228-32784-91
Old# 228-24311-01,
228-31828-00
PTFE Diaphragm, 2/pk



MODEL: LC-10ADvp,
10ATvp, 20AD/B, 2010HT

CTS-A11316
OEM 228-45707-91
**Suction Filter,
Stainless Improved**



MODEL: LC-6A, LC-10AS

CTS-A099
OEM 228-12353-91
Inlet Check Valve



MODEL: LC-10ATvp,
LC-10AT

CTS-A117
OEM 228-32166-91
Inlet Check Valve



MODEL: LC-6A, LC-8A,
LC-10AS

CTS-A100
OEM 228-09054-93
Outlet Check Valve



MODEL: LC-10ADvp

CTS-A145
OEM 228-39093-92
Inlet Check Valve



MODEL: LC-6A, LC-10AS

CTS-10595 KIT
OEM 228-11200-91
Check Valve Rebuild Kit
(for 2 valves)

- Ruby Balls 2
- Seats 2
- Thin Washers 2
- Thick Washers 2
- Retainer Washers 2



MODEL: LC-10ATvp,
LC-10ADvp

CTS-A121
OEM 228-34976-91
Outlet Check Valve

PUMP PARTS



MODEL: LC-20AD/AB XR,
LC-20AT, LC-30ADSF

CTS-A11186

OEM 228-48249-91
Old# 228-45704-91,
228-45557-91

Inlet Check Valve



MODEL: LC-10AD,
LC-10ADvp, LC-10ATvp,
LC-20AD/B, LC-20AT

CTS-11305

OEM 228-48607-91
Old# 228-32744-00

PEEK™ 5um Filter Assy



MODEL: LC-20AD/AB XR,
LC-10ADvp, LC-10ATvp,
i-Series, LC-30ADSF

CTS-A11187

OEM 228-45705-91
Old# 228-45563-95

Outlet Check Valve



MODEL: LC-20AD/B, LC-10ADvp,
LC-10ATvp, LC-2010

CTS-A11421

OEM 228-35871-96

Line Filter Assembly



MODEL: LC-20AT

CTS-11-3807

OEM 228-45704-92

Inlet Check Valve, Primary

AUTOSAMPLER



MODEL: SIL-10A,
SIL-10AXL, SIL-10Ai

CTS-10578
OEM 228-25237-04
Syringe, 500 μ L



MODEL: SIL-10AXL,
SIL-10ADvp, LC-2010

CTS-10686
OEM 228-33355-04
Needle Seal, Vespel[®]



MODEL: SIL-10AXL/vp,
SIL-10ADvp, SIL-HT

CTS-11435
OEM 228-33355-01
Seal, Needle PEEK[™]



MODEL: SIL-HT, SIL-20 A/C,
SIL-20HT A/C

CTS-11312
OEM 228-41310-92
PEEK[™] Rotor, HPV



MODEL: LC-2010,
LC-2010 HT, SIL-20 A/C,
SIL-20HT A/C, SIL-HT,
SIL-20 A/C XR

CTS-11313
OEM 228-36923-00
PEEK[™] Rotor, LPV



MODEL: SIL-20A/AC

CTS-11315
OEM 228-42325-01
**Needle Seal,
PEEK[™] SIL-20**



MODEL: SIL-10ADvp

CTS-10942
OEM 228-21217-97
Rotor Seal



MODEL: SIL-10A,
SIL-10AXL, SIL-10Ai

CTS-A106
OEM 228-21217-91
Rotor Seal



MODEL: SIL-10A,
SIL-10AXL, SIL-10Ai

CTS-A067
OEM 228-21220-91
Stator Assembly



MODEL: SIL-10ADvp,
SIL-20A/AC, LC-2010 A/C HT,
LC-HT, SIL-30AC,
SIL-30ACMP

CTS-11434
OEM 228-35010-91
Sapphire Plunger



MODEL: LC-10ATvp, LC-20AT,
LC-2010 A/C HT, SIL-20A/AC,
SIL-2010 A/C HT, SIL-10ADvp,
LC-HT, SIL-30AC, SIL-30ACMP

CTS-10694
OEM 228-35145-00
Plunger Seal, GFP

DETECTOR LAMPS



MODEL: SPD-10A,
SPD-10AV, SPD-10AVvp,
SPD-20AV

CTS-00135
OEM 228-34016-02
**Long Life Deuterium Lamp
(2000 hr)**