

Microsyringes for Agilent GC Autosampler

L-MARK® Syringes for Agilent GC Autosampler - Removable Needle



made
 in
 Germany

Microsyringes 0.5µl-10µl - Needle Tip Style AS, Needle Length 42mm

Part No.	Volume	Plunger Type	Needle	OD (Gauge)	Qty.	Agilent	ILS
3030151	0.5µl	SS (In Needle)	RN	0.64mm (G23s)	ea.	5188-5246	2204052
3030152	1µl	SS (In Needle)	RN	0.64mm (G23s)	ea.	5188-5247	2204103
3030129	5µl	SS	RN	0.47mm (G26s)	ea.	5182-0834	2104307
3030130	5µl	SS	RN	0.64-0.47mm (G23s-G26s)	ea.	5182-0835	2104308
3030138	10µl	SS	RN	0.64mm (G23s)	ea.	5181-8806	2104409
3030137	10µl	SS	RN	0.64-0.47mm (G23s-G26s)	ea.	5181-3321	2104408
3030166	10µl	PTFE	RN	0.64mm (G23s)	ea.	5181-8813	2604409
3030165	10µl	PTFE	RN	0.64-0.47mm (G23s-G26s)	ea.	5181-3356	2604408

L-MARK® Syringes for Agilent GC Autosampler - Fixed Needle



Microsyringes 5µl-100µl - Needle Tip Style AS, Needle Length 42mm

Part No.	Volume	Plunger Type	Needle	OD (Gauge)	Qty.	Agilent	ILS
3030128	5µl	SS	FN	0.64mm (G23s)	ea.	9301-0892	2104305
3030133	5µl	SS	FN	0.64mm (G23s)	pk.6	5182-0875	2104315
3030127	5µl	SS	FN	0.64-0.47mm (G23s-G26s)	ea.	5181-1273	2104304
3030132	5µl	SS	FN	0.64-0.47mm (G23s-G26s)	pk.6	5181-8810	2104314
3030126	5µl	SS	FN	0.47mm (G26s)	ea.	9301-0891	2104303
3030131	5µl	SS	FN	0.47mm (G26s)	pk.6	5183-4728	2104313
3030136	10µl	SS	FN	0.64mm (G23s)	ea.	9301-0713	2104405
3030141	10µl	SS	FN	0.64mm (G23s)	pk.6	9301-0725	2104415
3030135	10µl	SS	FN	0.64-0.47mm (G23s-G26s)	ea.	5181-1267	2104404
3030140	10µl	SS	FN	0.64-0.47mm (G23s-G26s)	pk.6	5181-3360	2104414
3030134	10µl	SS	FN	0.47mm (G26s)	ea.	9301-0714	2104403
3030139	10µl	SS	FN	0.47mm (G26s)	pk.6	5183-4729	2104413
3030163	10µl	PTFE	FN	0.64mm (G23s)	ea.	5181-8809	2604403
3030167	10µl	PTFE	FN	0.64mm (G23s)	pk.6	5183-4730	2604413
3030164	10µl	PTFE	FN	0.64-0.47mm (G23s-G26s)	ea.	5181-3354	2604404
3030168	10µl	PTFE	FN	0.64-0.47mm (G23s-G26s)	pk.6	5181-3361	2604414
3030169	25µl	PTFE	FN	0.64mm (G23s)	ea.	5183-0316	2604505
3030174	25µl	PTFE	FN	0.64-0.47mm (G23s-G26s)	ea.	5183-0316	2604504
3030170	50µl	PTFE	FN	0.64mm (G23s)	ea.	5183-0318	2604605
3030175	50µl	PTFE	FN	0.64-0.47mm (G23s-G26s)	ea.	5183-0314	2604604
3030171	100µl	PTFE	FN	0.64mm (G23s)	ea.	5183-2058	2604705
3030176	100µl	PTFE	FN	0.64-0.47mm (G23s-G26s)	ea.	5183-2058	2604704

Syringes 5µl-10µl - Needle Tip Style a, Needle Length 42mm

Part No.	Volume	Plunger Type	Needle	OD (Gauge)	Qty.	Hamilton	ILS
3031126	5µl	SS	FN	0.47mm (G26s)	ea.	87992	
3031127	5µl	SS	FN	0.64mm (G23s)	ea.	87991	
3031128	10µl	SS	FN ¹⁾	0.47mm (G26s)	ea.	80399	
3031129	10µl	PTFE	FN	0.64mm (G23s)	ea.	80398	2104410